

Емкости лабораторные CARL ROTH ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47

Bowls and tubs

Dish ROTILABO®

ROTH SELECTION. Material: PP.

With grooved base. Robust design. Colour: White. Suitable for foodstuffs.



Not a medical device / Not an IVD product

Vol. (l)	Upper length (mm)	Upper width (mm)	Lower length (mm)	Lower width (mm)	Height (mm)	Art. No.	Pack Qty.	
0.5	230	180	190	140	42	1722.1	1 unit(s)	
1.5	310	250	250	190	65	1724.1	1 unit(s)	
3.0	370	310	310	250	75	0772.1	1 unit(s)	
10.0	520	420	410	310	120	0778.1	1 unit(s)	
21.0	634	534	530	430	140	0782.1	1 unit(s)	
39.0	846	648	720	520	160	0799.1	1 unit(s)	

Range of dishes ROTILABO®

ROTH SELECTION. Material: PP.

Consists of six dishes, volumes: 0.5, 1.5, 3, 10, 21 and 39 l.

Not a medical device / Not an IVD product

Version	Art. No.	Pack Qty.	
ROTILABO® range of laboratory dishes	0938.1	1 set	

Dish ROTILABO®

ROTH SELECTION. Material: PVC, hard.

With grooved base. Colour: White.



Not a medical device / Not an IVD product

Upper length (mm)	Upper width (mm)	Lower length (mm)	Lower width (mm)	Height (mm)	Art. No.	Pack Qty.	
245	195	200	150	45	7032.1	1 unit(s)	
305	250	250	200	60	7034.1	1 unit(s)	
385	325	320	260	70	7035.1	1 unit(s)	
420	375	350	300	80	7038.1	1 unit(s)	
505	405	430	330	90	7041.1	1 unit(s)	
600	495	520	420	90	7043.1	1 unit(s)	
635	515	540	420	180	7047.1	1 unit(s)	
765	620	675	530	100	7045.1	1 unit(s)	



Bowls and tubs

Dish ROTILABO®

ROTH SELECTION. Material: PVC, break-proof.
Temperature stable up to +70 °C.



Not a medical device / Not an IVD product

With grooved base

Suitable as development dishes or for staining electrophoresis gels, the raised grooves in the base stop the gels from sticking to the bottom of the dish.

Upper length (mm)	Upper width (mm)	Lower length (mm)	Lower width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
180	230	140	190	40	white	C030.1	1 unit(s)	
180	230	140	190	40	red	C037.1	1 unit(s)	
240	320	190	260	50	white	C031.1	1 unit(s)	
240	320	190	260	50	red	C038.1	1 unit(s)	
240	320	190	260	50	black	C045.1	1 unit(s)	
320	370	260	320	50	white	C032.1	1 unit(s)	
320	370	260	320	50	red	C039.1	1 unit(s)	
320	370	260	320	50	black	C046.1	1 unit(s)	
430	530	310	410	105	white	C033.1	1 unit(s)	
430	530	310	410	105	red	C040.1	1 unit(s)	
520	610	420	510	110	white	C034.1	1 unit(s)	
520	610	420	510	110	red	C041.1	1 unit(s)	
520	610	420	510	110	black	C048.1	1 unit(s)	
630	730	510	610	110	white	C035.1	1 unit(s)	
630	730	510	610	110	black	C049.1	1 unit(s)	

Without grooved base

Upper length (mm)	Upper width (mm)	Lower length (mm)	Lower width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
150	200	110	160	30	white	C050.1	1 unit(s)	
150	200	110	160	30	red	C057.1	1 unit(s)	
150	200	110	160	30	black	C064.1	1 unit(s)	
180	230	140	190	40	white	C051.1	1 unit(s)	
180	230	140	190	40	red	C058.1	1 unit(s)	
180	230	140	190	40	black	C065.1	1 unit(s)	
240	320	190	260	50	white	C052.1	1 unit(s)	
240	320	190	260	50	red	C059.1	1 unit(s)	
240	320	190	260	50	black	C066.1	1 unit(s)	
320	370	260	320	50	white	C053.1	1 unit(s)	
320	370	260	320	50	red	C060.1	1 unit(s)	
320	370	260	320	50	black	C067.1	1 unit(s)	
430	530	310	410	105	white	C054.1	1 unit(s)	
430	530	310	410	105	red	C061.1	1 unit(s)	
430	530	310	410	105	black	C068.1	1 unit(s)	
520	610	420	510	110	white	C055.1	1 unit(s)	
520	610	420	510	110	red	C062.1	1 unit(s)	
520	610	420	510	110	black	C069.1	1 unit(s)	
630	730	510	610	110	white	C056.1	1 unit(s)	
630	730	510	610	110	red	C063.1	1 unit(s)	
630	730	510	610	110	black	C070.1	1 unit(s)	

Bowls and tubs

Dish ROTILABO®

ROTH SELECTION. Material: MF.

Temperature stable from -40 to +110 °C.



Delivery without lid, please order separately. **Note:** There is no lid available for Art. No. 8455.1.

Not a medical device / Not an IVD product

Vol. (l)	Upper length (mm)	Upper width (mm)	Height (mm)	Art. No.	Pack Qty.	
0.7	190	150	40	8449.1	1 unit(s)	
1	290	160	35	8453.1	1 unit(s)	
1.9	290	160	60	8454.1	1 unit(s)	
2.2	350	250	40	8455.1	1 unit(s)	
6.0	340	245	100	8457.1	1 unit(s)	

Accessories lid for laboratory dishes

ROTH SELECTION. Material: PS.

For ROTILABO® dishes made from melamine.



Not a medical device / Not an IVD product

Upper length (mm)	Upper width (mm)	Suitable for	Version	Art. No.	Pack Qty.	
190	150	Dish, Art. No. 8449.1	Button cover	P163.1	1 unit(s)	
290	160	Dishes, Art. No. 8453.1 and 8454.1	Button cover	P164.1	1 unit(s)	
340	245	Dish, Art. No. 8457.1	Grip dish lid, stackable	P165.1	1 unit(s)	

Dish ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10, polished.

Autoclavable: Yes.

Deep-drawn, with flat edge. Suitable for foodstuffs.



A 121 °C

Not a medical device / Not an IVD product

Pic.	Vol. (l)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
(1)	0.6	190	120	40	YH70.1	1 unit(s)	
(1)	1.2	240	170	40	YH71.1	1 unit(s)	
(1)	2.0	300	205	55	YH72.1	1 unit(s)	
(2)	3.0	340	210	60	YH73.1	1 unit(s)	



Bowls and tubs

Dish

Nalgene. Material: PP. Autoclavable: Yes.

Suitable for storing and sterilising instruments and small laboratory items. Heavy-duty design with reinforced edges for maximum stability.



A 121 °C



Not a medical device / Not an IVD product

Vol. (l)	Outer length (mm)	Outer width (mm)	Height (mm)	Art. No.	Pack Qty.	
2	264	161	66	ATE7.1	1 unit(s)	
2.8	262	161	104	ATE8.1	1 unit(s)	
3.1	325	264	66	ATE9.1	1 unit(s)	
5	323	264	107	ATH0.1	1 unit(s)	

Rinsing tray

Nalgene. Material: PP. Autoclavable: Yes.

Universally applicable tray for a variety of applications. Suitable for collecting, soaking and rinsing soiled laboratory items and for drying and sterilising instruments and small laboratory items.

- With integrated handles
- Robust and stackable



NEW / +

A 121 °C



Not a medical device / Not an IVD product

Vol. (l)	Outer length (mm)	Outer width (mm)	Height (mm)	Art. No.	Pack Qty.	
9.5	361	311	143	1L19.1	1 unit(s)	
20	543	435	130	1L1A.1	1 unit(s)	

Dish ROTILABO® stainless steel

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes.

Suitable for foodstuffs.



A 121 °C

Not a medical device / Not an IVD product

Upper length (mm)	Upper width (mm)	Lower length (mm)	Lower width (mm)	Height (mm)	Weight (g)	Art. No.	Pack Qty.	
370	275	350	250	50	1150	0375.1	1 unit(s)	
470	320	450	300	50	1450	0376.1	1 unit(s)	
575	375	550	350	50	1750	Y756.1	1 unit(s)	
670	430	650	400	50	2485	YH48.1	1 unit(s)	

Bowls and tubs

Dish ROTILABO®

ROTH SELECTION. Material: Chrome-nickel steel 18/10.
Autoclavable: Yes.
Rectangular with two drop handles on the sides. Suitable for foodstuffs.



A 121 °C

Not a medical device / Not an IVD product

Outer length (mm)	Outer width (mm)	Height (mm)	Art. No.	Pack Qty.	
365	270	75	8103.1	1 unit(s)	
470	325	90	8104.1	1 unit(s)	
520	325	90	8105.1	1 unit(s)	
630	450	90	8106.1	1 unit(s)	

Dish ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes.
Double-thickness walls 0.8 mm. Suitable for foodstuffs.



A 121 °C

Not a medical device / Not an IVD product

Outer length (mm)	Outer width (mm)	Height (mm)	Art. No.	Pack Qty.	
210	145	20	9492.1	1 unit(s)	
290	210	20	9493.1	1 unit(s)	
580	210	20	9494.1	1 unit(s)	
580	400	20	9495.1	1 unit(s)	

Work support ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes.
Thickness: 1.1 mm. Easy to clean and disinfect. Suitable for foodstuffs.



A 121 °C

Not a medical device / Not an IVD product

Outer length (mm)	Outer width (mm)	Height (mm)	Art. No.	Pack Qty.	
600	400	7	YH28.1	1 unit(s)	



Bowls and tubs

Dish ROTILABO® with perforated sheet

ROTH SELECTION. Material: Stainless steel 18/8. Autoclavable: Yes.
Removable grate. Perforation inside micrometer 19 mm, Hole Ø 12 mm.
Central hole Ø 19 mm.



A121 °C

Not a medical device / Not an IVD product

Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
310	125	30	X103.1	1 unit(s)	
320	200	27	X104.1	1 unit(s)	
440	270	27	X111.1	1 unit(s)	
560	200	27	X105.1	1 unit(s)	

Dish ROTILABO® with perforated sheet

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes.
With angled walls and inserted perforated sheet on 7 mm pegs. Suitable
for foodstuffs.



A121 °C

Not a medical device / Not an IVD product

Upper length (mm)	Upper width (mm)	Lower length (mm)	Lower width (mm)	Height (mm)	Inner height (mm)	Weight (g)	Art. No.	Pack Qty.	
500	350	435	285	60	50	1800	0377.1	1 unit(s)	

Bowls and tubs

Water bath dish ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10, matt polished.
Heavy-duty, with two side handles. Suitable for foodstuffs.



Delivery without lid, please order separately.

Not a medical device / Not an IVD product

Vol. (l)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.
0.5	155	105	55	YY42.1	1 unit(s)
1.0	155	105	75	YY43.1	1 unit(s)
2.5	155	105	160	YY44.1	1 unit(s)
5.0	155	155	235	YY45.1	1 unit(s)
13.0	240	240	235	YY46.1	1 unit(s)

Accessories lid for water bath trays

ROTH SELECTION. Material: Stainless steel 18/10, matt polished.

Not a medical device / Not an IVD product

Length (mm)	Width (mm)	Suitable for	Art. No.	Pack Qty.
155	105	Water bath trays, Art. No. YY42.1, YY43.1 and YY44.1	YY47.1	1 unit(s)
155	155	Water bath tray, Art. No. YY45.1	YY48.1	1 unit(s)
240	240	Water bath tray, Art. No. YY46.1	YY49.1	1 unit(s)

Sterilisation tray

Nalgene. Material: PP. Autoclavable: Yes.

The flat design of the sterilisation tray allows pipettes and instruments to be fully immersed in the disinfectant solution. The plastic construction reduces damage and breakage of glass pipettes. The tray is stackable, lightweight and quiet to use.

- Suitable for sterilising pipettes, instruments and small laboratory items up to a max. length of 40 cm
- Resistant to aggressive disinfectant solutions
- The lid that can be fitted into the tray groove prevents the spillage of hazardous liquids
- Lid with recessed grip



Delivery incl. lid.

Not a medical device / Not an IVD product

Outer length (mm)	Outer width (mm)	Height (mm)	Art. No.	Pack Qty.
460	150	67	1L1C.1	1 unit(s)

Алматы (7273)495-231	Иваново (4932)77-34-06	Магнитогорск (3519)55-03-13	Ростов-на-Дону (863)308-18-15	Тольятти (8482)63-91-07
Ангарск (3955)60-70-56	Ижевск (3412)26-03-58	Москва (495)268-04-70	Рязань (4912)46-61-64	Томск (3822)98-41-53
Архангельск (8182)63-90-72	Иркутск (395)279-98-46	Мурманск (8152)59-64-93	Самара (846)206-03-16	Тула (4872)74-02-29
Астрахань (8512)99-46-04	Казань (843)206-01-48	Набережные Челны (8552)20-53-41	Саранск (8342)22-96-24	Тюмень (3452)66-21-18
Барнаул (3852)73-04-60	Калининград (4012)72-03-81	Нижний Новгород (831)429-08-12	Санкт-Петербург (812)309-46-40	Ульяновск (8422)24-23-59
Белгород (4722)40-23-64	Калуга (4842)92-23-67	Новокузнецк (3843)20-46-81	Саратов (845)249-38-78	Улан-Удэ (3012)59-97-51
Благовещенск (4162)22-76-07	Кемерово (3842)65-04-62	Ноябрьск (3496)41-32-12	Севастополь (8692)22-31-93	Уфа (347)229-48-12
Брянск (4832)59-03-52	Киров (8332)68-02-04	Новосибирск (383)227-86-73	Симферополь (3652)67-13-56	Хабаровск (4212)92-98-04
Владивосток (423)249-28-31	Коломна (4966)23-41-49	Омск (3812)21-46-40	Ставрополь (8652)20-65-13	Чебоксары (8352)28-53-07
Владикавказ (8672)28-90-48	Кострома (4942)77-07-48	Орел (4862)44-53-42	Сочи (862)225-72-31	Челябинск (351)202-03-61
Владимир (4922)49-43-18	Краснодар (861)203-40-90	Оренбург (3532)37-68-04	Сыктывкар (8212)25-95-17	Череповец (8202)49-02-64
Волгоград (844)278-03-48	Красноярск (391)204-63-61	Пенза (8412)22-31-16	Тамбов (4752)50-40-97	Чита (3022)38-34-83
Вологда (8172)26-41-59	Курск (4712)77-13-04	Петрозаводск (8142)55-98-37	Сургут (3462)77-98-35	Якутск (4112)23-90-97
Воронеж (473)204-51-73	Курган (3522)50-90-47	Псков (8112)59-10-37	Тверь (4822)63-31-35	Ярославль (4852)69-52-93
Екатеринбург (343)384-55-89	Липецк (4742)52-20-81	Пермь (342)205-81-47		
Россия (495)268-04-70	Казахстан (772)734-952-31	Киргизия (996)312-96-26-47		