

Стойки, боксы, коробка для флаконов, пробирок CARL ROTH ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47



Combi racks

Test tube stands combined

ROTH SELECTION. Material: POM, rubber feet. Autoclavable: Yes.

Versatile rack with replaceable grid- inserts.

- Grid inserts in 5 colours for different test tube sizes
- 2 different test tube sizes are possible per rack
- Fast replacement of individual grids
- 3 levels for stability
- Stackable
- Temperature stability -80 to +121 °C



A121 °C



Delivery incl. rack and a set of 4 grids.

For glass Ø 10–13 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
84	14 x 6	318 x 111 x 67	orange	YA95.1	1 unit(s)	

For glass Ø 13–16 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
60	12 x 5	318 x 111 x 67	blue	YA96.1	1 unit(s)	

For glass Ø 16–20 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
40	10 x 4	318 x 111 x 67	green	YA97.1	1 unit(s)	

For glass Ø 20–25 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
24	8 x 3	318 x 111 x 67	pink	YA98.1	1 unit(s)	

For glass Ø 25–30 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
12	6 x 2	318 x 111 x 67	yellow	YA99.1	1 unit(s)	

Sample stands Pop-up™ racks

ROTH SELECTION. Material: PP. Autoclavable: No.

- Folds up for space-saving storage
- Height when folded up 21 mm
- Easy to unfold for quick use
- Numbered slots for identification of samples
- Robust design



For 15/50 ml centrifuge tubes

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
33	21 x Ø 17 / 12 x Ø 30 mm	255 x 137 x 72	blue	PC54.1	2 unit(s)	
33	21 x Ø 17 / 12 x Ø 30 mm	255 x 137 x 72	pink	PC55.1	2 unit(s)	

For 15 ml centrifuge tubes

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
45	45 x Ø 17 mm	255 x 137 x 72	violet	TP22.1	2 unit(s)	

For 50 ml centrifuge tubes

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
18	18 x Ø 30 mm	255 x 137 x 72	green	TP23.1	2 unit(s)	

Combi racks



Sample stands universal

ROTH SELECTION. Material: PP. Autoclavable: No.

Slots with adjustable size, can be opened by inserting the vessels or by turning a knob.

- Suitable for: 15/50 ml centrifuge tubes

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
12	194 x 140 x 51	NY31.1	1 unit(s)	



Sample stands for 50 ml centrifuge tubes

ROTH SELECTION.

For tubes with conical or round base, Ø 27.5 to 29 mm.

(1) Plastic

Material: PS.

No. of slots	Ø x H (mm)	Art. No.	Pack Qty.	
1	92 x 30	TP29.1	5 unit(s)	

(2) Metal

Material: Steel epoxy-coated. Autoclavable: Yes.

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
1	90 x 67 x 67	TP30.1	1 unit(s)	



A 121 °C



Sample stands ROTILABO® combined (4-in-1)

ROTH SELECTION. Material: PP. Autoclavable: Yes.

Space-saving use on a laboratory bench, beneath the vent, in a water bath and in a freezer. For centrifuge tubes, test tubes, reaction vials, PCR-vials. Temperature stability -90 to +121 °C.

With four sides which can be equipped as required:

- Hole Ø 30 mm – 4 x 50 ml centrifuge tubes
- Hole Ø 17 mm – 12 x 15 ml centrifuge tubes
- Hole Ø 11 mm – 32 x 1.5/2 ml reaction vials
- Hole Ø 8 mm – 32 x 0.5 ml reaction vials

L x W x H (mm)	Colour	Art. No.	Pack Qty.	
174 x 95 x 52	colourless	HYN6.1	5 unit(s)	
174 x 95 x 52	neon blue	K082.1	1 unit(s)	
174 x 95 x 52	neon green	K083.1	1 unit(s)	
174 x 95 x 52	neon orange	K084.1	1 unit(s)	
174 x 95 x 52	neon pink	K085.1	1 unit(s)	
174 x 95 x 52	neon yellow	K086.1	1 unit(s)	

Set

Assorted colours (neon blue, neon green, neon orange, neon pink, neon yellow).

L x W x H (mm)	Colour	Art. No.	Pack Qty.	
174 x 95 x 52	assorted colours	HYN7.1	5 unit(s)	



A 121 °C

Sample stands ROTILABO® combined (4-in-1) cubic

ROTH SELECTION. Material: PP. Autoclavable: Yes.

Can be connected in any way at the open sides. For centrifuge tubes, test tubes and reaction vials.

With four sides which can be equipped as required:

- Hole Ø 30 mm – 4 x 50 ml centrifuge tubes
- Hole Ø 17 mm – 10 x 15 ml centrifuge tubes
- Hole Ø 13 mm – 12 x test tubes Ø 12 x L 75/100 mm
- Hole Ø 12 mm – 16 x 1.5/2 ml reaction vials

Assorted colours (blue, yellow, green, orange, pink).

L x W x H (mm)	Colour	Art. No.	Pack Qty.	
108 x 108 x 108	assorted colours	AK71.1	5 unit(s)	



Combi racks



A 121 °C

Sample stands ROTILABO® combined (4-in-1)

ROTH SELECTION. Material: PP. Autoclavable: Yes.

For round, flat and conical blood collection and culture tubes from 3.5 to 10 ml, and for 15 ml centrifuge tubes.

Suitable for: Vessel Ø 13–16 mm.

Recommended tube size:

13 x 75/100 mm, 16 x 100/125 mm, 17 x 120 mm.

- Can be combined individually by fitting them together
- Alphanumerically coded
- Suitable for use in a water bath

With four sides which can be equipped as required:

- Side A: 12 x hole Ø 13 mm
- Side B: 18 x hole Ø 13 mm
- Side C: 12 x hole Ø 17 mm
- Side D: 18 x hole Ø 17 mm

L x W x H (mm)	Colour	Art. No.	Pack Qty.	
232 x 95 x 70	blue	KKY8.1	5 unit(s)	
232 x 95 x 70	violet	KKY9.1	5 unit(s)	
232 x 95 x 70	green	KLA0.1	5 unit(s)	
232 x 95 x 70	colourless	KLA1.1	5 unit(s)	

Set

Assorted colours (blue, purple, green, red neutral).

L x W x H (mm)	Colour	Art. No.	Pack Qty.	
232 x 95 x 70	assorted colours	KLA2.1	1 set	



A 121 °C

Sample stands Multi-tube rack

ROTH SELECTION. Material: POM. Autoclavable: Yes.

- Suitable for vessel Ø 12–30 mm
- Variable use thanks to 4 flexible slots
- Can be connected together

No. of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
4	181 x 56 x 76	blue/green	AXH9.1	2 unit(s)	
4	181 x 56 x 76	blue/yellow	AXK0.1	2 unit(s)	
4	181 x 56 x 76	blue	AXK1.1	2 unit(s)	
4	181 x 56 x 76	green	AXK2.1	2 unit(s)	
4	181 x 56 x 76	yellow	AXK3.1	2 unit(s)	
4	181 x 56 x 76	white	AXK4.1	2 unit(s)	
4	181 x 56 x 76	pink	AXK5.1	2 unit(s)	



A 121 °C

Sample stands ROTILABO® combined (4-in-1) rotatable

ROTH SELECTION. Material: PP. Autoclavable: Yes.

Four coloured modules are securely connected to each other but can each be turned individually so that reaction vials and tubes of different sizes can be stored and transported in a single rack.

For reaction vials

- Hole Ø 16.4 mm: 6 x 15 ml centrifuge tubes
- Hole Ø 11.2 mm: 9 x 1.5 to 2 ml reaction vials
- Hole Ø 8.0 mm: 12 x 0.5 ml reaction vials
- Hole Ø 6.5 mm: 32 x PCR single vials/4 x 8-tube PCR strips

L x W x H (mm)	Colour	Art. No.	Pack Qty.	
100 x 100 x 208	violet, blue, green, yellow	HC26.1	1 unit(s)	

For centrifuge tubes

- Hole Ø 29.6 mm: 2 x 50 ml centrifuge tubes
- Hole Ø 16.4 mm: 8 x 15 ml centrifuge tubes
- Hole Ø 20.5 mm: 6 x Ø 16 to 20 mm tubes
- Hole Ø 13.5 mm: 10 x Ø 10 to 13 mm tubes

L x W x H (mm)	Colour	Art. No.	Pack Qty.	
112 x 112 x 224	violet, blue, green, yellow	HC25.1	1 unit(s)	



A 121 °C

Reaction vial stands for 5 ml tubes

ROTH SELECTION. Material: PP. Autoclavable: Yes.

With side handles for transporting and stacking.

- Suitable for: 5 ml reaction vials/tube Ø 13–16 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
40	4 x 10	260 x 113 x 105	CXP6.1	1 unit(s)	

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Россия (495)268-04-70

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Казахстан (772)734-952-31

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Киргизия (996)312-96-26-47

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93