

Криобоксы, боксы для хранения CARL ROTH ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47



Cryo boxes, storage boxes

Cryogenic box ROTILABO® cardboard

ROTH SELECTION. Material: Cardboard (resistant to low temperatures).

Flexible storage of cryogenic vials and other sample containers through variable configurations using different dividers. Also suitable for larger containers without grid dividers.



Fitting table

Sample vessel type	Cryogenic box	Divider insert
0.5 ml reaction vials	Y054.2/CLE8.1	Y050.2
1.5/2 ml reaction vials	P895.3/CLE9.1	P900.2
Cryogenic vials up to 2 ml, internal thread	P895.3/CLE9.1	P900.2
Cryogenic vials up to 2 ml, external thread	P895.3/CLE9.1	P899.2
5 ml cryogenic vials	P897.2/CLH5.1	Y046.2

Delivery without divider inserts, please order separately.

With water-repellent standard coating

Version	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	136	136	50	white	P895.2	50 unit(s)	
Slip lid	136	136	50	white	P895.3	10 unit(s)	
Slip lid	136	136	75	white	P896.2	10 unit(s)	
Slip lid	136	136	100	white	P897.2	10 unit(s)	
Slip lid	136	136	50	blue	XP98.1	10 unit(s)	
Slip lid	136	136	50	red	XP99.1	10 unit(s)	
Slip lid	136	136	50	green	XT00.1	10 unit(s)	
Slip lid	136	136	50	yellow	XT01.1	10 unit(s)	
Slip lid	136	136	32	white	Y054.2	10 unit(s)	
Slip lid	136	136	130	white	Y055.2	10 unit(s)	

With waterproof plastic coating

Version	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	136	136	32	white	CLE8.1	10 unit(s)	
Slip lid	136	136	50	white	CLE9.1	10 unit(s)	
Slip lid	136	136	50	blue	CLH0.1	10 unit(s)	
Slip lid	136	136	50	red	CLH1.1	10 unit(s)	
Slip lid	136	136	50	green	CLH2.1	10 unit(s)	
Slip lid	136	136	50	yellow	CLH3.1	10 unit(s)	
Slip lid	136	136	75	white	CLH4.1	10 unit(s)	
Slip lid	136	136	100	white	CLH5.1	10 unit(s)	
Slip lid	136	136	130	white	CLH6.1	10 unit(s)	

Cryo boxes, storage boxes

Accessories divider inserts cardboard

ROTH SELECTION. Material: Cardboard.

For dividing the cryogenic boxes.

Compartment size (mm)	Slots	Height (mm)	Art. No.	Pack Qty.	
31 x 31	16 (4 x 4)	30	Y037.2	10 unit(s)	
31 x 31	16 (4 x 4)	65	Y038.2	10 unit(s)	
24 x 24	25 (5 x 5)	30	Y039.2	10 unit(s)	
24 x 24	25 (5 x 5)	65	Y040.2	10 unit(s)	
20 x 20	36 (6 x 6)	30	Y041.2	10 unit(s)	
20 x 20	36 (6 x 6)	65	Y042.2	10 unit(s)	
17 x 17	49 (7 x 7)	30	Y043.2	10 unit(s)	
17 x 17	49 (7 x 7)	65	Y044.2	10 unit(s)	
15 x 15	64 (8 x 8)	30	P898.2	10 unit(s)	
15 x 15	64 (8 x 8)	65	Y045.2	10 unit(s)	
13,5 x 13,5	81 (9 x 9)	30	P899.2	10 unit(s)	
13,5 x 13,5	81 (9 x 9)	65	Y046.2	10 unit(s)	
12,5 x 12,5	100 (10 x 10)	30	P900.2	10 unit(s)	
12,5 x 12,5	100 (10 x 10)	65	Y047.2	10 unit(s)	
10 x 10	144 (12 x 12)	25	Y048.2	10 unit(s)	
10 x 10	144 (12 x 12)	30	Y049.2	10 unit(s)	
9 x 9	169 (13 x 13)	25	Y050.2	10 unit(s)	
9 x 9	169 (13 x 13)	30	Y051.2	10 unit(s)	
8 x 8	196 (14 x 14)	25	Y052.2	10 unit(s)	
8 x 8	196 (14 x 14)	30	Y053.2	10 unit(s)	

Cryogenic box ROTILABO® cardboard maxi

ROTH SELECTION. Material: Cardboard (plastic-coated), waterproof and resistant to low temperatures.

For storing samples in reaction vials and centrifuge tubes. Also suitable for larger containers without grid dividers.

Temperature stability up to -80 °C



Delivery incl. divider inserts.

For 5 ml reaction vials

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	18 x 18	36 (6 x 6)	136	136	75	HNY0.1	10 unit(s)	

For 15 ml centrifuge tubes

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	18 x 18	36 (6 x 6)	148	148	130	CLE7.1	5 unit(s)	

For 50 ml centrifuge tubes

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	33 x 33	16 (4 x 4)	148	148	130	CLE6.1	5 unit(s)	



Cryo boxes, storage boxes

Storage box for centrifuge tubes

ROTH SELECTION. Material: PP. Autoclavable: Yes.

Replaceable inserts that can be inserted quickly offering ample depth for stable, upright storage of vials. Translucent lid and numbered inserts for fast identification of samples.

- Suitable for 15/50 ml centrifuge tubes
- Temperature stability -80 to +121 °C




Version	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	132	132	129	PC60.1	5 unit(s)	

Storage box ROTILABO®

ROTH SELECTION. Material: PP. Autoclavable: Yes.

With transparent lid, labelling area on two sides.

- Height without lid: 102 mm
- Temperature stable from -90 to +121 °C



(1) For 15 ml centrifuge tubes

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	18 x 18	36 (6 x 6)	144	144	127	green	HX01.1	2 unit(s)	

(2) For 50 ml centrifuge tubes

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	30 x 30	16 (4 x 4)	144	144	127	blue	HX02.1	2 unit(s)	

Storage box Work2Store™

ROTH SELECTION. Autoclavable: Yes.

Variable all-in-one rack and storage boxes: Can be pulled apart for easy sample handling with large distances between rows. Can be stored flat to save space.

Length when open: 310 mm.

- Alphanumeric numbering
- Stackable
- Temperature stable from -80 to +121 °C.



For reaction vials and cryogenic vials

- Slots: 64 x Ø 12 mm
- Holes with latching devices for cryogenic vials

Version	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	133	133	53	NL88.1	1 unit(s)	

Cryo boxes, storage boxes

Cryogenic box ROTILABO® without division

ROTH SELECTION. Material: PP. Autoclavable: Yes.

Flexible storage of cryogenic vials and other sample containers through variable configurations using different dividers. Also suitable for larger containers without grid dividers.

- Outlets in the base
- Temperature stability -90 to +121 °C



A 121 °C



Delivery without grid inserts, please order separately.

Version	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Hinged lid	133	133	52	natural	CLH7.1	5 unit(s)	
Hinged lid	133	133	52	yellow	CLH8.1	5 unit(s)	
Hinged lid	133	133	52	red	CLH9.1	5 unit(s)	
Hinged lid	133	133	52	green	CLK0.1	5 unit(s)	
Hinged lid	133	133	52	blue	CLK1.1	5 unit(s)	
Hinged lid	133	133	52	black	CLK2.1	5 unit(s)	
Hinged lid	133	133	75	natural	CLK3.1	5 unit(s)	
Hinged lid	133	133	75	yellow	CLK4.1	5 unit(s)	
Hinged lid	133	133	75	red	CLK5.1	5 unit(s)	
Hinged lid	133	133	75	green	CLK6.1	5 unit(s)	
Hinged lid	133	133	75	blue	CLK7.1	5 unit(s)	
Hinged lid	133	133	75	black	CLK8.1	5 unit(s)	

Accessories grid insert for cryogenic boxes

ROTH SELECTION. Material: PP. Autoclavable: Yes.

A 121 °C



Compartment size (mm)	Slots	Height (mm)	Art. No.	Pack Qty.	
11,2 x 11,2	100 (10 x 10)	30	CLL3.1	5 unit(s)	
12,8 x 12,8	81 (9 x 9)	30	CLL2.1	5 unit(s)	
15,8 x 15,8	64 (8 x 8)	30	CLL1.1	5 unit(s)	
17,3 x 17,3	49 (7 x 7)	30	CLL0.1	5 unit(s)	
20,5 x 20,5	36 (6 x 6)	30	CLK9.1	5 unit(s)	

Cryogenic box foldable

ROTH SELECTION. Material: PP.

- Suitable for 1.5/2 ml reaction vials, 2 ml cryogenic vials
- Thin-walled
- Moisture-resistant
- Temperature stability -196 to +50 °C
- Suitable for storing samples in the gas phase of liquid nitrogen



Delivery incl. divider inserts (pre-installed), boxes simply need to be assembled.

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	12 x 12	81 (9 x 9)	124	124	49	white	HYN2.1	10 unit(s)	
Slip lid	12 x 12	81 (9 x 9)	124	124	49	blue	HYN3.1	10 unit(s)	
Slip lid	12 x 12	81 (9 x 9)	124	124	49	green	HYN4.1	10 unit(s)	
Slip lid	12 x 12	81 (9 x 9)	124	124	49	violet	HYN5.1	10 unit(s)	



Cryo boxes, storage boxes

Cryogenic box foldable

ROTH SELECTION. Material: PP corrugated, moisture-resistant.

- Temperature stability -80 to +50 °C



Delivery incl. divider inserts.

For 0.2 ml PCR reaction vials

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	7,9 x 7,9	144 (12 x 12)	133	130	30	ENL0.1	10 unit(s)	

For 0.5 ml reaction vials

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	12 x 12	81 (9 x 9)	133	130	42	ENL1.1	10 unit(s)	

For 1.5/2 ml reaction vials

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	12 x 12	81 (9 x 9)	133	130	51	AXE4.1	10 unit(s)	

For 5 ml reaction vials

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	22,6 x 22,6	25 (5 x 5)	133	130	75	ENL2.1	10 unit(s)	

For 15 ml centrifuge tubes

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	19,8 x 19,8	36 (6 x 6)	147	145	120	ENL3.1	10 unit(s)	

For 50 ml centrifuge tubes


Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	31,6 x 31,6	16 (4 x 4)	147	145	120	ENL4.1	10 unit(s)	

Storage box with inserts

ROTH SELECTION. Material: PP.

Box with 8 inserts. Individually removable inserts for flexible working, alphanumeric marking for sample identification.

- Temperature stability -80 to +70 °C
- Autoclavable inserts



For 1.5/2 ml vessels

- Suitable for 1.5/2 ml reaction vials/2 ml cryogenic vials
- Latching device at bottom for cryogenic vials enables single-handed opening and closing

Version	Compartment size (mm)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	Ø 12,3	129	130	35	ENL8.1	1 unit(s)	

Cryo boxes, storage boxes

Cryogenic box without division

ROTH SELECTION. Material: PP, transparent. Autoclavable: Yes.
For storage of reaction vials, cryogenic vials and other sample containers.

- Height-adjustable lid, locks at 43 mm, 50 mm and 55 mm
- Temperature stability -196 to +121 °C



Delivery without divider inserts, please order separately.

Version	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	130	130	43	ACL8.1	5 unit(s)	

Divider inserts

ROTH SELECTION. Material: PP. Autoclavable: Yes.

Individual dividers for flexible division of ACL8.1 cryogenic boxes.

Variable division by using between 2 and 16 divider inserts:

2 divider inserts: 4 slots (50 x 55 to 65 x 65 mm)

4 divider inserts: 9 slots (40 x 40 mm)

8 divider inserts: 16 slots (25 x 25 mm)

16 divider inserts: 81 slots (12 x 12 mm)




L x W x H (mm)	Art. No.	Pack Qty.	
121 x 1 x 30	ACL9.1	16 unit(s)	

Cryogenic box ROTILABO®

ROTH SELECTION. Material: PC. Autoclavable: Yes.

- Suitable for 1–2 ml cryogenic vials
- Lid with alphanumeric labelling
- Stackable
- With outlets in the base
- Temperature stability -196 to +121 °C



Slip lid slots, round

Pic.	Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
(1)	Slip lid	Ø 12,5	25 (5 x 5)	76	76	53	red	HE00.1	1 unit(s)	

Hinged lid slots, round

Assorted colours: Blue, red, green, purple.

Delivery incl. removal aid.

Pic.	Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
(2)	Hinged lid	Ø 12,5	81 (9 x 9)	133	133	53	assorted colours	HE01.1	4 unit(s)	

Hinged lid pin pattern

Pic.	Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
(3)	Hinged lid	12,5 x 12,5	100 (10 x 10)	133	133	53	blue	HE03.1	1 unit(s)	



Cryo boxes, storage boxes

Cryogenic box ROTILABO® for 5 ml cryogenic vials

ROTH SELECTION. Material: PC, transparent. Autoclavable: Yes.

- Suitable for 5 ml cryogenic vials
- With embossed numeric labelling
- Removable insert
- Stackable
- Temperature stability from -196 °C




Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Hinged lid	13 x 13	81 (9 x 9)	145	142	97	ACL7.1	1 unit(s)	

Cryogenic box ROTILABO®

ROTH SELECTION. Material: PC. Autoclavable: Yes.

With transparent lid and removable coloured insert.

- Slots and lid are labelled numerically
- Labelling area on box and lid
- Outlets in the base
- Stackable
- Temperature stability -196 to +121 °C



For 1–2 ml cryogenic vials

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	12,2 x 12,2	81 (9 x 9)	133	133	52	blue	K585.1	1 unit(s)	
Slip lid	12,2 x 12,2	81 (9 x 9)	133	133	52	green	K586.1	1 unit(s)	
Slip lid	12,2 x 12,2	81 (9 x 9)	133	133	52	red	K587.1	1 unit(s)	
Slip lid	12,2 x 12,2	81 (9 x 9)	133	133	52	yellow	K588.1	1 unit(s)	

For 3–5 ml cryogenic vials

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	12,2 x 12,2	81 (9 x 9)	133	133	95	blue	K589.1	1 unit(s)	
Slip lid	12,2 x 12,2	81 (9 x 9)	133	133	95	green	K590.1	1 unit(s)	
Slip lid	12,2 x 12,2	81 (9 x 9)	133	133	95	red	K591.1	1 unit(s)	
Slip lid	12,2 x 12,2	81 (9 x 9)	133	133	95	yellow	K592.1	1 unit(s)	

Cryogenic box ROTILABO® 100 slots

ROTH SELECTION. Material: PC. Autoclavable: Yes.

- Suitable for 1/2 ml cryogenic vials, internal thread
- Slots: Pin pattern
- Labelled numerically (lid and base)
- Outlets in the base
- Labelling area on the side
- Stackable
- Temperature stability -196 to +121 °C



Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	12,5 x 12,5	100 (10 x 10)	133	133	52	blue	K593.1	1 unit(s)	
Slip lid	12,5 x 12,5	100 (10 x 10)	133	133	52	green	K594.1	1 unit(s)	
Slip lid	12,5 x 12,5	100 (10 x 10)	133	133	52	red	K595.1	1 unit(s)	
Slip lid	12,5 x 12,5	100 (10 x 10)	133	133	52	yellow	K596.1	1 unit(s)	



Cryo boxes, storage boxes

Storage box 81 slots

ROTH SELECTION. Material: PP. Autoclavable: Yes.

- Suitable for 1.5/2 ml reaction vials
- Slots are labelled alphanumerically
- Outlets in the base
- Black box with black lid for storing light-sensitive samples
- Temperature stability -90 to +121 °C



Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	13 x 13	81 (9 x 9)	130	130	47	colourless	NL80.1	5 unit(s)	
Slip lid	13 x 13	81 (9 x 9)	130	130	47	blue	NL81.1	5 unit(s)	
Slip lid	13 x 13	81 (9 x 9)	130	130	47	green	NL82.1	5 unit(s)	
Slip lid	13 x 13	81 (9 x 9)	130	130	47	pink	NL83.1	5 unit(s)	
Slip lid	13 x 13	81 (9 x 9)	130	130	47	yellow	NL84.1	5 unit(s)	
Slip lid	13 x 13	81 (9 x 9)	130	130	47	orange	NL85.1	5 unit(s)	
Slip lid	13 x 13	81 (9 x 9)	130	130	47	black	NL86.1	5 unit(s)	

Set

Assorted colours (blue, green, pink, yellow, orange).

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	13 x 13	81 (9 x 9)	130	130	47	assorted colours	NL87.1	5 unit(s)	

Storage box 50 slots

ROTH SELECTION. Material: PP. Autoclavable: Yes.

- Suitable for 1.5/2 ml reaction vials, 2 ml cryogenic vials



(1) Hinged lid set colourless

Slots are numbered alphanumerically.

Version	Compartment size (mm)	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Hinged lid	13 x 13	142	92	56	colourless	NY15.1	5 unit(s)	

(1) Hinged lid set assorted colours

Slots are numbered alphanumerically.

Assorted colours: Blue, green, purple, yellow, orange.

Version	Compartment size (mm)	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Hinged lid	13 x 13	142	92	56	assorted colours	NY16.1	5 unit(s)	

(2) Slip lid

Tube height variable up to 55 mm due to height-adjustable plug-in lid.

Version	Compartment size (mm)	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	12,5 x 12,5	144	75	35	colourless	N201.1	1 unit(s)	
Slip lid	12,5 x 12,5	144	75	35	pink	P184.1	1 unit(s)	
Slip lid	12,5 x 12,5	144	75	35	yellow	P185.1	1 unit(s)	
Slip lid	12,5 x 12,5	144	75	35	blue	P186.1	1 unit(s)	
Slip lid	12,5 x 12,5	144	75	35	green	P187.1	1 unit(s)	

Cryo boxes, storage boxes

Storage box ROTILABO® 200 slots

ROTH SELECTION. Material: PP. Autoclavable: Yes.

With two tabs on the front for inserting tags.

- Suitable for 0.5–2 ml reaction vials
- With alphanumeric labelling
- Stackable
- Temperature stability -80 to +135 °C



Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Hinged lid	12 x 12	200 (10 x 20)	275	151	64	black	ACN1.1	1 unit(s)	
Hinged lid	12 x 12	200 (10 x 20)	275	151	64	colourless	P786.1	1 unit(s)	

Storage box ROTILABO® slimline shape

ROTH SELECTION. Material: PP. Autoclavable: Yes.

Temperature stable from -80 to +135 °C.

- Suitable for 0.5–2 ml reaction vials
- Slip lid with alphanumeric labelling, loosely attached, adjustable height up to 52 mm
- Stackable



Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	12 x 12	100 (5 x 20)	287	76	38	black	ACN0.1	1 unit(s)	

Cryogenic rack drawer rack

ROTH SELECTION. Material: Stainless steel, polished, drawers made of aluminium.

For the storage of cryogenic boxes in deep freezers, less effort is required for removal and insertion. Cabinet modules without pull-out stop. Suitable for cryogenic boxes in sizes from 130 x 130 to 136 x 136 mm.



Suitable for	Width (mm)	Height (mm)	Depth (mm)	Arrangement of compartments	Art. No.	Pack Qty.	
Box height: 50 mm	144	224	570	4 x 4	CLL9.1	1 unit(s)	
Box height: 50 mm	144	280	570	5 x 4	CLN0.1	1 unit(s)	
Box height: 50 mm	144	336	570	6 x 4	CLN1.1	1 unit(s)	
Box height: 50 mm	144	280	709	5 x 5	HHX5.1	1 unit(s)	
Box height: 75 mm	144	250	570	3 x 4	CLN2.1	1 unit(s)	
Box height: 75 mm	144	335	570	4 x 4	CLN3.1	1 unit(s)	



Cryo racks

Cryogenic rack stainless steel

ROTH SELECTION. Material: Stainless steel, polished.

Extremely stable and easy to clean. For storing cryogenic boxes in deep freezers. With carry handles. Suitable for cryogenic boxes in sizes from 130 x 130 to 136 x 136 mm.



Wide (for upright freezers)

Pic.	Suitable for	Width (mm)	Height (mm)	Depth (mm)	Arrangement of compartments	Art. No.	Pack Qty.	
(1)	Box height: 50 mm	140	274	424	3 x 5	LH44.1	1 unit(s)	
(1)	Box height: 50 mm	140	219	566	4 x 4	LH43.1	1 unit(s)	
(1)	Box height: 50 mm	140	274	566	4 x 5	LH42.1	1 unit(s)	
(1)	Box height: 50 mm	140	328	566	4 x 6	LH41.1	1 unit(s)	


Tall (for chest freezers)

Pic.	Suitable for	Width (mm)	Height (mm)	Depth (mm)	Arrangement of compartments	Art. No.	Pack Qty.	
(2)	Box height: 50 mm	140	495	140	1 x 9	LH49.1	1 unit(s)	
(2)	Box height: 50 mm	140	550	140	1 x 10	T440.1	1 unit(s)	
(2)	Box height: 50 mm	140	604	140	1 x 11	LH47.1	1 unit(s)	
(2)	Box height: 50 mm	140	659	140	1 x 12	LH48.1	1 unit(s)	
(2)	Box height: 50 mm	140	713	140	1 x 13	LH46.1	1 unit(s)	
(2)	Box height: 75–77 mm	140	641	140	1 x 8	LH50.1	1 unit(s)	
(2)	Box height: 75–77 mm	140	720	140	1 x 9	LH51.1	1 unit(s)	

Cryogenic rack Aluminium

ROTH SELECTION. Material: Aluminium.

With carry handle. At least 30 % lighter than stainless steel racks. For the storage of cryogenic boxes in chest freezers or upright freezers, less effort is required for removal and insertion. Improved design for faster cooling rate. Suitable for cryogenic boxes in sizes from 130 x 130 to 136 x 136 mm.



Wide (for upright freezers)

Pic.	Suitable for	Width (mm)	Height (mm)	Depth (mm)	Arrangement of compartments	Art. No.	Pack Qty.	
(1)	Box height: 50 mm	140	274	566	4 x 5	CLL6.1	1 unit(s)	
(1)	Box height: 50 mm	140	328	566	4 x 6	CLL7.1	1 unit(s)	

Tall (for chest freezers)

Pic.	Suitable for	Width (mm)	Height (mm)	Depth (mm)	Arrangement of compartments	Art. No.	Pack Qty.	
(2)	Box height: 50 mm	140	550	140	1 x 10	CLL8.1	1 unit(s)	

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47