

Бутили-капельницы с пипетками CARL ROTH ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47

Dropper bottles



Dropper bottle with pipette

DWK Life Sciences. Material: Soda-lime glass.

Complete with suction cups and replaceable ground glass joint pipette made of clear glass.

Clear glass

Vol. (ml)	Height (mm)	Standard ground glass joint	Art. No.	Pack Qty.	
50	79	14/15	T930.1	10 unit(s)	
100	105	14/15	T931.1	10 unit(s)	

Brown glass

Vol. (ml)	Height (mm)	Standard ground glass joint	Art. No.	Pack Qty.	
50	79	14/15	T932.1	10 unit(s)	
100	105	14/15	T933.1	10 unit(s)	



Dropper bottle with pipette

Marlenfeld. Material: Glass.

Complete with clear-glass pipette, suction cups and ground glass stopper made of PP.

Clear glass

Vol. (ml)	Height (mm)	Standard ground glass joint	Art. No.	Pack Qty.	
50	77	14/15	TT98.1	10 unit(s)	
100	103	14/23	TT99.1	6 unit(s)	
250	131	19/26	TX00.1	10 unit(s)	

Brown glass

Vol. (ml)	Height (mm)	Standard ground glass joint	Art. No.	Pack Qty.	
50	77	14/15	TX01.1	10 unit(s)	
100	103	14/23	TX02.1	6 unit(s)	
250	131	19/26	TX03.1	10 unit(s)	



Dropper bottle with glass stopper

DWK Life Sciences. Material: Soda-lime glass.

Turning the plug seals the bottle tight.

Clear glass

Vol. (ml)	Height (mm)	Ø ext. (mm)	Art. No.	Pack Qty.	
50	95	45	ETL2.1	1 unit(s)	
100	120	56	ETL3.1	1 unit(s)	

Brown glass

Vol. (ml)	Height (mm)	Ø ext. (mm)	Art. No.	Pack Qty.	
50	95	45	ETL4.1	1 unit(s)	
100	120	56	ETL5.1	1 unit(s)	



Dropper bottle with pipette

Kartell. Material: PE, translucent.

Delivery incl. red screw cap, Dropper pipette and suction cup.

Vol. (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
50	35	96	TT69.1	10 unit(s)	
100	43	115	TT70.1	10 unit(s)	
250	60	140	TT71.1	10 unit(s)	



Dropper bottles



Dropper bottle

Nalgene. Material: LDPE, cap made of PP, dropper tip made of LDPE. Good chemical resistance. Ideal for small culture volumes. Precise single drop release. Droplet size approx. 40 µl (Aqua dest). **Leak-proof.**

Type 2752, transparent bottles, colour-coded lids

Lid colours: white, red, green, blue, yellow.

Vol. (ml)	Height (mm)	Art. No.	Pack Qty.	
4	41	PP10.1	25 unit(s)	
8	52	PP11.1	25 unit(s)	
15	62	PP12.1	25 unit(s)	

Type 2751, white bottles and lids

Suitable for UV-sensitive samples.

Vol. (ml)	Height (mm)	Art. No.	Pack Qty.	
4	41	HX22.1	25 unit(s)	
8	52	HX23.1	25 unit(s)	
15	62	HX24.1	25 unit(s)	



A 121 °C

Dropper bottle

Nalgene. Material: PP copolymer, dropper seal and cap made of PP. Autoclavable: Yes.

Model DS2420. For dispensing aqueous media. Also autoclavable when full with tip open. Tip sealing cap is combined with the screw cap and therefore cannot be lost.

Vol. (ml)	Art. No.	Pack Qty.	
125	ATL3.1	6 unit(s)	
250	ATL4.1	6 unit(s)	



Dropper bottle

Kartell. Material: LDPE, transparent.

Delivery incl. dropper lids.

Vol. (ml)	Thread	Art. No.	Pack Qty.	
50	GL 18	6541.1	6 unit(s)	
100	GL 18	6542.1	6 unit(s)	
250	GL 25	6544.1	4 unit(s)	
500	GL 25	6545.1	4 unit(s)	
1,000	GL 32	6548.1	4 unit(s)	



A 121 °C

Dropper bottle ROTILABO® fluoroplastic

ROTH SELECTION. Material: PTFE. Autoclavable: Yes. Chemically inert. Temperature stability -200 to +260 °C.

Delivery incl. caps.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Ø Inner neck (mm)	Art. No.	Pack Qty.	
25	33	80	25	LT80.1	1 unit(s)	
50	43	100	25	LT81.1	1 unit(s)	

Dropper bottles



Dosing bottle

Kartell. Material: PE, measuring cup made of PMP.

Adjustable measuring cup, graduated from 5 to 50 ml. Accuracy approx. 2%. Press on the bottle to fill the measuring cup. Surplus amounts flow back into the bottle automatically up to the measuring marker. Sealable filling tube.

Vol. (ml)	Art. No.	Pack Qty.	
250	6996.1	1 unit(s)	
500	7009.1	1 unit(s)	
1,000	7024.1	1 unit(s)	



Dropper bottle square

Bürkle. Material: HDPE, dropper closure made of LDPE, caps made of PVC.

With labelling space.

Vol. (ml)	Length (mm)	Width (mm)	Height (mm)	Ø Inner neck (mm)	Thread	Art. No.	Pack Qty.	
100	41	34	113	13.0	18 mm	KA95.1	10 unit(s)	
250	51	44	160	17.5	25 mm	KA96.1	10 unit(s)	
500	68	53	182	17.5	25 mm	KA97.1	10 unit(s)	
1,000	90	64	230	17.5	25 mm	KA98.1	10 unit(s)	



Spray bottle with atomiser

Rixius. Material: Glass brown, sprayer with cover cap and riser made of PE.

Spray volume approx. 0.05 ml per pump operation.

Vol. (ml)	Ø (mm)	Height (mm)	Thread	Art. No.	Pack Qty.	
20	29	98	GL 18	N145.1	3 unit(s)	



Spray bottle Turn'n'Spray

Bürkle. Material: LDPE, part in contact with media PP, PE, silicone and stainless steel.

Allows uninterrupted spraying from all positions, including overhead. Nozzle is continuously adjustable, nozzle Ø 0.6 mm. Spray quantity per stroke: 1.2 ml ±0.1 ml.

Vol. (ml)	Height (mm)	Ø Inner neck (mm)	Art. No.	Pack Qty.	
250	220	18	PX90.1	1 unit(s)	
500	240	18	PX91.1	1 unit(s)	
1,000	295	21	PX92.1	1 unit(s)	



Spray bottle

Bürkle. Material: PE, pump atomiser made of PP, transparent.

With protective cap, prevents unintentional pressing of the pump and serves as a dust protection cap.

Vol. (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
20	28	85	1C8X.1	5 unit(s)	
50	35	110	EP66.1	5 unit(s)	
100	37	161	EP67.1	2 unit(s)	
250	50	202	EP68.1	2 unit(s)	



Spray bottle with pump atomiser

Vitlab. Material: PP.

With pistol pump level and adjustable spray jet.

Vol. (ml)	Art. No.	Pack Qty.	
400	0499.1	1 unit(s)	
850	0500.1	1 unit(s)	

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47