

# Воронки лабораторные CARL ROTH ROTILABO

## Технические характеристики

Алматы (7273)495-231  
Ангарск (3955)60-70-56  
Архангельск (8182)63-90-72  
Астрахань (8512)99-46-04  
Барнаул (3852)73-04-60  
Белгород (4722)40-23-64  
Благовещенск (4162)22-76-07  
Брянск (4832)59-03-52  
Владивосток (423)249-28-31  
Владикавказ (8672)28-90-48  
Владимир (4922)49-43-18  
Волгоград (844)278-03-48  
Вологда (8172)26-41-59  
Воронеж (473)204-51-73  
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06  
Ижевск (3412)26-03-58  
Иркутск (395)279-98-46  
Казань (843)206-01-48  
Калининград (4012)72-03-81  
Калуга (4842)92-23-67  
Кемерово (3842)65-04-62  
Киров (8332)68-02-04  
Коломна (4966)23-41-49  
Кострома (4942)77-07-48  
Краснодар (861)203-40-90  
Красноярск (391)204-63-61  
Курск (4712)77-13-04  
Курган (3522)50-90-47  
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13  
Москва (495)268-04-70  
Мурманск (8152)59-64-93  
Набережные Челны (8552)20-53-41  
Нижний Новгород (831)429-08-12  
Новокузнецк (3843)20-46-81  
Ноябрьск (3496)41-32-12  
Новосибирск (383)227-86-73  
Омск (3812)21-46-40  
Орел (4862)44-53-42  
Оренбург (3532)37-68-04  
Пенза (8412)22-31-16  
Петрозаводск (8142)55-98-37  
Псков (8112)59-10-37  
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15  
Рязань (4912)46-61-64  
Самара (846)206-03-16  
Саранск (8342)22-96-24  
Санкт-Петербург (812)309-46-40  
Саратов (845)249-38-78  
Севастополь (8692)22-31-93  
Симферополь (3652)67-13-56  
Смоленск (4812)29-41-54  
Сочи (862)225-72-31  
Ставрополь (8652)20-65-13  
Сыктывкар (8212)25-95-17  
Тамбов (4752)50-40-97  
Сургут (3462)77-98-35  
Тверь (4822)63-31-35

Тольятти (8482)63-91-07  
Томск (3822)98-41-53  
Тула (4872)74-02-29  
Тюмень (3452)66-21-18  
Ульяновск (8422)24-23-59  
Улан-Удэ (3012)59-97-51  
Уфа (347)229-48-12  
Хабаровск (4212)92-98-04  
Чебоксары (8352)28-53-07  
Челябинск (351)202-03-61  
Череповец (8202)49-02-64  
Чита (3022)38-34-83  
Якутск (4112)23-90-97  
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47

# Funnels

## Funnels ROTILABO®

ROTH SELECTION. Material: Borosilicate glass 3.3. Autoclavable: Yes.  
With short stem. 60° angle.



| Ø Inner upper (mm) | Ø Outer outlet (mm) | Tube length (mm) | Suitable for              | Art. No. | Pack Qty.  |  |
|--------------------|---------------------|------------------|---------------------------|----------|------------|--|
| 35                 | 6                   | 35               | paper filter Ø 45-55 mm   | ACC3.1   | 10 unit(s) |  |
| 45                 | 6                   | 45               | paper filter Ø 55-70 mm   | ACC4.1   | 10 unit(s) |  |
| 55                 | 8                   | 55               | paper filter Ø 70-90 mm   | ACC5.1   | 10 unit(s) |  |
| 75                 | 8                   | 75               | paper filter Ø 125-150 mm | ACC6.1   | 10 unit(s) |  |
| 85                 | 9                   | 85               | paper filter Ø 150 mm     | ACC7.1   | 10 unit(s) |  |
| 100                | 10                  | 100              | paper filter Ø 150-185 mm | ACC8.1   | 10 unit(s) |  |
| 125                | 16                  | 125              | paper filter Ø 185 mm     | ACC9.1   | 1 unit(s)  |  |
| 150                | 16                  | 150              | paper filter Ø 240-270 mm | ACE0.1   | 1 unit(s)  |  |
| 200                | 24                  | 150              | paper filter Ø 320-385 mm | ACE1.1   | 1 unit(s)  |  |
| 250                | 30                  | 175              | paper filter Ø 400-450 mm | ACE2.1   | 1 unit(s)  |  |
| 300                | 34                  | 175              | paper filter Ø 500 mm     | ACE3.1   | 1 unit(s)  |  |

## Funnels

DWK Life Sciences. Material: DURAN® glass. Autoclavable: Yes.  
With short stem.



| Ø Outer upper (mm) | Ø Outer outlet (mm) | Tube length (mm) | Suitable for              | Art. No. | Pack Qty.  |  |
|--------------------|---------------------|------------------|---------------------------|----------|------------|--|
| 35                 | 6                   | 35               | Paper filter-Ø 45-55 mm   | C153.1   | 10 unit(s) |  |
| 45                 | 6                   | 45               | Paper filter-Ø 55-70 mm   | C154.1   | 10 unit(s) |  |
| 55                 | 8                   | 55               | Paper filter-Ø 70-90 mm   | C155.1   | 10 unit(s) |  |
| 70                 | 8                   | 70               | Paper filter-Ø 110-125 mm | C156.1   | 10 unit(s) |  |
| 80                 | 10                  | 80               | Paper filter-Ø 125-150 mm | C157.1   | 10 unit(s) |  |
| 100                | 10                  | 100              | Paper filter-Ø 150-185 mm | C158.1   | 1 unit(s)  |  |
| 150                | 16                  | 150              | Paper filter-Ø 240-270 mm | C159.1   | 1 unit(s)  |  |
| 200                | 26                  | 175              | Paper filter-Ø 320-385 mm | C160.1   | 1 unit(s)  |  |
| 300                | 30                  | 175              | Paper filter-Ø 500 mm     | C161.1   | 1 unit(s)  |  |



# Funnels



## Funnels with ribs

Bohemia. Material: Borosilicate glass 3.3. Autoclavable: Yes.  
For quick filtration.

| Ø Outer upper (mm) | Ø Outer outlet (mm) | Tube length (mm) | Art. No.      | Pack Qty. |  |
|--------------------|---------------------|------------------|---------------|-----------|--|
| 105                | 18                  | 80               | <b>NY75.1</b> | 1 unit(s) |  |
| 138                | 22                  | 90               | <b>NY76.1</b> | 1 unit(s) |  |
| 170                | 26                  | 120              | <b>NY77.1</b> | 1 unit(s) |  |
| 210                | 38                  | 140              | <b>NY78.1</b> | 1 unit(s) |  |
| 250                | 40                  | 155              | <b>NY79.1</b> | 1 unit(s) |  |



## Funnels ROTILABO®

ROTH SELECTION. Material: Soda-lime glass.  
With short stem. 60° angle.

| Ø Outer upper (mm) | Ø Outer outlet (mm) | Tube length (mm) | Art. No.      | Pack Qty.  |  |
|--------------------|---------------------|------------------|---------------|------------|--|
| 30                 | 6                   | 30               | <b>EE58.1</b> | 10 unit(s) |  |
| 35                 | 6                   | 35               | <b>EE59.1</b> | 10 unit(s) |  |
| 40                 | 6                   | 40               | <b>EE60.1</b> | 10 unit(s) |  |
| 45                 | 7                   | 45               | <b>EE61.1</b> | 10 unit(s) |  |
| 50                 | 7                   | 50               | <b>EE62.1</b> | 10 unit(s) |  |
| 55                 | 8                   | 55               | <b>EE63.1</b> | 10 unit(s) |  |
| 60                 | 8                   | 60               | <b>EE64.1</b> | 10 unit(s) |  |
| 70                 | 9                   | 70               | <b>EE65.1</b> | 10 unit(s) |  |
| 80                 | 9                   | 80               | <b>EE66.1</b> | 10 unit(s) |  |
| 100                | 10                  | 100              | <b>EE67.1</b> | 10 unit(s) |  |



## Funnels

Vitlab. Material: PP, transparent.  
With handle and lifting eye. Suitable for foodstuffs.

| Ø Upper (mm) | Ø Outlet (mm) | Height (mm) | Tube length (mm) | Art. No.      | Pack Qty. |  |
|--------------|---------------|-------------|------------------|---------------|-----------|--|
| 30           | 1.5           | 45          | 25               | <b>YA44.1</b> | 1 unit(s) |  |
| 30           | 4             | 47          | 25               | <b>YA45.1</b> | 1 unit(s) |  |
| 40           | 4             | 65          | 35               | <b>YA46.1</b> | 1 unit(s) |  |
| 50           | 7             | 85          | 43               | <b>YA47.1</b> | 1 unit(s) |  |
| 75           | 7             | 108         | 55               | <b>YA48.1</b> | 1 unit(s) |  |
| 100          | 8             | 155         | 77               | <b>YA49.1</b> | 1 unit(s) |  |
| 120          | 11            | 180         | 90               | <b>YA50.1</b> | 1 unit(s) |  |
| 150          | 15            | 220         | 95               | <b>YA51.1</b> | 1 unit(s) |  |



## Funnels ROTILABO®

ROTH SELECTION. Material: HDPE.  
Natural colours with smooth surface. With lifting eye and 4 ventilation ribs (except Art. No. 0583.1).  
Art. No. A828.1 incl. sieve Ø 50 mm, mesh size 0.4 mm.

| Ø Upper (mm) | Ø Outlet (mm) | Height (mm) | Art. No.      | Pack Qty. |  |
|--------------|---------------|-------------|---------------|-----------|--|
| 50           | 3             | 60          | <b>0583.1</b> | 1 unit(s) |  |
| 60           | 8             | 80          | <b>0584.1</b> | 1 unit(s) |  |
| 80           | 11            | 105         | <b>0585.1</b> | 1 unit(s) |  |
| 100          | 12            | 125         | <b>0586.1</b> | 1 unit(s) |  |
| 120          | 12            | 145         | <b>0587.1</b> | 1 unit(s) |  |
| 140          | 15            | 165         | <b>0588.1</b> | 1 unit(s) |  |
| 164          | 17            | 195         | <b>0589.1</b> | 1 unit(s) |  |
| 208          | 20            | 220         | <b>A828.1</b> | 1 unit(s) |  |

## Funnel range ROTILABO®

ROTH SELECTION. Material: HDPE.  
7 pieces, one for each funnel, Art. No. 0583.1 to 0589.1.

| Version                | Art. No.      | Pack Qty. |  |
|------------------------|---------------|-----------|--|
| ROTILABO® funnel range | <b>0594.1</b> | 1 set     |  |

# Funnels

## Funnels ROTILABO® disposable

ROTH SELECTION.

With ventilation grooves. 60° angle. Inner outlet Ø 4.3 mm.



### Transparent

Material: PP.

| Ø Inner upper (mm) | Ø Inner outlet (mm) | Height (mm) | Tube length (mm) | Suitable for          | Art. No.      | Pack Qty.   |  |
|--------------------|---------------------|-------------|------------------|-----------------------|---------------|-------------|--|
| 55                 | 4.3                 | 69          | 25               | paper filter-Ø 110 mm | <b>HY47.1</b> | 100 unit(s) |  |
| 65                 | 4.3                 | 78          | 30               | paper filter-Ø 125 mm | <b>HY48.1</b> | 100 unit(s) |  |

### Clear

Material: PS.

| Ø Inner upper (mm) | Ø Inner outlet (mm) | Height (mm) | Tube length (mm) | Suitable for          | Art. No.      | Pack Qty.   |  |
|--------------------|---------------------|-------------|------------------|-----------------------|---------------|-------------|--|
| 55                 | 4.3                 | 69          | 25               | Paper filter-Ø 110 mm | <b>HY45.1</b> | 100 unit(s) |  |
| 65                 | 4.3                 | 78          | 30               | Paper filter-Ø 125 mm | <b>HY46.1</b> | 100 unit(s) |  |

## Funnels ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes.

With lifting eye and side groove to bleed off the air from the container (except for Art. No. 0162.1). Suitable for foodstuffs.



**A 121 °C**

| Ø Outer upper (mm) | Ø Outlet (mm) | Height (mm) | Version               | Art. No.      | Pack Qty. |  |
|--------------------|---------------|-------------|-----------------------|---------------|-----------|--|
| 55                 | 8             | 65          | without handle/eyelet | <b>0162.1</b> | 1 unit(s) |  |
| 94                 | 13            | 120         | with handle           | <b>0164.1</b> | 1 unit(s) |  |
| 127                | 15            | 132         | with eyelet           | <b>0166.1</b> | 1 unit(s) |  |
| 154                | 15            | 152         | with eyelet           | <b>0167.1</b> | 1 unit(s) |  |
| 210                | 23            | 227         | with eyelet           | <b>0168.1</b> | 1 unit(s) |  |
| 257                | 23            | 265         | with eyelet           | <b>0169.1</b> | 1 unit(s) |  |

## Funnel assortment ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes.

Includes six funnels (Art. No. 0162.1, 0164.1, 0166.1, 0167.1, 0168.1 and 0169.1).

**A 121 °C**

| Version                     | Art. No.      | Pack Qty. |  |
|-----------------------------|---------------|-----------|--|
| Funnel assortment ROTILABO® | <b>5838.1</b> | 1 set     |  |



# Funnels



A 121 °C



## Funnels ROTILABO® made of fluoroplastic

ROTH SELECTION. Material: PTFE. Autoclavable: Yes.

Chemically inert, temperature stable from -200 to + 250 °C. With smooth, hydrophobic inner surfaces.

| Ø Outer upper (mm) | Ø Outlet (mm) | Height (mm) | Art. No.      | Pack Qty. |  |
|--------------------|---------------|-------------|---------------|-----------|--|
| 30                 | 8             | 50          | <b>NX77.1</b> | 1 unit(s) |  |
| 50                 | 10            | 85          | <b>NX78.1</b> | 1 unit(s) |  |
| 70                 | 12            | 110         | <b>NX79.1</b> | 1 unit(s) |  |
| 100                | 14            | 170         | <b>NX80.1</b> | 1 unit(s) |  |
| 150                | 16            | 220         | <b>NX81.1</b> | 1 unit(s) |  |



A 121 °C

## Funnels

Bürtle. Material: Stainless steel V2A (1.4301). Autoclavable: Yes.

With an air channel at the spout for unhindered flow.

| Ø Upper (mm) | Ø Outlet (mm) | Version        | Art. No.      | Pack Qty. |  |
|--------------|---------------|----------------|---------------|-----------|--|
| 60           | 13            | without handle | <b>1C82.1</b> | 1 unit(s) |  |
| 100          | 13            | with handle    | <b>1C83.1</b> | 1 unit(s) |  |
| 150          | 19            | with handle    | <b>1C84.1</b> | 1 unit(s) |  |
| 300          | 20            | with handle    | <b>1C85.1</b> | 1 unit(s) |  |



A 121 °C

## Funnels Micro

Bel-Art Products. Material: PP. Autoclavable: Yes.

For filling burettes and vessels with small opening.

| Vol. (ml) | Ø Upper (mm) | Ø Outlet (mm) | Tube length (mm) | Art. No.      | Pack Qty.  |  |
|-----------|--------------|---------------|------------------|---------------|------------|--|
| 2.6       | 24           | 5,4           | 16               | <b>YA66.1</b> | 12 unit(s) |  |
| 5.2       | 28           | 6             | 19               | <b>YA67.1</b> | 12 unit(s) |  |
| 8.1       | 35           | 8             | 23               | <b>YA68.1</b> | 12 unit(s) |  |
| 15.7      | 43           | 8             | 26               | <b>YA69.1</b> | 12 unit(s) |  |



A 121 °C

## Funnels ROTILABO® Mini

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes.

Seamless, matt polished. Suitable for foodstuffs.

| Ø Upper (mm) | Ø Outlet (mm) | Height (mm) | Art. No.      | Pack Qty.  |  |
|--------------|---------------|-------------|---------------|------------|--|
| 27           | 4             | 50          | <b>YY23.1</b> | 10 unit(s) |  |
| 27           | 6             | 50          | <b>YY24.1</b> | 10 unit(s) |  |
| 28           | 7             | 50          | <b>YY25.1</b> | 10 unit(s) |  |
| 29           | 8             | 50          | <b>YY26.1</b> | 10 unit(s) |  |
| 31           | 10            | 50          | <b>YY27.1</b> | 10 unit(s) |  |
| 32           | 11            | 50          | <b>YY28.1</b> | 10 unit(s) |  |
| 33           | 12            | 50          | <b>YY29.1</b> | 10 unit(s) |  |
| 34           | 13            | 50          | <b>YY30.1</b> | 10 unit(s) |  |
| 35           | 14            | 50          | <b>YY31.1</b> | 10 unit(s) |  |

# Funnels



**A 121 °C**

## Funnels with core

Lenz. Material: DURAN® glass. Autoclavable: Yes.  
Filler funnel with 60° angle.

| Standard ground glass joint core | Ø Outer upper (mm) | Art. No.      | Pack Qty. |  |
|----------------------------------|--------------------|---------------|-----------|--|
| 14/23                            | 45                 | <b>KX83.1</b> | 1 unit(s) |  |
| 19/26                            | 55                 | <b>EHL5.1</b> | 1 unit(s) |  |
| 24/29                            | 70                 | <b>EHL6.1</b> | 1 unit(s) |  |
| 29/32                            | 80                 | <b>KY11.1</b> | 1 unit(s) |  |



**A 121 °C**

## Powder funnel

Bohemia. Material: Borosilicate glass 3.3. Autoclavable: Yes.  
With short stem.

| Ø Outer upper (mm) | Ø Outer outlet (mm) | Height (mm) | Tube length (mm) | Art. No.      | Pack Qty.  |  |
|--------------------|---------------------|-------------|------------------|---------------|------------|--|
| 50                 | 20                  | 50          | 20               | <b>X669.1</b> | 10 unit(s) |  |
| 60                 | 20                  | 60          | 20               | <b>X670.1</b> | 10 unit(s) |  |
| 80                 | 25                  | 75          | 25               | <b>X672.1</b> | 10 unit(s) |  |
| 100                | 25                  | 95          | 25               | <b>X673.1</b> | 10 unit(s) |  |
| 120                | 30                  | 110         | 25               | <b>X674.1</b> | 10 unit(s) |  |
| 160                | 35                  | 137         | 30               | <b>NX16.1</b> | 1 unit(s)  |  |
| 200                | 38                  | 180         | 30               | <b>NX17.1</b> | 1 unit(s)  |  |



## Funnels ROTILABO® with standard taper joint

ROTH SELECTION. Material: PP.

For clean filling of flasks and equipment. Suitable for foodstuffs.

| Standard ground glass joint | Ø Upper (mm) | Height (mm) | Tube length (mm) | Art. No.      | Pack Qty. |  |
|-----------------------------|--------------|-------------|------------------|---------------|-----------|--|
| 14/23                       | 40           | 75          | 17               | <b>2060.1</b> | 1 unit(s) |  |
| 19/26                       | 50           | 95          | 23               | <b>2073.1</b> | 1 unit(s) |  |
| 29/32                       | 75           | 135         | 30               | <b>2075.1</b> | 1 unit(s) |  |

## Funnel range ROTILABO® with standard ground joint

ROTH SELECTION. Material: PP.

Includes 3 funnels (Art. No. 2060.1, 2073.1 and 2075.1).

| Version                             | Art. No.      | Pack Qty. |  |
|-------------------------------------|---------------|-----------|--|
| ROTILABO® ground joint funnel range | <b>2080.1</b> | 1 set     |  |



**A 121 °C**

## Powder funnel

Kartell. Material: PP. Autoclavable: Yes.

With a ribbed exterior to avoid pressure build-up. 60° inclination angle.

| Vol. (ml) | Ø Outer upper (mm) | Ø Outer outlet (mm) | Tube length (mm) | Art. No.      | Pack Qty. |  |
|-----------|--------------------|---------------------|------------------|---------------|-----------|--|
| 45        | 60                 | 15                  | 18               | <b>TT41.1</b> | 1 unit(s) |  |
| 94        | 80                 | 15                  | 22               | <b>TT42.1</b> | 1 unit(s) |  |
| 200       | 100                | 25                  | 25               | <b>TT43.1</b> | 1 unit(s) |  |
| 360       | 120                | 30                  | 27               | <b>TT44.1</b> | 1 unit(s) |  |
| 730       | 150                | 36                  | 40               | <b>TT45.1</b> | 1 unit(s) |  |
| 1,290     | 180                | 43                  | 49               | <b>TT46.1</b> | 1 unit(s) |  |



# Funnels



## Powder funnel ROTILABO®

ROTH SELECTION. Material: PP or HDPE, transparent.

With smooth surface. Autoclavable, except for Art. No. 2057.1.

Art. No. 2057.1 incl. sieve Ø 105 mm.

| Ø Inner upper (mm) | Ø Inner outlet (mm) | Height (mm) | Version | Art. No. | Pack Qty. |  |
|--------------------|---------------------|-------------|---------|----------|-----------|--|
| 65                 | 16                  | 68          | PP      | 2041.1   | 1 unit(s) |  |
| 80                 | 22                  | 75          | PP      | 2046.1   | 1 unit(s) |  |
| 100                | 23                  | 94          | PP      | 2047.1   | 1 unit(s) |  |
| 120                | 27                  | 110         | PP      | 2049.1   | 1 unit(s) |  |
| 150                | 29                  | 138         | PP      | 2051.1   | 1 unit(s) |  |
| 253                | 50                  | 190         | HDPE    | 2057.1   | 1 unit(s) |  |

## Powder funnel range ROTILABO®

ROTH SELECTION.

Includes 6 funnels (Art. No. 2041.1, 2046.1, 2047.1, 2049.1, 2051.1 and 2057.1).

| Version                       | Art. No. | Pack Qty. |  |
|-------------------------------|----------|-----------|--|
| ROTILABO® powder funnel range | 2020.1   | 1 set     |  |



## Drum funnels ROTILABO®

ROTH SELECTION.

Made of PP or HDPE, transparent. Art. No. HY26.1 and HY27.1 with side handle.

| Ø Upper (mm) | Ø Outlet (mm) | Height (mm) | Version | Art. No. | Pack Qty. |  |
|--------------|---------------|-------------|---------|----------|-----------|--|
| 200          | 22            | 200         | PP      | HY26.1   | 1 unit(s) |  |
| 250          | 30            | 260         | PP      | HY27.1   | 1 unit(s) |  |
| 350          | 35            | 440         | PP      | HY28.1   | 1 unit(s) |  |
| 400          | 42            | 365         | HDPE    | HY29.1   | 1 unit(s) |  |

## Accessories sieve insert for ROTILABO® drum funnel

ROTH SELECTION. Material: Stainless steel, Aluminium.

Suitable for Art. No. HY26.1 and HY27.1.

| Ø (mm) | Mesh size (mm) | Art. No. | Pack Qty. |  |
|--------|----------------|----------|-----------|--|
| 50     | 0.45           | HY30.1   | 1 unit(s) |  |



**A** 121 °C

## Analytical funnels ROTILABO® Urbanti

ROTH SELECTION. Material: PMP. Autoclavable: Yes.

With spiral rib on the inside of the funnel to accelerate filtration speed.

With ventilation ribs. Temperature stable up to 170 °C.

| Ø Inner upper (mm) | Ø Inner outlet (mm) | Tube length (mm) | Art. No. | Pack Qty. |
|--------------------|---------------------|------------------|----------|-----------|
| 51                 | 3                   | 150              | EA42.1   | 1 unit(s) |
| 70                 | 3                   | 150              | EA43.1   | 1 unit(s) |
| 100                | 7                   | 108              | EA44.1   | 1 unit(s) |
| 140                | 10                  | 132              | EA45.1   | 1 unit(s) |
| 196                | 20                  | 160              | EA46.1   | 1 unit(s) |

# Funnels

## Analytical funnels ROTILABO®

ROTH SELECTION. Material: Borosilicate glass 3.3. Autoclavable: Yes.  
Acc. to ISO 4798. Stem Ø 8 mm. Length 150 mm. For quick filtration.



| Ø Inner upper (mm) | Ø Outlet (mm) | Height (mm) | Suitable for              | Art. No.      | Pack Qty. |  |
|--------------------|---------------|-------------|---------------------------|---------------|-----------|--|
| 55                 | 8             | 150         | paper filter-Ø 70-90 mm   | <b>NT85.1</b> | 1 unit(s) |  |
| 75                 | 8             | 150         | paper filter-Ø 110-125 mm | <b>NT86.1</b> | 1 unit(s) |  |
| 100                | 8             | 150         | paper filter-Ø 150-185 mm | <b>NT87.1</b> | 1 unit(s) |  |

## Analytical funnels

DWK Life Sciences. Material: DURAN® glass. Autoclavable: Yes.  
Ribbed.



| Ø Outer upper (mm) | Ø Outer outlet (mm) | Tube length (mm) | Suitable for              | Art. No.      | Pack Qty. |  |
|--------------------|---------------------|------------------|---------------------------|---------------|-----------|--|
| 70                 | 8                   | 70               | paper filter-Ø 110-125 mm | <b>K243.1</b> | 1 unit(s) |  |
| 80                 | 10                  | 80               | paper filter-Ø 125-150 mm | <b>K244.1</b> | 1 unit(s) |  |
| 100                | 10                  | 100              | paper filter-Ø 150-185 mm | <b>K245.1</b> | 1 unit(s) |  |
| 150                | 16                  | 150              | paper filter-Ø 240-270 mm | <b>K246.1</b> | 1 unit(s) |  |
| 200                | 26                  | 175              | paper filter-Ø 320-385 mm | <b>K247.1</b> | 1 unit(s) |  |

## Analytical funnels ROTILABO®

ROTH SELECTION. Material: PP, transparent. Autoclavable: Yes.



| Ø Outer upper (mm) | Ø Outlet (mm) | Tube length (mm) | Suitable for              | Art. No.      | Pack Qty. |  |
|--------------------|---------------|------------------|---------------------------|---------------|-----------|--|
| 60                 | 8             | 150              | paper filter-Ø 70-90 mm   | <b>9518.1</b> | 1 unit(s) |  |
| 80                 | 8             | 150              | paper filter-Ø 110-125 mm | <b>9519.1</b> | 1 unit(s) |  |
| 100                | 8             | 150              | paper filter-Ø 150-185 mm | <b>9520.1</b> | 1 unit(s) |  |



## Funnel holder

Bürkle. Material: PP, Steel chrome-plated.

For funnel Ø of 40 to 180 mm (for smaller funnels, the supplied reducing inserts can be used) or for sediment funnels from Ø 90 mm. Infinitely adjustable thanks to magnet concept. Height is adjustable up to 450 mm. With 2 or 4 spaces on the base.



| Version                   | Art. No. | Pack Qty. |  |
|---------------------------|----------|-----------|--|
| for 2 funnels             | HCC8.1   | 1 unit(s) |  |
| for 4 funnels             | HCC9.1   | 1 unit(s) |  |
| for 2 sedimentation cones | HCE0.1   | 1 unit(s) |  |

## Accessories

Bürkle.


| Type                                              | Art. No. | Pack Qty. |  |
|---------------------------------------------------|----------|-----------|--|
| Holder (2 holding rings) for 1 sedimentation cone | HCE2.1   | 1 unit(s) |  |
| Funnel holder with reduction insert               | HCE1.1   | 1 unit(s) |  |

## Funnel holder ROTILABO®

ROTH SELECTION. Material: PP.

With stand clamp for mounting on stands of Ø 8 to 14 mm. Suitable for funnel Ø of 50 to 120 mm and hose outer Ø of max. 39 mm.

With replaceable inner ring for varying funnel sizes.



| Version       | Art. No. | Pack Qty. |  |
|---------------|----------|-----------|--|
| for 1 funnel  | 9768.1   | 1 unit(s) |  |
| for 2 funnels | 9769.1   | 1 unit(s) |  |

**Алматы** (7273)495-231  
**Ангарск** (3955)60-70-56  
**Архангельск** (8182)63-90-72  
**Астрахань** (8512)99-46-04  
**Барнаул** (3852)73-04-60  
**Белгород** (4722)40-23-64  
**Благовещенск** (4162)22-76-07  
**Брянск** (4832)59-03-52  
**Владивосток** (423)249-28-31  
**Владикавказ** (8672)28-90-48  
**Владимир** (4922)49-43-18  
**Волгоград** (844)278-03-48  
**Вологда** (8172)26-41-59  
**Воронеж** (473)204-51-73  
**Екатеринбург** (343)384-55-89

**Иваново** (4932)77-34-06  
**Ижевск** (3412)26-03-58  
**Иркутск** (395)279-98-46  
**Казань** (843)206-01-48  
**Калининград** (4012)72-03-81  
**Калуга** (4842)92-23-67  
**Кемерово** (3842)65-04-62  
**Киров** (8332)68-02-04  
**Коломна** (4966)23-41-49  
**Кострома** (4942)77-07-48  
**Краснодар** (861)203-40-90  
**Красноярск** (391)204-63-61  
**Курск** (4712)77-13-04  
**Курган** (3522)50-90-47  
**Липецк** (4742)52-20-81

**Россия** (495)268-04-70

**Магнитогорск** (3519)55-03-13  
**Москва** (495)268-04-70  
**Мурманск** (8152)59-64-93  
**Набережные Челны** (8552)20-53-41  
**Нижний Новгород** (831)429-08-12  
**Новокузнецк** (3843)20-46-81  
**Ноябрьск** (3496)41-32-12  
**Новосибирск** (383)227-86-73  
**Омск** (3812)21-46-40  
**Орел** (4862)44-53-42  
**Оренбург** (3532)37-68-04  
**Пенза** (8412)22-31-16  
**Петрозаводск** (8142)55-98-37  
**Псков** (8112)59-10-37  
**Пермь** (342)205-81-47

**Казахстан** (772)734-952-31

**Ростов-на-Дону** (863)308-18-15  
**Рязань** (4912)46-61-64  
**Самара** (846)206-03-16  
**Саранск** (8342)22-96-24  
**Санкт-Петербург** (812)309-46-40  
**Саратов** (845)249-38-78  
**Севастополь** (8692)22-31-93  
**Симферополь** (3652)67-13-56  
**Смоленск** (4812)29-41-54  
**Сочи** (862)225-72-31  
**Ставрополь** (8652)20-65-13  
**Сыктывкар** (8212)25-95-17  
**Тамбов** (4752)50-40-97  
**Тургут** (3462)77-98-35  
**Тверь** (4822)63-31-35

**Киргизия** (996)312-96-26-47

**Тольятти** (8482)63-91-07  
**Томск** (3822)98-41-53  
**Тула** (4872)74-02-29  
**Тюмень** (3452)66-21-18  
**Ульяновск** (8422)24-23-59  
**Улан-Удэ** (3012)59-97-51  
**Уфа** (347)229-48-12  
**Хабаровск** (4212)92-98-04  
**Чебоксары** (8352)28-53-07  
**Челябинск** (351)202-03-61  
**Череповец** (8202)49-02-64  
**Чита** (3022)38-34-83  
**Якутск** (4112)23-90-97  
**Ярославль** (4852)69-52-93