

Стаканы лабораторные стеклянные CARL ROTH DURAN, DURAN Super Duty, ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47



Glass beakers

Beaker DURAN®

DWK Life Sciences. Material: DURAN® glass. Autoclavable: Yes.

Graduated, with a spout (Art. No. C108.1 and C109.1 not graduated).

According to DIN 12331, ISO 3819 (Art. No. C122.1 not according to DIN). With retrace code.



A121 °C



Low form

Vol. (ml)	Ø ext. (mm)	Height (mm)	Standards	Art. No.	Pack Qty.
5	22	30	DIN 12331, ISO 3819	C108.1	10 unit(s)
10	26	35	DIN 12331, ISO 3819	C109.1	10 unit(s)
25	34	50	DIN 12331, ISO 3819	C110.1	10 unit(s)
50	42	60	DIN 12331, ISO 3819	C111.1	10 unit(s)
100	50	70	DIN 12331, ISO 3819	C112.1	10 unit(s)
150	60	80	DIN 12331, ISO 3819	C113.1	10 unit(s)
250	70	95	DIN 12331, ISO 3819	C114.1	10 unit(s)
400	80	110	DIN 12331, ISO 3819	C115.1	10 unit(s)
600	90	125	DIN 12331, ISO 3819	C116.1	10 unit(s)
800	100	135	DIN 12331, ISO 3819	C117.1	10 unit(s)
1,000	105	145	DIN 12331, ISO 3819	C118.1	1 unit(s)
2,000	132	185	DIN 12331, ISO 3819	C119.1	1 unit(s)
3,000	152	210	DIN 12331, ISO 3819	C120.1	1 unit(s)
5,000	170	270	DIN 12331, ISO 3819	C121.1	1 unit(s)
10,000	217	350	not in acc. with DIN	C122.1	1 unit(s)

High form

Vol. (ml)	Ø ext. (mm)	Height (mm)	Standards	Art. No.	Pack Qty.
50	38	70	DIN 12331, ISO 3819	C123.1	10 unit(s)
100	48	80	DIN 12331, ISO 3819	C124.1	10 unit(s)
150	54	95	DIN 12331, ISO 3819	C125.1	10 unit(s)
250	60	120	DIN 12331, ISO 3819	C126.1	10 unit(s)
400	70	130	DIN 12331, ISO 3819	C127.1	10 unit(s)
600	80	150	DIN 12331, ISO 3819	C128.1	10 unit(s)
800	90	175	DIN 12331, ISO 3819	C129.1	10 unit(s)
1,000	95	180	DIN 12331, ISO 3819	C130.1	1 unit(s)
2,000	120	240	DIN 12331, ISO 3819	C131.1	1 unit(s)
3,000	135	280	DIN 12331, ISO 3819	C132.1	1 unit(s)

Beaker DURAN® Super Duty

DWK Life Sciences. Material: DURAN® glass. Autoclavable: Yes.

- Proven DURAN® properties
- Glass type I, Neutral glass according to USP, EP and JP
- High mechanical stability
- Reinforced, stable rim
- Regular wall-thickness
- Autoclavable, sterilisable
- With retrace code
- Easily readable scale
- Fired-on highly durable labelling



A121 °C



Notes for use:

To avoid thermal stresses in the glass, uniform and slow heating of Super Duty products is recommended. When working at very high temperatures, or if rapid temperature changes are expected, standard DURAN® beakers and Erlenmeyer flasks should be used, as these have excellent thermal shock resistance. However, their mechanical sturdiness is limited compared with the Super Duty range.

Glass beakers

Low form

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
150	60	80	PC33.1	10 unit(s)	
250	70	95	PC34.1	10 unit(s)	
400	80	110	PC35.1	10 unit(s)	
600	90	125	PC36.1	10 unit(s)	
1,000	105	145	PC37.1	10 unit(s)	
2,000	132	185	PC38.1	10 unit(s)	
5,000	170	270	PC39.1	1 unit(s)	

High form

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
150	54	95	CHL2.1	10 unit(s)	
250	60	120	CHL3.1	10 unit(s)	
600	80	150	CHL4.1	10 unit(s)	

Beaker ROTILABO®

ROTH SELECTION. Material: Borosilicate glass 3.3. Autoclavable: Yes. Graduated, with spout (Art. No. X685.1, X686.1, X687.1, and X700.1 not graduated). According to DIN 12331, ISO 3819 (Art. No. ACK7.1 not according to DIN).



A 121 °C

Low form

Vol. (ml)	Ø ext. (mm)	Height (mm)	Standards	Art. No.	Pack Qty.	
5	22	30	DIN 12331, ISO 3819	X685.1	10 unit(s)	
10	26	35	DIN 12331, ISO 3819	X686.1	10 unit(s)	
25	34	50	DIN 12331, ISO 3819	X687.1	10 unit(s)	
50	42	60	DIN 12331, ISO 3819	X688.1	10 unit(s)	
100	50	70	DIN 12331, ISO 3819	X689.1	10 unit(s)	
150	60	80	DIN 12331, ISO 3819	X690.1	10 unit(s)	
250	70	95	DIN 12331, ISO 3819	X691.1	10 unit(s)	
400	80	110	DIN 12331, ISO 3819	X692.1	10 unit(s)	
600	90	125	DIN 12331, ISO 3819	X693.1	10 unit(s)	
800	100	135	DIN 12331, ISO 3819	X694.1	10 unit(s)	
1,000	105	145	DIN 12331, ISO 3819	X695.1	10 unit(s)	
2,000	130	185	DIN 12331, ISO 3819	X696.2	1 unit(s)	
3,000	150	210	DIN 12331, ISO 3819	X697.1	1 unit(s)	
5,000	170	270	DIN 12331, ISO 3819	X698.1	1 unit(s)	
10,000	217	350	not in acc. with DIN	ACK7.1	1 unit(s)	

High form

Vol. (ml)	Ø ext. (mm)	Height (mm)	Standards	Art. No.	Pack Qty.	
25	30	55	DIN 12331, ISO 3819	X700.1	10 unit(s)	
50	38	70	DIN 12331, ISO 3819	X701.1	10 unit(s)	
100	48	80	DIN 12331, ISO 3819	X702.1	10 unit(s)	
150	54	95	DIN 12331, ISO 3819	X703.1	10 unit(s)	
250	60	120	DIN 12331, ISO 3819	X704.1	10 unit(s)	
400	70	130	DIN 12331, ISO 3819	X705.1	10 unit(s)	
600	80	150	DIN 12331, ISO 3819	X706.1	10 unit(s)	
800	90	175	DIN 12331, ISO 3819	X707.1	10 unit(s)	
1,000	95	180	DIN 12331, ISO 3819	X708.1	10 unit(s)	
2,000	120	240	DIN 12331, ISO 3819	X709.2	1 unit(s)	
3,000	135	280	DIN 12331, ISO 3819	X710.2	1 unit(s)	

Алматы (7273)495-231	Иваново (4932)77-34-06	Магнитогорск (3519)55-03-13	Ростов-на-Дону (863)308-18-15	Тольятти (8482)63-91-07
Ангарск (3955)60-70-56	Ижевск (3412)26-03-58	Москва (495)268-04-70	Рязань (4912)46-61-64	Томск (3822)98-41-53
Архангельск (8182)63-90-72	Иркутск (395)279-98-46	Мурманск (8152)59-64-93	Самара (846)206-03-16	Тула (4872)74-02-29
Астрахань (8512)99-46-04	Казань (843)206-01-48	Набережные Челны (8552)20-53-41	Саранск (8342)22-96-24	Тюмень (3452)66-21-18
Барнаул (3852)73-04-60	Калининград (4012)72-03-81	Нижний Новгород (831)429-08-12	Санкт-Петербург (812)309-46-40	Ульяновск (8422)24-23-59
Белгород (4722)40-23-64	Калуга (4842)92-23-67	Новокузнецк (3843)20-46-81	Саратов (845)249-38-78	Улан-Удэ (3012)59-97-51
Благовещенск (4162)22-76-07	Кемерово (3842)65-04-62	Ноябрьск (3496)41-32-12	Севастополь (8692)22-31-93	Уфа (347)229-48-12
Брянск (4832)59-03-52	Киров (8332)68-02-04	Новосибирск (383)227-86-73	Симферополь (3652)67-13-56	Хабаровск (4212)92-98-04
Владивосток (423)249-28-31	Коломна (4966)23-41-49	Омск (3812)21-46-40	Смоленск (4812)29-41-54	Чебоксары (8352)28-53-07
Владикавказ (8672)28-90-48	Кострома (4942)77-07-48	Орел (4862)44-53-42	Сочи (862)225-72-31	Челябинск (351)202-03-61
Владимир (4922)49-43-18	Краснодар (861)203-40-90	Оренбург (3532)37-68-04	Ставрополь (8652)20-65-13	Череповец (8202)49-02-64
Волгоград (844)278-03-48	Красноярск (391)204-63-61	Пенза (8412)22-31-16	Сыктывкар (8212)25-95-17	Чита (3022)38-34-83
Вологда (8172)26-41-59	Курск (4712)77-13-04	Петрозаводск (8142)55-98-37	Тамбов (4752)50-40-97	Якутск (4112)23-90-97
Воронеж (473)204-51-73	Курган (3522)50-90-47	Псков (8112)59-10-37	Сургут (3462)77-98-35	Ярославль (4852)69-52-93
Екатеринбург (343)384-55-89	Липецк (4742)52-20-81	Пермь (342)205-81-47		
	Россия (495)268-04-70	Казахстан (772)734-952-31	Киргизия (996)312-96-26-47	