

Коннекторы для шлангов, клапаны дренажные CARL ROTH ROTILABO, ROTILABO Mini, FLUXX, Keck

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93


Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47



Hose connectors



Hose connectors ROTILABO® Mini

ROTH SELECTION. Material: PP transparent.


Not a medical device / Not an IVD product

(1) Straight shape

Suitable for	Art. No.	Pack Qty.	
Hose inner Ø 1.6 mm	CT33.1	10 unit(s)	
Hose inner Ø 2.4 mm	CT34.1	10 unit(s)	
Hose inner Ø 3.2 mm	CT35.1	10 unit(s)	

(2) Y-shape

Suitable for	Art. No.	Pack Qty.	
Hose inner Ø 1.6 mm	CT36.1	10 unit(s)	
Hose inner Ø 2.4 mm	CT37.1	10 unit(s)	
Hose inner Ø 3.2 mm	CT38.1	10 unit(s)	



Hose reducers ROTILABO® Mini

ROTH SELECTION. Material: PP transparent.


Not a medical device / Not an IVD product

(1) Straight shape

Suitable for	Art. No.	Pack Qty.	
Hose inner Ø 1.6/4.0 mm	CT42.1	10 unit(s)	
Hose inner Ø 2.4/4.0 mm	CT43.1	10 unit(s)	
Hose inner Ø 3.2/4.0 mm	CT44.1	10 unit(s)	
Hose inner Ø 2.4/4.8 mm	CT46.1	10 unit(s)	
Hose inner Ø 3.2/4.8 mm	CT47.1	10 unit(s)	

T-shape

Pic.	Suitable for	Art. No.	Pack Qty.	
(3)	Hose inner Ø 1.6/2.4/1.6 mm	CT52.1	10 unit(s)	
(3)	Hose inner Ø 1.6/3.2/1.6 mm	CT53.1	10 unit(s)	
(3)	Hose inner Ø 2.4/3.2/2.4 mm	CT55.1	10 unit(s)	
(3)	Hose inner Ø 3.2/1.6/3.2 mm	CT56.1	10 unit(s)	




Luer hose connectors ROTILABO®

ROTH SELECTION. Material: PP transparent.

Not a medical device / Not an IVD product

Pic.	Suitable for	Art. No.	Pack Qty.	
(1)	Luer Lock Female / hose inner Ø 1.6 mm	CT58.1	10 unit(s)	
(1)	Luer Lock Female / hose inner Ø 2.4 mm	CT59.1	10 unit(s)	
(1)	Luer Lock Female / hose inner Ø 3.2 mm	CT60.1	10 unit(s)	
(1)	Luer Lock Female / hose inner Ø 4.0 mm	CT61.1	10 unit(s)	
(2)	Luer Lock Male / hose inner Ø 1.6 mm	CT62.1	10 unit(s)	
(2)	Luer Lock Male / hose inner Ø 2.4 mm	CT63.1	10 unit(s)	
(2)	Luer Lock Male / hose inner Ø 3.2 mm	CT64.1	10 unit(s)	
(2)	Luer Lock Male / hose inner Ø 4.0 mm	CT65.1	10 unit(s)	




Luer connectors ROTILABO®

ROTH SELECTION. Material: PP transparent.

Not a medical device / Not an IVD product

Pic.	Suitable for	Art. No.	Pack Qty.	
(1)	Luer Lock female / Luer Male	CT66.1	10 unit(s)	
(2)	Luer Lock Male / Luer Lock Male	CT67.1	10 unit(s)	
(3)	Luer Female / Luer Female	CT68.1	10 unit(s)	



Luer stopper ROTILABO®

ROTH SELECTION. Material: PP transparent.

Not a medical device / Not an IVD product

Pic.	Suitable for	Art. No.	Pack Qty.	
(1)	Luer Lock female	CT69.1	10 unit(s)	
(2)	Luer Lock male	CT70.1	10 unit(s)	
(3)	Luer female	CT71.1	10 unit(s)	

Hose connectors



A 121 °C

Hose connectors ROTILABO® T-shape

ROTH SELECTION. Material: PP natural, translucent.
Autoclavable: Yes.

Purpose: for any tubing system and vacuum line

In vitro diagnostic medical device according to IVDD

Ø ext. (mm)	Suitable for	Art. No.	Pack Qty.	
4	Hose inner Ø 3 mm	E763.1	10 unit(s)	
6	Hose inner Ø 4–5 mm	E764.1	10 unit(s)	
8	Hose inner Ø 6–7 mm	E765.1	10 unit(s)	
10	Hose inner Ø 8–9 mm	E766.1	10 unit(s)	
12	Hose inner Ø 10–11 mm	E767.1	10 unit(s)	
14	Hose inner Ø 12–13 mm	E768.1	10 unit(s)	
16	Hose inner Ø 14–15 mm	E769.1	10 unit(s)	



A 121 °C

Hose connectors ROTILABO® cross shape

ROTH SELECTION. Material: PP natural, translucent.
Autoclavable: Yes.

Purpose: for any tubing system and vacuum line

In vitro diagnostic medical device according to IVDD

Ø ext. (mm)	Suitable for	Art. No.	Pack Qty.	
4	Hose inner Ø 3 mm	E777.1	10 unit(s)	
6	Hose inner Ø 4–5 mm	E778.1	10 unit(s)	
8	Hose inner Ø 6–7 mm	E779.1	10 unit(s)	
10	Hose inner Ø 8–9 mm	E780.1	10 unit(s)	
12	Hose inner Ø 10–11 mm	E781.1	10 unit(s)	
14	Hose inner Ø 12–13 mm	E782.1	10 unit(s)	



A 121 °C

Hose connectors ROTILABO® angle shape

ROTH SELECTION. Material: PP natural, translucent.
Autoclavable: Yes.

Purpose: for any tubing system and vacuum line

In vitro diagnostic medical device according to IVDD

Ø ext. (mm)	Suitable for	Art. No.	Pack Qty.	
4	Hose inner Ø 3 mm	E784.1	10 unit(s)	
6	Hose inner Ø 4–5 mm	E785.1	10 unit(s)	
8	Hose inner Ø 6–7 mm	E786.1	10 unit(s)	
10	Hose inner Ø 8–9 mm	E787.1	10 unit(s)	
12	Hose inner Ø 10–11 mm	E788.1	10 unit(s)	
14	Hose inner Ø 12–13 mm	E789.1	10 unit(s)	
16	Hose inner Ø 14–15 mm	E790.1	10 unit(s)	



A 121 °C

Hose connectors ROTILABO® Y-shape

ROTH SELECTION. Autoclavable: Yes.

Purpose: for any tubing system and vacuum line

In vitro diagnostic medical device according to IVDD

(1) With straight ends

Material: PP natural, translucent.

Ø ext. (mm)	Suitable for	Art. No.	Pack Qty.	
4	Hose inner Ø 3 mm	E770.1	10 unit(s)	
6	Hose inner Ø 4–5 mm	E771.1	10 unit(s)	
8	Hose inner Ø 6–7 mm	E772.1	10 unit(s)	
10	Hose inner Ø 8–9 mm	E773.1	10 unit(s)	
12	Hose inner Ø 10–11 mm	E774.1	10 unit(s)	
14	Hose inner Ø 12–13 mm	E775.1	10 unit(s)	
16	Hose inner Ø 14–15 mm	E776.1	10 unit(s)	



Hose connectors

(2) With conical ends

Material: PP.

Ø ext. (mm)	Suitable for	Art. No.	Pack Qty.	
4 / 5	Hose inner Ø 3–5 mm	TT50.1	10 unit(s)	
7 / 9	Hose inner Ø 7–9 mm	TT52.1	10 unit(s)	
10 / 11,5	Hose inner Ø 9–11 mm	TT53.1	10 unit(s)	
12 / 14	Hose inner Ø 11–13 mm	TT54.1	10 unit(s)	
14 / 15,5	Hose inner Ø 13–15 mm	TT55.1	10 unit(s)	
5,5 / 7	Hose inner Ø 5–7 mm	TT51.1	10 unit(s)	



Hose connectors ROTILABO® straight shape

ROTH SELECTION. Material: HDPE white.

Not a medical device / Not an IVD product

Ø Internal (mm)	Suitable for	Art. No.	Pack Qty.	
7	Hose inner Ø 10 mm	E797.1	10 unit(s)	
9	Hose inner Ø 12 mm	E798.1	10 unit(s)	
15	Hose inner Ø 19 mm	E799.1	10 unit(s)	
20	Hose inner Ø 25 mm	E800.1	10 unit(s)	



Hose connectors made of glass

ROTH SELECTION. Material: DURAN® glass.

Not a medical device / Not an IVD product

(1) T-shape

With straight ends.

Ø ext. (mm)	Length (mm)	Suitable for	Art. No.	Pack Qty.	
8	115	Hose inner Ø 6–8 mm	EKH1.1	1 unit(s)	
10	115	Hose inner Ø 8–10 mm	EKH2.1	1 unit(s)	

(2) Y-shape

With straight ends.

Ø ext. (mm)	Length (mm)	Suitable for	Art. No.	Pack Qty.	
8	115	Hose inner Ø 6–8 mm	EKH3.1	1 unit(s)	
10	115	Hose inner Ø 8–10 mm	EKH4.1	1 unit(s)	

(3) Straight

With conical ends.

Ø ext. (mm)	Length (mm)	Suitable for	Art. No.	Pack Qty.	
6–14	110	Hose inner Ø 6–14 mm	EKH5.1	1 unit(s)	



A 121 °C

Hose connectors ROTILABO® straight shape with conical ends

ROTH SELECTION. Material: PP white. Autoclavable: Yes.

Purpose: for any tubing system and vacuum line

In vitro diagnostic medical device according to IVDD

Ø Internal (mm)	Suitable for	Art. No.	Pack Qty.	
2	Hose inner Ø 4–5 mm	E791.1	10 unit(s)	
3	Hose inner Ø 5–7 mm	E792.1	10 unit(s)	
4	Hose inner Ø 8–10 mm	E793.1	10 unit(s)	
6	Hose inner Ø 10–12 mm	E794.1	10 unit(s)	
8	Hose inner Ø 12–14 mm	E795.1	10 unit(s)	
10	Hose inner Ø 14–16 mm	E796.1	10 unit(s)	

Hose connectors



A 121 °C

Hose reducers ROTILABO®

ROTH SELECTION. Material: PP white. Autoclavable: Yes.

Purpose: for any tubing system and vacuum line

With conical ends. Temperature stable up to +130 °C. Length 87 mm.

In vitro diagnostic medical device according to IVDD

Suitable for	Art. No.	Pack Qty.	
Hose inner Ø 8–12/4–8 mm	E801.1	10 unit(s)	
Hose inner Ø 12–16/4–8 mm	E802.1	10 unit(s)	
Hose inner Ø 12–16/8–12 mm	E803.1	10 unit(s)	



A 121 °C

Hose connectors straight shape universal

ROTH SELECTION. Autoclavable: Yes.

With conical ends.

Not a medical device / Not an IVD product

Pic.	Materialspe- cification	Suitable for	Art. No.	Pack Qty.	
(1)	PP	Hose inner Ø 4–17/8–17 mm	Y523.1	10 unit(s)	
(2)	PVDF	Hose inner Ø 4–17/8–17 mm	Y524.1	10 unit(s)	



Hose couplings ROTILABO®

ROTH SELECTION. Material: PE natural, translucent.

Purpose: for joining hoses made of glass, rubber or plastic.

Excellent sealing properties and can be disassembled effortlessly at any time. Conical, temperature stable up to +100 °C.

In vitro diagnostic medical device according to IVDD

Ø Internal (mm)	Suitable for	Art. No.	Pack Qty.	
2	Hose inner Ø 4–5 mm	E806.1	10 unit(s)	
3	Hose inner Ø 6–7 mm	E807.1	10 unit(s)	
4	Hose inner Ø 8–9 mm	E808.1	10 unit(s)	
6	Hose inner Ø 9–12 mm	E809.1	10 unit(s)	
7	Hose inner Ø 11–14 mm	E810.1	10 unit(s)	



Non-return valves ROTILABO®

ROTH SELECTION. Material: PP natural colours, translucent.

Suitable for liquids. Consists of two parts. Prevents pressure drop and backflow. With conical ends. Temperature stable up to +100 °C.

Not a medical device / Not an IVD product

Length (mm)	Suitable for	Art. No.	Pack Qty.	
60	Hose inner Ø 7.8–10 mm	E804.1	5 unit(s)	
70	Hose inner Ø 9.5–12.5 mm	E805.1	5 unit(s)	



Hose connectors



Quick-release coupling

ROTH SELECTION. Material: PP, EPDM, Stainless steel.

For easy, leak-free connection and disconnection of hoses. With hose connector. Optionally with shut-off, stops the flow automatically upon disconnection. Can be sterilised with gamma rays.

- NW flow opening: 3.2 mm
- Temperature range: 0 to +71 °C
- Operating pressure: Vacuum up to 8.3 bar

Not a medical device / Not an IVD product

(1) Female connection straight

Suitable for	Shut-off	Art. No.	Pack Qty.	
Hose inner Ø 3.2 mm	Yes	HAH8.1	1 unit(s)	
Hose inner Ø 6.4 mm	Yes	HAK0.1	1 unit(s)	

(2) Male connection straight

Suitable for	Shut-off	Art. No.	Pack Qty.	
Hose inner Ø 1.6 mm	Yes	HAK2.1	1 unit(s)	
Hose inner Ø 3.2 mm	No	HAK3.1	1 unit(s)	
Hose inner Ø 3.2 mm	Yes	HAK4.1	1 unit(s)	
Hose inner Ø 6.4 mm	Yes	HAK6.1	1 unit(s)	

(3) Male connection angled

Suitable for	Shut-off	Art. No.	Pack Qty.	
Hose inner Ø 6.4 mm	No	KA53.1	1 unit(s)	
Hose inner Ø 6.4 mm	Yes	HAK9.1	1 unit(s)	



A121 °C

Luer hose connectors made of metal with conical ends

ROTH SELECTION. Material: Brass nickel-plated, bright chrome-plated. Autoclavable: Yes.

For connecting needles, syringes and taps to hoses.

LLM: Luer Lock Male, LLF: Luer Lock Female

Not a medical device / Not an IVD product

Pic.	Suitable for	Art. No.	Pack Qty.	
(1)	LLM/hose inner Ø 2–6 mm	1HCC.1	1 unit(s)	
(2)	LLM/hose inner Ø 5–9 mm	1HCE.1	1 unit(s)	
(3)	LLF/hose inner Ø 2–6 mm	1HCH.1	1 unit(s)	
(4)	LLF/hose inner Ø 5–9 mm	1HCK.1	1 unit(s)	



A121 °C

Luer hose connectors made of metal with straight ends

ROTH SELECTION. Material: Brass nickel-plated, bright chrome-plated. Autoclavable: Yes.

For connecting needles, syringes and taps to hoses.

LLM: Luer Lock Male, LLF: Luer Lock Female

Not a medical device / Not an IVD product

Pic.	Suitable for	Art. No.	Pack Qty.	
(1)	LLM/hose inner Ø 1.5 mm	1HC7.1	1 unit(s)	
(2)	LLM/hose inner Ø 3.0 mm	1HC8.1	1 unit(s)	
(3)	LLF/hose inner Ø 1.5 mm	1HC9.1	1 unit(s)	
(4)	LLF/hose inner Ø 3.0 mm	1HCA.1	1 unit(s)	

Hose connectors



A 121 °C

Luer valve made of metal

ROTH SELECTION. Material: Brass nickel-plated, bright chrome-plated. Autoclavable: Yes.

Suitable for pressures up to 3 bar.

LLM: Luer Lock Male, LLF: Luer Lock Female

Not a medical device / Not an IVD product

Pic.	Version	Suitable for	Art. No.	Pack Qty.	
(1)	2-way tap	LLF / LLM	1HCL.1	1 unit(s)	
(2)	3-way tap	LLF / LLM / LLM	1HCN.1	1 unit(s)	




Hose valve ROTILABO®

ROTH SELECTION. Material: LDPE natural colours.

Temperature range: -20 to +70 °C. Withstands pressures up to 0.4 bar.

Not a medical device / Not an IVD product

Ø ext. (mm)	Length (mm)	Suitable for	Art. No.	Pack Qty.	
10	70	Hose inner Ø 8–9 mm	1622.1	5 unit(s)	
13	95	Hose inner Ø 11–12 mm	1623.1	5 unit(s)	



Hose valve

Bürkle.

For the transfer of liquids or gases. Max. pressure 1 bar. Operating temperature +5 to +40 °C. Available in various materials: PVDF (natural colours) offers better chemical resistance than PP/PE (red/white).

Not a medical device / Not an IVD product

2-way tap

Materialspecification	Ø Internal (mm)	Suitable for	Art. No.	Pack Qty.	
PP/PE	4	Hose inner Ø 5–7 mm	1015.1	1 unit(s)	
PP/PE	6	Hose inner Ø 7–9 mm	1016.1	1 unit(s)	
PP/PE	8	Tubing of internal Ø 9–11 mm	1017.1	1 unit(s)	
PP/PE	8	Hose inner Ø 11–13 mm	CPH6.1	1 unit(s)	
PVDF	4	Hose inner Ø 5–7 mm	1761.1	1 unit(s)	
PVDF	6	Hose inner Ø 7–9 mm	1762.1	1 unit(s)	
PVDF	8	Hose inner Ø 9–11 mm	1763.1	1 unit(s)	
PVDF	8	Hose inner Ø 11–13 mm	CPH7.1	1 unit(s)	

3-way tap

Materialspecification	Ø Internal (mm)	Suitable for	Art. No.	Pack Qty.	
PP/PE	4	Hose inner Ø 5–7 mm	1018.1	1 unit(s)	
PP/PE	6	Hose inner Ø 7–9 mm	1019.1	1 unit(s)	
PP/PE	8	Hose inner Ø 9–11 mm	1020.1	1 unit(s)	
PP/PE	8	Hose inner Ø 11–13 mm	CPH8.1	1 unit(s)	
PVDF	4	Hose inner Ø 5–7 mm	1766.1	1 unit(s)	
PVDF	6	Hose inner Ø 7–9 mm	1767.1	1 unit(s)	
PVDF	8	Hose inner Ø 9–11 mm	1768.1	1 unit(s)	
PVDF	8	Hose inner Ø 11–13 mm	CPH9.1	1 unit(s)	



Hose connectors



Hose connections ROTILABO® system Keck

ROTH SELECTION. Material: PBT, aeals made of EPDM.

For quickly establishing secure and leak-tight connections between hoses and glass nozzles/glass threads. Excellent chemical and temperature resistance up to +130 °C.

Not a medical device / Not an IVD product

For glass threads

Delivery incl. screw cap, hose part and seals.

Ø ext. (mm)	Suitable for	Art. No.	Pack Qty.
9,0	GL 14 (Schott), yellow	H075.1	10 unit(s)

For glass nozzles

Delivery incl. clamping ring, screw cap, hose part and seals.

Ø ext. (mm)	Suitable for	Art. No.	Pack Qty.
9,0	GL 14 (Schott), yellow	H072.1	10 unit(s)



A 121 °C

Hose couplings ROTILABO® system Keck

ROTH SELECTION. Autoclavable: Yes.

Not a medical device / Not an IVD product

Hose nozzles

Made of PBT, seal made of EPDM. Outer Ø 9 mm, length 43 mm. Light blue colour.

Delivery without connecting pieces, please order separately.

Version	Art. No.	Pack Qty.
Hose nozzles	0354.1	10 unit(s)

Connecting pieces

White colour, 8 mm outlet.

Delivery without hose nozzles, please order separately.

Version	Art. No.	Pack Qty.
Connecting pieces	0343.1	5 unit(s)



Draincock red/natural

ROTH SELECTION. Material: PE.

With separate ventilation channel. Continuously variable regulation of the discharge speed.

Not suitable for glass threads.

Thread	Art. No.	Pack Qty.
DIN 45	X530.1	3 unit(s)
DIN 50	X531.1	3 unit(s)



Draincock FLUXX®

ROTH SELECTION. Material: PE.

The patented ventilation system allows easy, clean and fast filling and transfer of liquids without gurgling.

Not suitable for glass threads.

Thread	Art. No.	Pack Qty.
DIN 51 / DIN 40	T358.1	3 unit(s)
DIN 61 / DIN 40	T359.1	3 unit(s)



Spout

ROTH SELECTION. Material: PE, Alkozell seal.

Flexible spout tube with snap-on lid, screwed into a screw cap. With riser for smooth pouring.

Not suitable for glass threads

Thread	Art. No.	Pack Qty.
DIN 45	TE90.1	2 unit(s)
DIN 51	TE91.1	2 unit(s)
DIN 61	TE92.1	2 unit(s)

Draincocks



Draincock Compact

Bürkle. Material: PP, seals made of PE/PTFE.

Closes tight. Easy to open and close. Can be operated single-handed. Space-saving design: The tap does not protrude beyond the base of the container and can remain on the container in storage. Optimised drainage spout for an even, splash-free flow.

Thread	Art. No.	Pack Qty.	
G $\frac{3}{4}$ " (internal thread)	PE42.1	1 unit(s)	



Draincock

Bürkle.

Closes tight. Easy to open and close. Can be operated single-handed. With rotating non-drip nozzle and metering scale. With reducing outlet spout (outside Ø 6–8 mm) for narrow bottle openings.

Non-self-ventilating draincock.

Stopcock

Material: PP Seal made of PE/PTFE.

Delivery incl. reducing outlet spout.

Thread	Conductive	Art. No.	Pack Qty.	
G $\frac{3}{4}$ " (internal thread)	No	2683.1	1 unit(s)	

PA stopcock

Material: PA electrically conductive, Seal made of PTFE.

Note: The drum and the tap must be grounded using the Roth antistatic set (Art. No. 0213.1) during filling with flammable liquids.

Delivery incl. reducing outlet spout and electrically conductive PA adapter for drum with G $\frac{3}{4}$ " internal thread.

Thread	Conductive	Art. No.	Pack Qty.	
G $\frac{3}{4}$ " (internal thread)	Yes	LP75.1	1 unit(s)	



Draincock Air

Bürkle. Material: PE, seal made of natural rubber.

Self-ventilating.

Thread	Art. No.	Pack Qty.	
G $\frac{3}{4}$ " external thread	Y496.1	3 unit(s)	



Draincock PP

Bürkle. Material: PP.

Ball valve suitable for crystallising media. Pressure max. 4 bar at 20 °C.

Delivery without hose connector, please order separately.

Thread	Art. No.	Pack Qty.	
G $\frac{3}{4}$ " external thread	P943.1	1 unit(s)	

Accessories hose connector for PP draincock

Bürkle.

Version	Art. No.	Pack Qty.	
Bent 90°	P946.1	1 unit(s)	
straight	P947.1	1 unit(s)	

Алматы (7273)495-231	Иваново (4932)77-34-06	Магнитогорск (3519)55-03-13	Ростов-на-Дону (863)308-18-15	Тольятти (8482)63-91-07
Ангарск (3955)60-70-56	Ижевск (3412)26-03-58	Москва (495)268-04-70	Рязань (4912)46-61-64	Томск (3822)98-41-53
Архангельск (8182)63-90-72	Иркутск (395)279-98-46	Мурманск (8152)59-64-93	Самара (846)206-03-16	Тула (4872)74-02-29
Астрахань (8512)99-46-04	Казань (843)206-01-48	Набережные Челны (8552)20-53-41	Саранск (8342)22-96-24	Тюмень (3452)66-21-18
Барнаул (3852)73-04-60	Калининград (4012)72-03-81	Нижний Новгород (831)429-08-12	Санкт-Петербург (812)309-46-40	Ульяновск (8422)24-23-59
Белгород (4722)40-23-64	Калуга (4842)92-23-67	Новокузнецк (3843)20-46-81	Саратов (845)249-38-78	Улан-Удэ (3012)59-97-51
Благовещенск (4162)22-76-07	Кемерово (3842)65-04-62	Ноябрьск (3496)41-32-12	Севастополь (8692)22-31-93	Уфа (347)229-48-12
Брянск (4832)59-03-52	Киров (8332)68-02-04	Новосибирск (383)227-86-73	Симферополь (3652)67-13-56	Хабаровск (4212)92-98-04
Владивосток (423)249-28-31	Коломна (4966)23-41-49	Омск (3812)21-46-40	Ставрополь (8652)20-65-13	Чебоксары (8352)28-53-07
Владикавказ (8672)28-90-48	Кострома (4942)77-07-48	Орел (4862)44-53-42	Сочи (862)225-72-31	Челябинск (351)202-03-61
Владимир (4922)49-43-18	Краснодар (861)203-40-90	Оренбург (3532)37-68-04	Сыктывкар (8212)25-95-17	Череповец (8202)49-02-64
Волгоград (844)278-03-48	Красноярск (391)204-63-61	Пенза (8412)22-31-16	Тамбов (4752)50-40-97	Чита (3022)38-34-83
Вологда (8172)26-41-59	Курск (4712)77-13-04	Петрозаводск (8142)55-98-37	Сургут (3462)77-98-35	Якутск (4112)23-90-97
Воронеж (473)204-51-73	Курган (3522)50-90-47	Псков (8112)59-10-37	Тверь (4822)63-31-35	Ярославль (4852)69-52-93
Екатеринбург (343)384-55-89	Липецк (4742)52-20-81	Пермь (342)205-81-47		
Россия (495)268-04-70	Казахстан (772)734-952-31	Киргизия (996)312-96-26-47		