

Инкубаторы CARL ROTH DipSlide EM4, Cultura M, Heratherm IMC18, TIN-IN16B, TIN-IN35B, BD-S Solid.Line, KT 53, KT 115

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47



Incubators



DipSlide incubator EM4 EM4

Laboratories Humeau.

Lightweight, small incubator specially designed for DipSlides. With transparent, easily removable cover for quick insertion and removal of DipSlides. Easy to operate, even for inexperienced users. A stand insert and a thermometer are also supplied. Also usable without a stand for incubation of Petri dishes and contact plates.

Technical specifications:

Art. No.	KK43.1
Type	DipSlide incubator EM4
Suitable for	23 DipSlides
Operating temperature	Room temperature up to +44 °C
Vol.	4 l
Inner width	180 mm
Inner depth	180 mm
Inner height	125 mm
Outer width	230 mm
Outer depth	230 mm
Outer height	180 mm
Weight	1,5 kg
Mains connection	200 V

Type	Suitable for	Vol. (l)	Art. No.	Pack Qty.	
DipSlide incubator EM4	23 DipSlides	4	KK43.1	1 unit(s)	



Small incubator Cultura® M

For incubating 18 dip slides (e.g. ROTI® DipSlides culture medium) or 12 Petri dishes and other small vessels for verifying bacteria, fungi and yeasts.


- **Operating temperature range from +5 °C above room temperature to 45 °C; temperature adjustment range from 25 to 45 °C**
- Thermometer in tube with special petroleum filling. Measuring range 7 to 60 °C. Graduations 1 °C. Accuracy ±1 °C. Dimensions: H 95 x Ø 25 mm.
- Transparent full-vision door for visual inspection
- Minimal footprint, i.e. space-saving

Technical specifications:

Art. No.	AHH9.1
Type	Cultura® M
Operating temperature	From +5 °C above room temperature to +45 °C
Set temperature from	+25 °C
Set temperature up to	+45 °C
Vol.	4 l
Inner width	220 mm
Inner depth	150 mm
Inner height	120 mm
Outer width	310 mm
Outer depth	168 mm
Outer height	155 mm
Temperature variance	±1 K
Weight	1.1 kg
Mains connection	230 V, 50 Hz

Delivery incl. control thermometer in the tube and shelf.

Type	Vol. (l)	Art. No.	Pack Qty.	
Cultura® M	4	AHH9.1	1 unit(s)	



Multirack for Cultura® M small incubator

Suitable for Cultura® M mini incubator. Made of ABS. For 18 tubes (Ø 25 to 30 mm). Dimensions: W 213 mm x D 142 mm x H 50 mm.

Designation	Art. No.	Pack Qty.	
Multirack	AHK0.1	1 unit(s)	

Incubators



Small incubator

Labnet International.

Space-saving unit for incubating culture bottles, petri dishes, micro-titration plates and other small-volume vessels.

- **Operating temperature range from + 5 °C above ambient temperature to 60 °C; temperature adjustment range from 20 to 60 °C**
- Cable conduit (Ø 55 mm) at the rear
- Loading options, e.g. 48 petri dishes or 28 microtitration plates in 96 well format
- Thermometer holder (Ø 7 mm) on the top allows the inner temperature to be monitored

Technical specifications:

Art. No.	YC53.1
Type	Small incubator
Operating temperature	From +5 °C above room temperature to +60 °C
Set temperature from	+20 °C
Set temperature up to	+60 °C
Vol.	9.2 l
Inner width	230 mm
Inner depth	200 mm
Inner height	200 mm
Outer width	285 mm
Outer depth	280 mm
Outer height	335 mm
Max. removable shelves	3
Temperature variance	±0,7 K
Weight	8.3 kg
Mains connection	230 V, 50/60 Hz

Delivery incl. 1 shelf (L 220 x W 180 mm). Please order thermometer separately.

Type	Vol. (l)	Art. No.	Pack Qty.
Small incubator	9.2	YC53.1	1 unit(s)



Small incubator Heratherm™ IMC18

Thermo.

- Operating temperature range from 17 to 40 °C
- With Peltier cooling
- Minimal footprint for space-restricted lab areas
- Microprocessor control
- Round edges make it easy to clean
- Interior lighting facilitates sample monitoring
- Fan for air circulation; runs permanently and is not adjustable

Technical specifications:

Art. No.	HAP4.1	PL74.1
Type	Darkened front panel	Standard model – glass front
Operating temperature	Operating temperature range from 17 to 40 °C	
Vol.	18 l	
Resolution	1 °C	
Inner width	180 mm	
Inner depth	310 mm	
Inner height	290 mm	
Outer width	260 mm	
Outer depth	470 mm	
Outer height	415 mm	
Max. removable shelves	3	
Temporal temp. variance	±0,2 K	
Spatial temp. deviation	±1,2 K	
Weight	7.2 kg	
Mains connection	100–240 V, 50/60 Hz	

Delivery incl. two removable shelves.

With glass front – standard model

Type	Vol. (l)	Art. No.	Pack Qty.
Standard model – glass front	18	PL74.1	1 unit(s)

With darkened front panel – IMC18 bio model

For incubation without the influence of light, cool long-term storage in the dark, light-sensitive micro-organisms and much more.

Type	Vol. (l)	Art. No.	Pack Qty.
Darkened front panel	18	HAP4.1	1 unit(s)

Accessories additional tray for IMC18

Thermo.

Designation	Art. No.	Pack Qty.
Additional shelf for IMC18	PL75.1	1 unit(s)



Incubators



Small incubator TIN series

Phönix Instrument.

- Operating temperature range from +5 °C above room temperature to 70 °C; temperature adjustment range from 20 to 70 °C
- Digital timer (adjustable from 1 min to 99 hrs or continuous run)
- Glass window in the door for monitoring the samples
- Natural air flow (convection)
- Stainless steel interior
- Test opening (Ø 5 mm) on top, e.g. for introducing sensors

Technical specifications:

Art. No.	EEH0.1	EEH1.1
Type	TIN-IN16B	TIN-IN35B
Operating temperature	From +5 °C above room temperature to +70 °C	
Set temperature from	+20 °C	
Set temperature up to	+70 °C	
Vol.	16 l	35 l
Inner width	270 mm	360 mm
Inner depth	230 mm	300 mm
Inner height	255 mm	320 mm
Outer width	505 mm	595 mm
Outer depth	370 mm	440 mm
Outer height	400 mm	460 mm
Max. removable shelves	6	
Rated power	85 W	125 W
Warm-up period	18 min (37 °C)	22 min (37 °C)
Temperature variance	±0,4 K	
Weight	23 kg	33 kg
Mains connection	230 V, 50/60 Hz	

Delivery incl. 2 removable shelves.

Type	Vol. (l)	Art. No.	Pack Qty.
TIN-IN16B	16	EEH0.1	1 unit(s)
TIN-IN35B	35	EEH1.1	1 unit(s)

Accessories grate made of stainless steel for TIN series incubators

Phönix Instrument.

Suitable for	Art. No.	Pack Qty.
for TIN-IN16B	EEH2.1	1 unit(s)
for TIN-IN35B	EEH3.1	1 unit(s)



Incubator BD-S Solid.Line

Binder.

- Operating temperature range from +5 °C above room temperature to 70 °C
- APT.line™ pre-heating chamber technology
- Natural convection
- Adjustable exhaust air flap
- Controller with timer function (1 min to 99 h)
- Inner door made of safety glass
- Integrated Class 3.1 temperature safety device (DIN 12880) with visual alarm

Technical specifications:

Art. No.	KCT5.1	KCT6.1
Type	BD-S 056	BD-S 115
Operating temperature	From +5 °C above room temperature to 70 °C	
Vol.	62 l	118 l
Resolution	1 °C	
Inner width	400 mm	550 mm
Inner depth	350 mm	390 mm
Inner height	440 mm	550 mm
Outer width	614 mm	764 mm
Outer depth	580 mm	662 mm
Outer height	625 mm	735 mm
Max. removable shelves	4	5
Temporal temp. variance	±0,1 °C	±0,2 °C
Spatial temp. deviation	±0,7 °C	±0,5 °C
Weight	41 kg	57 kg
Mains connection	230 V, 50 Hz	

Delivery incl. two chrome-plated insertable grates.

Type	Vol. (l)	Art. No.	Pack Qty.
BD-S 056	62	KCT5.1	1 unit(s)
BD-S 115	118	KCT6.1	1 unit(s)

Accessories racks, chrome-plated for Solid.Line incubators and drying ovens

Binder.

Suitable for	Art. No.	Pack Qty.
for BD-S 056/ED-S 056/FD-S 056	KCT7.1	1 unit(s)
for BD-S 115/ED-S 115/FD-S 115	KCT8.1	1 unit(s)

Incubators



Incubator Heratherm™ Protocol series

Thermo.

General features:

- Housing interior with round edges make it easy to clean
- With additional inner glass door
- RS 232 port
- Natural convection minimises desiccation of samples; air circulation achieves more accurate temperatures
- Handle and display (65 mm) not included in depth; adjustable feet (35 mm) not included in height
- Mains connection 230 V, 50 Hz

Technical specifications:

Art. No.	PL96.1	PL97.1	PL98.1	PL99.1	PP00.1	PP01.1
Type	IGS60	IGS100	IGS180	IMH60	IMH100	IMH180
Operating temperature	From +5 °C above room temperature to +75 °C			From +5 °C above room temperature to +105 °C		
Vol.	75 l	117 l	194 l	66 l	104 l	178 l
Inner width	354 mm	464 mm		354 mm	464 mm	
Inner depth	414 mm		589 mm	368 mm		543 mm
Inner height	508 mm	608 mm	708 mm	508 mm	608 mm	708 mm
Outer width	530 mm	640 mm		530 mm	640 mm	
Outer depth	565 mm		738 mm	565 mm		738 mm
Outer height	720 mm	820 mm	920 mm	720 mm	820 mm	920 mm
Max. removable shelves	13	16	19	13	16	19
Rated power	300 W	520 W	710 W	850 W	1,100 W	1,300 W
Temporal temp. variance	±0,2 °C			±0,1 °C		
Spatial temp. deviation	±0,6 °C			±0,2 °C	±0,3 °C	±0,4 °C
Weight	40 kg	51 kg	65 kg	45 kg	56 kg	70 kg
Mains connection	230 V, 50 Hz					

Delivery incl. two removable shelves. Please note: Additional transport costs apply for delivery to the place of use. We will be happy to inform you of these on request.

General Protocol with natural convection

- **Operating temperature range: From +5 °C above room temperature to 75 °C**
- Stainless steel interior (1.4016)
- Integrated timer (1 min to 99 h, 59 min)
- Microprocessor control
- Actual temperature display
- Setting accuracy ±0.1 K
- Automatic overheating alarm

Type	Vol. (l)	Art. No.	Pack Qty.
IGS60	75	PL96.1	1 unit(s)
IGS100	117	PL97.1	1 unit(s)
IGS180	194	PL98.1	1 unit(s)

Advanced Protocol with air-circulation fan

- **Operating temperature range: From +5 °C above room temperature to 105 °C**
- Stainless steel interior (1.4301)
- Digital timer for daily/weekly on/off cycles
- Setting accuracy ±0.1 K
- Automatic overheating alarm
- Cable conduit (Ø 50 mm) at the rear
- With varnished housing
- With air-circulation fan: Dual convection with adjustable fan speed from 0 (corresponds to natural convection) to 100%

Type	Vol. (l)	Art. No.	Pack Qty.
IMH60	66	PL99.1	1 unit(s)
IMH100	104	PP00.1	1 unit(s)
IMH180	178	PP01.1	1 unit(s)

Accessories shelf kit for Heratherm® incubators

Thermo.

Shelf kit: One shelf with two mounting clips.

Designation	Art. No.	Pack Qty.
for IGS60	PP02.1	1 unit(s)
for IGS100	PP03.1	1 unit(s)
for IGS180	PP04.1	1 unit(s)
for IMH60	PP05.1	1 unit(s)
for IMH100	PP06.1	1 unit(s)
Insert kit for IMH180 and IMP180	PP07.1	1 unit(s)



Incubators

Incubator IN/INplus series

Memmert.

General features:

- **Operating temperature range from 5 °C above room temperature to 80 °C**
Set temperatures from 20 to 80 °C
Adjustment accuracy 0.1 K
- Stainless steel interior (1.4301)
- Fully insulated stainless steel door (two-wing from size 450 upwards)
- **With additional inner glass door**
- Fresh air supply by electronically adjustable exhaust air flap
- Textured stainless steel housing
- Model sizes 450 and 750, with moveable and lockable rollers
- Depth plus door handle: 56 mm
- Incl. factory calibration certificate valid for +37 °C



Other features:

- **With natural air circulation (convection)**
- German, English, French and Spanish language options
- Digital timer (adjustable from 1 min to 99 d, 23 h) or continuous operation
- "SetpointWAIT" function – process time does not start until the target temperature is reached

Delivery incl. stainless steel grating (models 30 and 55: One grate, models 75 to 750: Two grates). Please note: Additional transport costs apply for delivery to the place of use. We will be happy to inform you of these on request.

IN series models

- Graphical single TFT display with touch screen for operation and displaying temperature, program runtime and air flap position
- Integrated data logger with 10-year memory
- With Ethernet communication port for reading out log data
- Protection against overheating

Technical specifications:

Art. No.	TY30.1	TY31.1	TY32.1	TY33.1	TY34.1	TY35.1	TY36.1	TY37.1
Type	IN 30	IN 55	IN 75	IN 110	IN 160	IN 260	IN 450	IN 750
Operating temperature	From +5 °C above room temperature to +80 °C							
Vol.	32 l	53 l	74 l	108 l	161 l	256 l	449 l	749 l
Inner width	400 mm			560 mm		640 mm	1,040 mm	
Inner depth	250 mm	330 mm		400 mm		500 mm	600 mm	
Inner height	320 mm	400 mm	560 mm	480 mm	720 mm	800 mm	720 mm	1,200 mm
Outer width	585 mm			745 mm		824 mm	1,224 mm	
Outer depth	434 mm	514 mm		584 mm		684 mm	784 mm	
Outer height	707 mm	787 mm	947 mm	867 mm	1,107 mm	1,186 mm	1,247 mm	1,726 mm
Max. removable shelves	3	4	6	5	8	9	8	14
Rated power	800 W	1,000 W	1,250 W	1,400 W	1,600 W	1,700 W	1,800 W	2,000 W
Temporal temp. variance	<±0,2 °C							
Spatial temp. deviation	<±0,7 °C							
Weight	44 kg	55 kg	65 kg	74 kg	96 kg	110 kg	170 kg	215 kg
Mains connection	230 V, 50/60 Hz							

Type	Vol. (l)	Art. No.	Pack Qty.
IN 30	32	TY30.1	1 unit(s)
IN 55	53	TY31.1	1 unit(s)
IN 75	74	TY32.1	1 unit(s)
IN 110	108	TY33.1	1 unit(s)
IN 160	161	TY34.1	1 unit(s)
IN 260	256	TY35.1	1 unit(s)
IN 450	449	TY36.1	1 unit(s)
IN 750	749	TY37.1	1 unit(s)

Incubators

BF series with forced air circulation (fan)

- Fan speed cannot be regulated or switched off
- Timer function selection:
 - Heater is switched off when the set time expires
 - The set timer does not start until the actual value reaches the selected set value. The heater is switched off when the time expires.
 - Heater is switched on after set time expires and stays in continuous operation.
- Models BF 56, BF 115 and BF 260 also have a USB port for recording data (time, actual temperature value, set temperature value, etc.)

Technical specifications:

Art. No.	CKH2.1	CKH3.1	CKH4.1	CKH5.1	CKH6.1
Type	BF 56	BF 115	BF 260	BF 400	BF 720
Operating temperature	From +7 °C above room temperature to +100 °C				
Set temperature from	+30 °C				
Set temperature up to	+100 °C				
Vol.	59 l	114 l	257 l	400 l	734 l
Inner width	400 mm	550 mm	650 mm	1,000 mm	
Inner depth	340 mm	380 mm	510 mm	500 mm	565 mm
Inner height	440 mm	550 mm	780 mm	800 mm	1,300 mm
Outer width	560 mm	710 mm	810 mm	1,235 mm	1,165 mm
Outer depth	640 mm	680 mm	815 mm	875 mm	910 mm
Outer height	625 mm	735 mm	965 mm	1,025 mm	1,590 mm
Max. removable shelves	4	5	7	10	15
Rated power	400 W		900 W	850 W	1,250 W
Temporal temp. variance	±0,1 °C				
Spatial temp. deviation	±0,2 °C	±0,3 °C	±0,2 °C	±0,4 °C	
Weight	38 kg	54 kg	85 kg	145 kg	165 kg
Mains connection	230 V, 50/60 Hz				

Type	Vol. (l)	Art. No.	Pack Qty.
BF 56	59	CKH2.1	1 unit(s)
BF 115	114	CKH3.1	1 unit(s)
BF 260	257	CKH4.1	1 unit(s)
BF 400	400	CKH5.1	1 unit(s)
BF 720	734	CKH6.1	1 unit(s)

Accessories grate for incubators and drying ovens for the BD, BF, ED, FD and FED series

Binder.

Grids, chrome plated

Suitable for	Art. No.	Pack Qty.
for BD 23, ED 23 and FD 23	CKE0.1	1 unit(s)
for BD 56 and ED 56	CKE1.1	1 unit(s)
for BF 56, FD 56 and FED 56	ECY1.1	1 unit(s)
for BD 115 and ED 115	CKE2.1	1 unit(s)
for BF 115, FD 115 and FED 115	ECY2.1	1 unit(s)
for BD 260 and ED 260	CKE3.1	1 unit(s)
for BF 260, FD 260 and FED 260	ECY3.1	1 unit(s)
for BD 400, BF 400, ED 400 and FED 400	CKE4.1	1 unit(s)
for BD 720 and ED 720	CKE5.1	1 unit(s)
for BF 720 and FED 720	HPA7.1	1 unit(s)



Incubators

High-capacity incubators Heratherm™

Thermo.

General features:

- Free-standing devices with fixed rollers
- With internal glass door
- Stainless steel interior
- RS 232 port
- Not including depth of handle and display (65 mm) and height of rollers (108 mm)
- Mains connection 230 V, 50 Hz



Technical specifications:

Art. No.	XH26.1	XH27.1	XH28.1	XH29.1
Type	IGS400	IGS750	IMH400-S	IMH750-S
Operating temperature	From +5 °C above room temperature to +75 °C		From +5 °C above room temperature to +105 °C	
Vol.	405 l	747 l	381 l	702 l
Inner width	544 mm	1,004 mm	544 mm	1,004 mm
Inner depth	569 mm		524 mm	
Inner height	1,307 mm		1,335 mm	
Outer width	778 mm	1,261 mm	778 mm	1,261 mm
Outer depth	770 mm			
Outer height	1,545 mm			
Max. removable shelves	39			
Rated power	1,200 W	1,500 W	1,400 W	1,600 W
Temporal temp. variance	±0,4 °C		±0,2 °C	
Spatial temp. deviation	±0,5 °C	±1,3 °C	±0,2 °C	±0,3 °C
Weight	145 kg	201 kg	144 kg	205 kg
Mains connection	230 V, 50 Hz			

Delivery incl. two stainless steel perforated sheets. Please note: Additional transport costs apply for delivery to the place of use. We will be happy to inform you of these on request.

General Protocol

- **Operating temperature range from +5 °C above room temperature to 75 °C**
- With natural convection
- Integrated timer (1 min to 99 h, 59 min)
- Actual temperature display; setting accuracy ±0.1 K

Type	Vol. (l)	Art. No.	Pack Qty.
IGS400	405	XH26.1	1 unit(s)
IGS750	747	XH27.1	1 unit(s)

Advanced Protocol

- **Operating temperature range from +5 °C above room temperature to 105 °C**
- Setting accuracy ±0.1 K
- Mechanical convection – two-stage fan: Can run fast or slow, but cannot be switched off
- Cable conduit (Ø 19 mm), e.g. for fitting sensors
- An audible alarm sounds if the door is accidentally left open
- Switch on/off after a preset time
- Weekly time setting or real-time control
- Programmable control for temperature ramps and dwell times (memory for 10 programs each with 10 individual steps)

Type	Vol. (l)	Art. No.	Pack Qty.
IMH400-S	381	XH28.1	1 unit(s)
IMH750-S	702	XH29.1	1 unit(s)

Accessories perforated sheet kit for Heratherm® high-capacity incubators

Thermo.

Perforated shelf kit: One shelf with two supports.

Designation	Art. No.	Pack Qty.
Insert kit for IGS400, IMH400-S and IMP400	XH30.1	1 unit(s)
Insert kit for IGS750, IMH750-S	XH31.1	1 unit(s)

Incubators



Incubator inserts

Made of PP, stackable. For easy organisation of Petri dishes and convenient transfer to other incubation or storage systems.
With 16 ventilation holes and two labelling areas.
Outer dimensions: W 251 mm x D 237 mm x H 35 mm.
Loading options: E.g. Five Petri dishes Ø 100 mm
or 15 Petri dishes Ø 60 mm or 21 culture bottles of 25 cm².



Peltier cooling incubator KT series

Binder.

General features:

- **Operating temperature range: From 4 to 100 °C**
Setting accuracy: ± 0.1 °C
- 5.7-inch LCD monitor and turn/push button for input
- Variable time and weekday programming
 - Storage capacity for up to 60 programs (for KT 53) or 100 programs (for KT 115)
 - Measurement data can be saved to internal memory
 - Readout of programs and data via USB port
 - Adjustable fan speed (min. speed 30%)
- With additional inner door of safety glass
- Ethernet interface
- Class 3.1 temperature safety device (DIN 12880) with visual and audible alarms

Technical specifications:

Art. No.	CKK0.1	CKK1.1
Type	KT 53	KT 115
Vol.	53 l	102 l
Inner width	400 mm	600 mm
Inner depth	334 mm	355 mm
Inner height	400 mm	455 mm
Outer width	660 mm	860 mm
Outer depth	695 mm	815 mm
Outer height	635 mm	715 mm
Max. removable shelves	5	6
Rated power	400 W	700 W
Temporal temp. variance	± 0.1 K	
Spatial temp. deviation	± 0.3 K	
Weight	63 kg	83 kg
Mains connection	200–240 V, 50/60 Hz	

Delivery incl. two stainless steel grids.

Type	Vol. (l)	Art. No.	Pack Qty.	
KT 53	53	CKK0.1	1 unit(s)	
KT 115	102	CKK1.1	1 unit(s)	

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Россия (495)268-04-70

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Казахстан (772)734-952-31

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Киргизия (996)312-96-26-47

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93