

Стаканы лабораторные мерные, колбы Эрленмейера CARL ROTH ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47



Measuring beakers



Griffin beakers

Vitlab. Material: PP, transparent. Autoclavable: Limited.

With spout and blue scale. Cleaning at max. 60 °C recommended to protect the printing. Suitable for foodstuffs.

Vol. (ml)	Gradua- tions (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
10	2	30	36	2870.1	1 unit(s)	
25	5	38	50	2871.1	1 unit(s)	
50	10	47	60	2872.1	1 unit(s)	
100	10	55	70	2873.1	1 unit(s)	
150	20	66	80	2874.1	1 unit(s)	
250	25	77	95	2875.1	1 unit(s)	
400	50	87	112	2876.1	1 unit(s)	
500	50	94	118	HX26.1	1 unit(s)	
600	50	100	127	2877.1	1 unit(s)	
1,000	100	120	147	2879.1	1 unit(s)	
2,000	200	149	187	P804.1	1 unit(s)	
3,000	250	170	212	P805.1	1 unit(s)	
5,000	500	203	247	P806.1	1 unit(s)	

Griffin beakers assortment

Vitlab. Material: PP.

10 pieces (one beaker each from 10 to 1000 ml).

Version	Art. No.	Pack Qty.	
Griffin beakers assortment	2880.1	1 set	



Griffin beakers

Vitlab. Material: PFA, transparent. Autoclavable: Yes.

With raised scale. Highly resistant to chemicals. Temperature-stable from -200 to +250 °C.

Vol. (ml)	Gradua- tions (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
25	5	32	50	A862.1	1 unit(s)	
50	10	39	59	A863.1	1 unit(s)	
100	20	50	72	A864.1	1 unit(s)	
250	50	67	96	A865.1	1 unit(s)	
500	100	88	122	A866.1	1 unit(s)	
1,000	100	105	143	A867.1	1 unit(s)	



Griffin beakers

Vitlab. Material: PMP, transparent. Autoclavable: Limited.

With spout and printed red scale. Cleaning at max. 60 °C, recommended to protect the printing. Good chemical resistance.

Vol. (ml)	Gradua- tions (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
10	2	30	36	2890.1	1 unit(s)	
25	5	38	50	2891.1	1 unit(s)	
50	10	47	60	2892.1	1 unit(s)	
100	10	55	70	2893.1	1 unit(s)	
150	20	66	80	2894.1	1 unit(s)	
250	25	77	95	2895.1	1 unit(s)	
400	50	87	112	2896.1	1 unit(s)	
500	50	94	118	HX27.1	1 unit(s)	
600	50	100	127	2897.1	1 unit(s)	
1,000	100	120	147	2899.1	1 unit(s)	
2,000	200	149	187	P807.1	1 unit(s)	
3,000	250	170	212	P808.1	1 unit(s)	
5,000	500	203	247	P809.1	1 unit(s)	

Griffin beakers assortment

Vitlab. Material: PMP.

10 pieces (one beaker each from 10 to 1000 ml).

Version	Art. No.	Pack Qty.	
Griffin beakers assortment	2900.1	1 set	



Griffin beakers

BRAND. Material: ETFE, transparent. Autoclavable: Limited.

With spout and printed black scale. Cleaning at max. 60 °C recommended to protect the printing. Very good chemical resistance. Temperature-stable from -100 to +150 °C.

Vol. (ml)	Gradua- tions (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
25	5	32	50	2881.1	1 unit(s)	
50	10	39	59	2882.1	1 unit(s)	
100	20	50	72	2883.1	1 unit(s)	
250	50	67	96	2884.1	1 unit(s)	
400	50	77	109	2885.1	1 unit(s)	
500	100	88	122	HX28.1	1 unit(s)	
600	100	91	125	2886.1	1 unit(s)	
1,000	100	105	143	2887.1	1 unit(s)	

Measuring beakers



A 121 °C

Griffin beakers

Vitlab. Material: PP, transparent. Autoclavable: Yes.

With spout and embossed scale. Temperature-stable from 0 to +125 °C. Suitable for foodstuffs.

Vol. (ml)	Gradua- tions (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
25	5	38	50	CKK4.1	1 unit(s)	
50	10	47	60	CKK5.1	1 unit(s)	
100	10	55	70	CKK6.1	1 unit(s)	
150	20	66	80	CKK7.1	1 unit(s)	
250	25	77	95	CKK8.1	1 unit(s)	
400	50	87	111.5	CKK9.1	1 unit(s)	
500	50	94	117.5	CKL0.1	1 unit(s)	
600	50	100	126.5	CKL1.1	1 unit(s)	
1,000	100	120	147	CKL2.1	1 unit(s)	
2,000	200	149	187	CKL3.1	1 unit(s)	
3,000	250	170	212	CKL4.1	1 unit(s)	
5,000	500	203	247	CKL5.1	1 unit(s)	



Measuring beakers ROTILABO®

ROTH SELECTION. Material: PP, Transparent. Autoclavable: Limited.

With smooth surface and raised, red embossed scale. Cleaning at max. 60 °C recommended to protect the printing. Suitable for foodstuffs.

Vol. (ml)	Gradua- tions (ml)	Ø Upper (mm)	Height (mm)	Art. No.	Pack Qty.	
50	2	40	70	YA29.1	1 unit(s)	
100	2	50	80	YA30.1	1 unit(s)	
250	5	70	120	YA31.1	1 unit(s)	
500	10	91	133	0757.1	1 unit(s)	
1,000	10	116	170	0758.1	1 unit(s)	
2,000	20	150	215	0768.1	1 unit(s)	
3,000	50	170	242	0779.1	1 unit(s)	
5,000	100	210	270	0780.1	1 unit(s)	



A 121 °C

Measuring beakers

Vitlab. Material: PP, Highly transparent. Autoclavable: Yes.

With embossed scale and sturdy handle with good grip. Temperature-stable from 0 to +125 °C. Suitable for foodstuffs.

Vol. (ml)	Gradua- tions (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
50	2	40	70	CKL6.1	1 unit(s)	
100	2	50	80	CKL7.1	1 unit(s)	
250	5	70	120	CKL8.1	1 unit(s)	
500	10	91	133	CKL9.1	1 unit(s)	
1,000	20	116	170	CKN0.1	1 unit(s)	
2,000	50	150	215	CKN1.1	1 unit(s)	
3,000	50	170	242	CKN2.1	1 unit(s)	
5,000	100	210	270	CKN3.1	1 unit(s)	



Measuring beakers

ROTH SELECTION. Material: PP, Highly transparent.

Autoclavable: Limited.

With black scale printed on two sides. Stackable due to open handle. Cleaning at max. 60 °C recommended to protect the printing.

Vol. (ml)	Gradua- tions (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
250	5	75	115	XY21.1	1 unit(s)	
500	10	100	140	P800.1	1 unit(s)	
1,000	10	125	167	P801.1	1 unit(s)	
2,000	20	148	212	P802.1	1 unit(s)	
3,000	50	170	242	P803.1	1 unit(s)	



Measuring beakers



Measuring beakers

Vitlab. Material: PP, transparent.

With clearly legible scale printed on both sides and stable handle. Stackable. To preserve the imprint, cleaning at temperatures not exceeding 60 °C is recommended. Suitable for foodstuffs.

Vol. (ml)	Colour	Graduations (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.
500	blue	10	100	140	APH6.1	1 unit(s)
500	yellow	10	100	140	APH7.1	1 unit(s)
500	red	10	100	140	APH8.1	1 unit(s)
500	green	10	100	140	APH9.1	1 unit(s)



Measuring beakers

Vitlab. Material: SAN, transparent.

With embossed scale and sturdy handle with good grip. Temperature-stable between -40 and +70 °C. Suitable for foodstuffs.

Vol. (ml)	Graduations (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.
250	5	70	120	CKN4.1	1 unit(s)
500	10	91	133	CKN5.1	1 unit(s)
1,000	10	116	170	CKN6.1	1 unit(s)
2,000	20	150	215	CKN7.1	1 unit(s)
3,000	50	170	242	CKN8.1	1 unit(s)



A 121 °C



Measuring jugs

Hünersdorff. Material: PP, transparent. Autoclavable: Yes.

With injected, transparent measuring scale and additional blue stamped scale. Closed handle. Temperature-stable between -20 and +100 °C.

Delivery without lid, please order separately.

Vol. (ml)	Graduations (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.
1,000	10	117	165	EKE7.1	1 unit(s)
2,000	20	151	206	EKE8.1	1 unit(s)

Accessories lid for measuring jugs

Hünersdorff. Material: PP.

Designation	Art. No.	Pack Qty.
Lid for 1000-ml jug EKE7.1	EKE9.1	1 unit(s)
Lid for 2000-ml jug EKE8.1	EKH0.1	1 unit(s)



A 121 °C

Measuring beakers ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10, polished.

Autoclavable: Yes.

Graduated inside. With open handle, stackable. Suitable for foodstuffs.

Vol. (ml)	Graduations (ml)	Height (mm)	Art. No.	Pack Qty.
100	50	50	YY00.1	1 unit(s)
250	50	70	YY01.1	1 unit(s)
500	100	110	YY02.1	1 unit(s)
1,000	500	135	YY03.1	1 unit(s)
2,000	500	195	YY04.1	1 unit(s)

Measuring beakers



A 121 °C

Measuring beakers ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes. Cylindrical. Graduated inside and outside.

Vol. (ml)	Gradua- tions (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
700	100	100	100	9472.1	1 unit(s)	
1,300	100	120	120	9473.1	1 unit(s)	
2,100	200	140	140	9475.1	1 unit(s)	
3,200	200	160	160	9476.1	1 unit(s)	
4,000	500	180	180	9478.1	1 unit(s)	
5,000	1,000	200	200	9479.1	1 unit(s)	



A 121 °C

Measuring beakers ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes. Inside graduated. Suitable for foodstuffs.

Vol. (ml)	Gradua- tions (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
30	10	50	25	LT24.1	5 unit(s)	
50	10	53	45	LT25.1	5 unit(s)	



A 121 °C

Measuring beakers 51

Haldenwanger. Material: Porcelain, glazed (outside surface non-glazed). Autoclavable: Yes.

Graduated. Temperature stable up to +400 °C.

Vol. (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
250	80	100	XL85.1	1 unit(s)	
500	95	125	XL86.1	1 unit(s)	
1,000	110	165	XL87.1	1 unit(s)	
2,000	140	193	XL88.1	1 unit(s)	



A 121 °C

Erlenmeyer flasks disposable

Simport. Material: PP, transparent. Autoclavable: Yes.

Graduated from 0 to 170 ml. Specially shaped edge enables drip-free pouring. Can also be used for cultivating drosophila.

Delivery without snap-on lid, please order separately.

Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
55	55	105	AYP4.1	500 unit(s)	

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47