

Пипетки Пастера CARL ROTH CELLSTAR, ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47



Pasteur/disposable pipettes

Pasteur pipettes

DWK Life Sciences. Material: Borosilicate glass 5.1. Autoclavable: Yes.
For single use in microbiological and serological laboratories. Ideal for distributing aliquot quantities of liquid.

- In accordance with ISO 7712
- Practically inert borosilicate glass (USP type I neutral glass)
- With neck for cotton stoppers
- Good resistance to temperature variation, suitable for transferring hot liquids
- Excellent chemical resistance to acid and alkaline solutions



A121 °C

DURAN
WHEATON
KIMBLE

Vol. (ml)	Total length (mm)	Ø ext. (mm)	Tip length (mm)	Outer tip Ø (mm)	Pack.	Art. No.	Pack Qty.	
2	147.5	6.95	55	1.1	4 x 250	HXH7.1	1,000 unit(s)	
2	227.5	6.95	135	1.1	4 x 250	HXH8.1	1,000 unit(s)	

Pasteur pipettes plastic sterile

ROTH SELECTION. Material: PS. **Sterile**.

For use with pipette bulbs or as aspiration pipettes in combination with suction systems. Safe alternative to Pasteur pipettes made of glass.



S

Vol. (ml)	Total length (mm)	Ø ext. (mm)	Outer tip Ø (mm)	Sterility	Pack.	Art. No.	Pack Qty.	
2	145	6	1.5	Sterile	4 x 50, individual	EKN7.1	200 unit(s)	
3	230	6	1.5	Sterile	4 x 50, individual	EKN8.1	200 unit(s)	

Aspiration pipettors sterile

ROTH SELECTION. Material: PS. **Sterile**.

Pipettes without cotton plug, without graduation.

Suitable for vacuum applications and liquid handling where no high precision is required.

Safe alternative to Pasteur pipettes made of glass.



S

Vol. (ml)	Total length (mm)	Sterility	Pack.	Art. No.	Pack Qty.	
2	271	Sterile	individually	NA41.2	200 unit(s)	


Pasteur/disposable pipettes

Pasteur pipettes ungraduated

ROTH SELECTION. Material: LDPE.

Single-use pipettes with integrated suction ball for straightforward transfer of liquids where precision volume is not required. Suitable for media aliquoting thanks to reproducible droplet size.

- No slipping with latex gloves
- Precise, fine tip for reproducible results



Not a medical device / Not an IVD product

1 ml, universal

For routine work and quick liquid transfers.

Vol. (ml)	Filling capacity (ml)	Total volume (l)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
1	0.9	1.8	87	4.8	40	26	Non-sterile	1 x 500	1CC3.1	500 unit(s)	

1 ml, fine tip

Manageable pipette for applications involving limited space.

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
1	1	1.6	51	1.2	30	34	Non-sterile	1 x 500	1CC4.1	500 unit(s)	

2 ml, extra long

For use in sampling, coating of columns in chromatography or for transferring into columns.

Vol. (ml)	Filling capacity (ml)	Total volume (l)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
2	2.3	5.6	225	5	40	25	Non-sterile	1 x 400	1CC5.1	400 unit(s)	

3 ml, fine tip

For precision liquid application.

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
3	3.3	5.8	147	4.6	19	53	Non-sterile	1 x 500	1CC6.1	500 unit(s)	

3 ml, narrow stem

Narrow stem for special applications, e.g. pipetting in semi-micro cuvettes. Suitable for filling and heat sealing for the purpose of storing liquids.

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
3	3.2	6.9	155	8	32	31	Non-sterile	1 x 500	1CC7.1	500 unit(s)	

3 ml, narrow stem, extra long

Extra long, slimline body for pipetting in narrow vessels.

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
3	3.3	4.5	155	2.5	38	26	Non-sterile	1 x 500	1CC8.1	500 unit(s)	

4 ml, fine tip (extended)

Suitable for loading electrophoresis gels.

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
4	3.7	5.8	157	4.6	20	50	Non-sterile	1 x 500	1CC9.1	500 unit(s)	

6 ml, wide opening

For viscous samples, e.g. creams, oils, sludges, paints. Can be used as a funnel if the suction ball is cut open centrally and a filter is placed. Suitable for use as a spatula if the suction ball is cut at an angle.

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
6	6.2	9.5	137	6.1	50	20	Non-sterile	1 x 250	1CCA.1	250 unit(s)	



Pasteur/disposable pipettes

10 ml, extra long

For use in sampling, pipetting of media from tall vessels, coating of columns in chromatography or for transferring into columns.

Vol. (ml)	Filling capacity (ml)	Total length (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
10	10.7	300	50	20	Non-sterile	1 x 100	E251.1	100 unit(s)	

10 ml, narrow stem


For dispensing larger liquid volumes. Slimline stem for narrow vessels.

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
10	11.6	15.3	155	3	52	19	Non-sterile	1 x 250	1CCC.1	250 unit(s)	

Pasteur pipettes ungraduated

ROTH SELECTION. Material: LDPE.

Precision-formed, unbreakable, non toxic and with controlled, reproducible drop size.




Not a medical device / Not an IVD product

Vol. (ml)	Filling capacity (ml)	Total length (mm)	Drop size (µl)	Drops per ml	Version	Sterility	Pack.	Art. No.	Pack Qty.	
7	7.3	300	50	20	very long	Sterile	indiv.	AEA7.1	100 unit(s)	
5	4.8	137	62	16	wide opening, for viscose media	Non-sterile	1 x 250	E305.1	250 unit(s)	
5	4.8	137	62	16	Wide opening, for cell culture	Sterile	indiv.	EA54.1	250 unit(s)	
3	3.0	155	34	29	long tip, for cells	Non-sterile	1 x 500	E306.1	500 unit(s)	
3	3.0	155	34	29	long tip, for cells	Sterile	50 x 10	EA55.1	500 unit(s)	
2	2.3	225	40	25	long	Sterile	indiv.	AEA8.1	400 unit(s)	
2	2.3	225	40	25	long	Non-sterile	1 x 400	EA56.1	400 unit(s)	
1	0.9	87	40	25	short	Non-sterile	1 x 500	EA57.1	500 unit(s)	

Mixing pipettes ungraduated

ROTH SELECTION. Material: LDPE.

Precision-formed, unbreakable, non toxic and with controlled, reproducible drop size.



Not a medical device / Not an IVD product


Vol. (ml)	Total length (mm)	Drop size (µl)	Drops per ml	Version	Sterility	Art. No.	Pack Qty.	
0.2	124	50	20	standard	Non-sterile	EA52.1	500 unit(s)	
0.2	131	25	40	fine	Non-sterile	EA53.1	500 unit(s)	

Pasteur/disposable pipettes

Pasteur pipettes ungraduated

ROTH SELECTION. Material: LDPE.

Precision-formed, unbreakable, non toxic and with controlled, reproducible drop size.




Not a medical device / Not an IVD product

Vol. (ml)	Filling capacity (ml)	Total length (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
3	3.3	155	40	25	Non-sterile	1 x 500	EA58.1	500 unit(s)	
8	8.0	155	52	19	Non-sterile	1 x 250	E303.1	250 unit(s)	

Pasteur pipettes ultra-fine tip ungraduated

ROTH SELECTION. Material: LDPE.

Precision-formed, unbreakable, non toxic and with controlled, reproducible drop size.



Not a medical device / Not an IVD product

Vol. (ml)	Filling capacity (ml)	Total length (mm)	Outer tip Ø (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
3.5	3.5	157	1.1	20	50	Sterile	individually	AEA5.1	400 unit(s)	
3.5	3.5	157	1.1	20	50	Non-sterile	1 x 500	E304.1	500 unit(s)	
3	3.3	147	1.1	18	55	Non-sterile	1 x 500	EA94.1	500 unit(s)	
1	1.0	51	1.5	22	45	Non-sterile	1 x 500	E250.1	500 unit(s)	
0.5	0.5	36	1.1	25	40	Non-sterile	1 x 500	AEA6.1	500 unit(s)	




Pasteur/disposable pipettes

Pasteur pipettes graduated

ROTH SELECTION. Material: LDPE.

Precision-formed, unbreakable, non toxic and with controlled, reproducible drop size.



Not a medical device / Not an IVD product

Filling volume 0.8 ml

Vol. (ml)	Filling capacity (ml)	Graduation range (ml)	Total length (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
0.3	0.8	0,1-0,3	116	47	21	Non-sterile	1 x 500	EA95.1	500 unit(s)	
0.3	0.8	0,1-0,3	116	47	21	Sterile	25 x 20	EA96.1	500 unit(s)	

Filling volume 1.8 ml

Filling volume 1.8 ml is achieved by pressing the suction ball. To achieve a filling capacity 2 ml, the stem must also be squeezed.

Vol. (ml)	Filling capacity (ml)	Graduation range (ml)	Total length (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
2	1.8	0,5-2	152	41	24	Non-sterile	1 x 500	EA70.1	500 unit(s)	
2	1.8	0,5-2	152	41	24	Sterile	25 x 20	EA71.1	500 unit(s)	

Filling volume 2.2 ml

Vol. (ml)	Filling capacity (ml)	Graduation range (ml)	Total length (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
1	2.2	0,25-1	140	40	25	Non-sterile	1 x 500	EA68.1	500 unit(s)	
1	2.2	0,25-1	140	40	25	Sterile	25 x 20	EA69.1	500 unit(s)	

Filling volume 3.2 ml

Vol. (ml)	Filling capacity (ml)	Graduation range (ml)	Total length (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
3	3.2	0,5-3	155	50	20	Non-sterile	1 x 500	EA65.1	500 unit(s)	
3	3.2	0,5-3	155	50	20	Sterile	individually	EA66.1	500 unit(s)	
3	3.2	0,5-3	155	50	20	Sterile	25 x 20	EA67.1	500 unit(s)	

Filling volume 3.4 ml

Vol. (ml)	Filling capacity (ml)	Graduation range (ml)	Total length (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
1	3.4	0,25-1	155	40	25	Non-sterile	1 x 500	EA62.1	500 unit(s)	
1	3.4	0,25-1	155	40	25	Sterile	individually	EA63.1	500 unit(s)	
1	3.4	0,25-1	155	40	25	Sterile	25 x 20	EA64.1	500 unit(s)	

Filling volume 6.2 ml

Vol. (ml)	Filling capacity (ml)	Graduation range (ml)	Total length (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
5	6.2	1-5	217	54	19	Sterile	12 x 20	AEA4.1	240 unit(s)	
5	6.2	1-5	217	54	19	Non-sterile	1 x 250	EA61.1	250 unit(s)	

Pasteur/disposable pipettes

Pasteur pipettes graduated

ROTH SELECTION. Material: LDPE.

Single-use pipettes with integrated suction ball for straightforward transfer of liquids. With graduation for approximation of dispensed volume. Suitable for media aliquoting thanks to reproducible droplet size.

- Clear graduation
- No slipping with latex gloves
- Precise, fine tip for reproducible results



Not a medical device / Not an IVD product

0.3 ml

For drawing and dispensing small volumes, for transferring into and out of multiwell plates, and for pH measuring instruments (drawing and redispensing of filling solutions for electrodes).

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Graduation range (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
0.3	0.8	1.5	0,1 / 0,2 / 0,3	116	3	45	22	Non-sterile	1 x 500	1CAC.1	500 unit(s)	
0.3	0.8	1.5	0,1 / 0,2 / 0,3	116	3	45	22	Sterile	25 x 20	1CAE.1	500 unit(s)	

1 ml, small suction ball

For drawing and dispensing small volumes, for sample preparation in chromatography (autosampler vials 1 ml), and for pH measuring instruments (addition of acids and bases).

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Graduation range (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
1	2.1	3.9	0,25 / 0,5 / 0,75 / 1	140	4.7	40	25	Non-sterile	1 x 500	1CAH.1	500 unit(s)	
1	2.1	3.9	0,25 / 0,5 / 0,75 / 1	140	4.7	40	25	Sterile	25 x 20	1CAK.1	500 unit(s)	

1 ml, large suction ball

For sample preparation in chromatography (autosampler vials 1 ml).

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Graduation range (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
1	3.5	5.8	0,25 / 0,5 / 0,75 / 1	155	4.6	40	25	Non-sterile	1 x 500	1CAL.1	500 unit(s)	
1	3.5	5.8	0,25 / 0,5 / 0,75 / 1	155	4.6	40	25	Sterile	25 x 20	1CAN.1	500 unit(s)	

1 ml, large suction ball, wide opening

For handling reagents, for microscope slide examinations and for quickly transferring media.

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Graduation range (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
1	4.3	5.8	0,25 / 0,5 / 0,75 / 1	152	2.8	50	20	Non-sterile	1 x 500	1CAP.1	500 unit(s)	

2 ml

For handling reagents, for microscope slide examinations and for quickly transferring media.

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Graduation range (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
2	1.8	4.6	1 / 2	152	7.1	40	25	Non-sterile	1 x 500	1CAT.1	500 unit(s)	
2	1.8	4.6	1 / 2	152	7.1	40	25	Sterile	25 x 20	1CAX.1	500 unit(s)	



Pasteur/disposable pipettes

3 ml

For sample preparation in chromatography (autosampler vials 2 ml). For handling reagents, for microscope slide examinations. Suitable for filling and heat sealing for the purpose of storing liquids.

Vol. (ml)	Filling capacity (ml)	Total volume (ml)	Graduation range (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
3	3.2	7.7	1 / 2 / 3	155	7.9	47	21	Non-sterile	1 x 500	1CAY.1	500 unit(s)	
3	3.2	7.7	1 / 2 / 3	155	7.9	47	21	Sterile	500 x 1	1CC0.1	500 unit(s)	
3	3.2	7.7	1 / 2 / 3	155	7.9	47	21	Sterile	25 x 20	1CC1.1	500 unit(s)	

5 ml

For dispensing larger liquid volumes and quickly transferring media. Suitable for filling and heat sealing for the purpose of storing liquids.

Vol. (ml)	Filling capacity (ml)	Total volume (l)	Graduation range (ml)	Total length (mm)	Ø ext. (mm)	Drop size (µl)	Drops per ml	Sterility	Pack.	Art. No.	Pack Qty.	
5	7.4	13.7	1 / 2 / 3 / 4 / 5	217	3.0	52	19	Non-sterile	1 x 250	1CC2.1	250 unit(s)	

Pasteur pipettes with folding bellows graduated

Kartell. Material: PE.

Dropping pipettes for use in laboratory analytics.



In vitro diagnostic medical device according to IVDD

Vol. (ml)	Graduation range (ml)	Total length (mm)	Ø ext. (mm)	Ø Inner outlet (mm)	Sterility	Pack.	Art. No.	Pack Qty.	
1.5	0,5-1,5	134	8	2	Non-sterile	1 x 100	TT22.1	100 unit(s)	
5	1-5	195	5.5	1.5	Non-sterile	1 x 100	TT23.1	100 unit(s)	

Serological pipettes ROTILABO®

ROTH SELECTION. Material: PS, Plug made of polyester fibre. **Sterile**.

- Easy-to-read graduation
- With additional reverse graduation
- RNase/DNase-free and not pyrogenic



S

Not a medical device / Not an IVD product

Vol. (ml)	Filling capacity (ml)	Graduations (ml)	Total length (mm)	Ø ext. (mm)	Outer tip Ø (mm)	Sterility	Pack.	Art. No.	Pack Qty.	
1	1.04	0.01	279	4.7	3.4	Sterile	10 x (100 x 1)	N231.2	1,000 unit(s)	
1	1.04	0.01	279	4.7	3.4	Sterile	40 x 25	N232.1	1,000 unit(s)	
2	2.08	0.01	270	5.5	3.7	Sterile	10 x (100 x 1)	N236.2	1,000 unit(s)	
2	2.08	0.01	270	5.5	3.7	Sterile	40 x 25	N237.1	1,000 unit(s)	
5	7.5	0.1	348	8	3.4	Sterile	4 x (50 x 1)	N239.2	200 unit(s)	
10	13	0.1	347	10	3.35	Sterile	4 x (50 x 1)	N242.1	200 unit(s)	
10	13	0.1	347	10	3.35	Sterile	20 x 25	N243.2	500 unit(s)	
25	32	0.25	311	15	4.8	Sterile	4 x (50 x 1)	N245.1	200 unit(s)	
50	60	1	368	18	5.5	Sterile	4 x (25 x 1)	ET25.2	100 unit(s)	

Pasteur/disposable pipettes

Serological pipettes CELLSTAR®

Greiner Bio-One. Material: PS, Polyester plugs. **Sterile.**

- Maximum accuracy
- High visual transparency
- Easy-to-read graduation
- Pipette tip design guarantees pipetting without droplets
- Increased volume capacity thanks to the extended negative graduation
- Fitted with a synthetic filter to prevent liquids from being drawn into the pipette
- Pipette colour code according to international standards for easy identification
- Expiration date and lot number printed on each individual piece of packaging
- Individually packed in plastic/plastic packaging with a peel-off function
- **Free of detectable DNase, RNase and human DNA**
- **Free of detectable endotoxins, free of cytotoxic substances**



S

Not a medical device / Not an IVD product

Vol. (ml)	Filling capacity (ml)	Graduations (ml)	Total length (mm)	Ø ext. (mm)	Outer tip Ø (mm)	Colour coding	Sterility	Pack.	Art. No.	Pack Qty.	
1	1.25	0.01	277	4.55	2.75	yellow	Sterile	10 x (100 x 1)	1A74.1	1,000 unit(s)	
2	2.5	0.01	277	5.65	4.05	green	Sterile	10 x (100 x 1)	1A75.1	1,000 unit(s)	
5	7.5	0.1	336.5	8	3.05	blue	Sterile	4 x (50 x 1)	1A76.1	200 unit(s)	
10	13	0.1	340	9.5	3.2	orange	Sterile	4 x (50 x 1)	1A77.1	200 unit(s)	
25	35	0.2	343	14.3	4.76	red	Sterile	4 x (50 x 1)	1A78.1	200 unit(s)	
50	60	0.5	343	18.1	5.3	violet	Sterile	5 x (20 x 1)	1A79.1	100 unit(s)	

Short shape


Short shape pipettes enable users to work comfortably and without straining their back. Practical for use at low heights, e.g. under the flow cabinet.

Vol. (ml)	Filling capacity (ml)	Graduations (ml)	Total length (mm)	Ø ext. (mm)	Outer tip Ø (mm)	Colour coding	Sterility	Pack.	Art. No.	Pack Qty.	
5	8	0.1	221	9.5	3.2	blue	Sterile	4 x (50 x 1)	1A7A.1	200 unit(s)	
10	20	0.2	221	14.3	4.76	orange	Sterile	4 x (50 x 1)	1A7C.1	200 unit(s)	

Serological pipettes

Bürkle. Material: PS. **Sterile.**

For aqueous solutions.



S

Not a medical device / Not an IVD product

Vol. (ml)	Sterility	Pack.	Art. No.	Pack Qty.	
25	Sterile	50 x 1	CPP2.1	50 unit(s)	
50	Sterile	25 x 1	CPP3.1	25 unit(s)	
100	Sterile	10 x 1	CPP4.1	10 unit(s)	



Pasteur/disposable pipettes

Serological pipettes ROTILABO® large opening

ROTH SELECTION. Material: PS, Plug made of polyester fibre. **Sterile.**

For pepping viscose liquids, for the cell culture.

- Easy-to-read graduation
- RNase/DNase-free, non-pyrogenic, non-cytotoxic and non-haemolytic



S

Not a medical device / Not an IVD product

Vol. (ml)	Graduations (ml)	Sterility	Pack.	Art. No.	Pack Qty.	
10	0.1	Sterile	individually	ET30.1	200 unit(s)	

Demeter pipettes sterile

ROTH SELECTION. Material: PS, Plug made of polyester fibre. **Sterile.**

Dilution pipettes for bacteriological examinations in food laboratories and in the dairy products industry.

- RNase/DNase-free, not pyrogenic



S

Not a medical device / Not an IVD product

Vol. (ml)	Graduation range (ml)	Total length (mm)	Ø ext. (mm)	Outer tip Ø (mm)	Sterility	Pack.	Art. No.	Pack Qty.	
2.2	1,0 / 2,0 / 2,1 / 2,2	282.5	6.3	3.7	Sterile	40 x 25	NA44.1	1,000 unit(s)	

Graduated pipettes plastic

Vitlab. Material: PP.

- Blue graduated
- EX calibrated
- Shatter-proof



Cleaning with low alkaline cleaning agents is recommended up to a max. temperature of 60 °C in order to avoid changes in volume.

Vol. (ml)	Graduation (ml)	Limit of error (ml)	Total length (mm)	Ø Outer upper (mm)	Art. No.	Pack Qty.	
1	0,1	±0,02	300	8	XY57.1	12 unit(s)	
2	0,1	±0,02	300	8	XY58.1	12 unit(s)	
5	0,1	±0,05	330	8	XY59.1	12 unit(s)	
10	0,1	±0,10	330	10	XY60.1	12 unit(s)	
10	0,1	±0,10	320	8	XY61.1	12 unit(s)	

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47