

Стойки, боксы, коробка для флаконов, пробирок CARL ROTH ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35


Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47

Racks and boxes



A 121 °C

Test tube stands

Nalgene. Material: PP. Autoclavable: Yes.


Secure hold for test tubes and reaction vials, e.g. when used as a holder in shakers. The contents of the tubes are always clearly visible, and the racks do not float in water baths. Can also be used as a drip rack for tubes, electrophoresis and chromatography plates.

For Ø 10–13 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.
96	8 x 12	179 x 127 x 64	white	E088.1	1 unit(s)
96	8 x 12	179 x 127 x 64	blue	E091.1	1 unit(s)

For Ø 14–17 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.
50	5 x 10	187 x 105 x 70	white	E094.1	1 unit(s)
50	5 x 10	187 x 105 x 70	blue	E097.1	1 unit(s)



Sample stands ROTILABO®

Kartell. Material: PS, foamed.

Ideal for large quantities of test tubes or reaction vials. The holes are deep enough to keep test tubes, vials, etc. straight. Stackable. Extremely lightweight. With finger grips for easy lifting.

For 1.5 ml reaction vials

- Compartment size (hole) Ø 8,5 mm, depth of 15 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.
250	218 x 218 x 33	K828.1	5 unit(s)

For vessel Ø 14 mm


- Compartment size (hole) Ø 14,5 mm, depth of 18 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.
119	218 x 218 x 33	K829.1	5 unit(s)

For vessel Ø 16 mm

- Compartment size (hole) Ø 16,5 mm, depth of 22 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.
96	218 x 218 x 33	K830.1	5 unit(s)



A 121 °C

Test tube stands

Kartell. Material: PP. Autoclavable: Yes.

Temperature stability: 0 to +120 °C

For vessel Ø 10 mm

- Compartment size: Ø 12 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.
12	190 x 60 x 80	EPC0.1	1 unit(s)

For vessel Ø 16 mm

- Compartment size: Ø 18 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.
10	250 x 60 x 74	EPA9.1	1 unit(s)

For vessel Ø 18 mm

- Compartment size: Ø 20 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.
12	190 x 60 x 80	EPC1.1	1 unit(s)
24	375 x 65 x 95	EPC2.1	1 unit(s)



Sample stands foldable for vessel Ø 56 mm, plastic

ROTH SELECTION. Material: PP corrugated, moisture-resistant.

Removable partitions, so can also be used as a transport container or tray. Stackable.

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.
24	4 x 6	406 x 260 x 89	1LEH.1	10 unit(s)



Racks and boxes

Test tube stands

Nalgene. Material: Plastic Resmer®. Autoclavable: Yes.

Moulded from a single piece, with three tiers, contents always clearly visible.

- Alphanumeric labelling on the top level
- Large labelling area on the front
- Does not float in a water bath
- Stackable
- Temperature stability -70 to +121 °C



A121 °C



For glass Ø 13 mm

- Compartment size: 14 x 14 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
36	6 x 6	102 x 102 x 56	white	H529.1	1 unit(s)	
36	6 x 6	102 x 102 x 56	blue	H530.1	1 unit(s)	
36	6 x 6	102 x 102 x 56	red	H532.1	1 unit(s)	
72	6 x 12	200 x 100 x 57	white	7390.1	1 unit(s)	
72	6 x 12	200 x 100 x 57	orange	7391.1	1 unit(s)	
72	6 x 12	200 x 100 x 57	yellow	7392.1	1 unit(s)	
72	6 x 12	200 x 100 x 57	blue	7393.1	1 unit(s)	
72	6 x 12	200 x 100 x 57	green	7394.1	1 unit(s)	
72	6 x 12	200 x 100 x 57	red	H524.1	1 unit(s)	

For glass Ø 16 mm

- Compartment size: 18 x 18 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
36	6 x 6	126 x 126 x 68	white	H533.1	1 unit(s)	
36	6 x 6	126 x 126 x 68	blue	H534.1	1 unit(s)	
36	6 x 6	126 x 126 x 68	green	H535.1	1 unit(s)	
36	6 x 6	126 x 126 x 68	red	H536.1	1 unit(s)	
72	6 x 12	250 x 130 x 70	white	7395.1	1 unit(s)	
72	6 x 12	250 x 130 x 70	orange	7401.1	1 unit(s)	
72	6 x 12	250 x 130 x 70	yellow	7407.1	1 unit(s)	
72	6 x 12	250 x 130 x 70	blue	7431.1	1 unit(s)	
72	6 x 12	250 x 130 x 70	green	7435.1	1 unit(s)	
72	6 x 12	250 x 130 x 70	red	H525.1	1 unit(s)	

For glass Ø 20 mm

- Compartment size: 22 x 22 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
20	4 x 5	128 x 103 x 83	white	H537.1	1 unit(s)	
20	4 x 5	128 x 103 x 83	blue	H538.1	1 unit(s)	
20	4 x 5	128 x 103 x 83	green	H539.1	1 unit(s)	
20	4 x 5	128 x 103 x 83	red	H540.1	1 unit(s)	
40	4 x 10	250 x 100 x 83	white	7438.1	1 unit(s)	
40	4 x 10	250 x 100 x 83	orange	7441.1	1 unit(s)	
40	4 x 10	250 x 100 x 83	yellow	7475.1	1 unit(s)	
40	4 x 10	250 x 100 x 83	blue	7477.1	1 unit(s)	
40	4 x 10	250 x 100 x 83	green	7478.1	1 unit(s)	
40	4 x 10	250 x 100 x 83	red	H526.1	1 unit(s)	

For glass Ø 25 mm

- Compartment size: 26 x 26 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
16	4 x 4	122 x 122 x 75	white	H541.1	1 unit(s)	
16	4 x 4	122 x 122 x 75	blue	H542.1	1 unit(s)	
40	4 x 10	300 x 120 x 92	white	7479.1	1 unit(s)	
40	4 x 10	300 x 120 x 92	orange	7484.1	1 unit(s)	
40	4 x 10	300 x 120 x 92	yellow	7487.1	1 unit(s)	
40	4 x 10	300 x 120 x 92	blue	7488.1	1 unit(s)	
40	4 x 10	300 x 120 x 92	green	7489.1	1 unit(s)	
40	4 x 10	300 x 120 x 92	red	H527.1	1 unit(s)	

Racks and boxes

For glass Ø 30 mm

- Compartment size: 32 x 32 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
24	3 x 8	285 x 110 x 83	white	7492.1	1 unit(s)	
24	3 x 8	285 x 110 x 83	orange	7493.1	1 unit(s)	
24	3 x 8	285 x 110 x 83	yellow	7495.1	1 unit(s)	
24	3 x 8	285 x 110 x 83	blue	7497.1	1 unit(s)	
24	3 x 8	285 x 110 x 83	green	7508.1	1 unit(s)	
24	3 x 8	285 x 110 x 83	red	H528.1	1 unit(s)	
9	3 x 3	109 x 109 x 84	white	H543.1	1 unit(s)	
9	3 x 3	109 x 109 x 84	blue	H544.1	1 unit(s)	
9	3 x 3	109 x 109 x 84	green	H545.1	1 unit(s)	
9	3 x 3	109 x 109 x 84	red	H546.1	1 unit(s)	

Test tube stands ROTILABO®

Kartell. Material: PP. Autoclavable: Yes.

- Does not float in water
- Alphanumerically labelled slots
- Temperature stability -10 to +120 °C



For glass Ø 13 mm

- Compartment size: 13 x 13 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
90	6 x 15	250 x 105 x 60	white	4279.1	1 unit(s)	
90	6 x 15	250 x 105 x 60	red	4283.1	1 unit(s)	
90	6 x 15	250 x 105 x 60	blue	4287.1	1 unit(s)	
90	6 x 15	250 x 105 x 60	yellow	4288.1	1 unit(s)	

For glass Ø 16 mm

- Compartment size: 16 x 16 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
60	5 x 12	250 x 105 x 70	white	4299.1	1 unit(s)	
60	5 x 12	250 x 105 x 70	red	4327.1	1 unit(s)	
60	5 x 12	250 x 105 x 70	blue	4367.1	1 unit(s)	
60	5 x 12	250 x 105 x 70	yellow	4383.1	1 unit(s)	

For glass Ø 20 mm

- Compartment size: 21 x 21 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
40	4 x 10	250 x 105 x 70	white	4388.1	1 unit(s)	
40	4 x 10	250 x 105 x 70	red	4398.1	1 unit(s)	
40	4 x 10	250 x 105 x 70	blue	4399.1	1 unit(s)	
40	4 x 10	250 x 105 x 70	yellow	4430.1	1 unit(s)	

For glass Ø 25 mm

- Compartment size: 26 x 26 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
40	4 x 10	300 x 125 x 85	white	4448.1	1 unit(s)	
40	4 x 10	300 x 125 x 85	red	4449.1	1 unit(s)	
40	4 x 10	300 x 125 x 85	blue	4450.1	1 unit(s)	
40	4 x 10	300 x 125 x 85	yellow	4466.1	1 unit(s)	



Racks and boxes

For glass Ø 30 mm

- Compartment size: 31 x 31 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
24	3 x 8	300 x 110 x 85	white	K833.1	1 unit(s)	
24	3 x 8	300 x 110 x 85	red	K834.1	1 unit(s)	
24	3 x 8	300 x 110 x 85	blue	K835.1	1 unit(s)	
24	3 x 8	300 x 110 x 85	yellow	K836.1	1 unit(s)	

Sample stands

ROTH SELECTION. Material: Steel, epoxy-coated. Autoclavable: Yes.

Temperature stability: -40 to +121 °C



A121 °C

For tube Ø 10–13 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
72	6 x 12	205 x 105 x 65	blue	AXH0.1	1 unit(s)	
108	9 x 12	204 x 161 x 65	blue	1LHX.1	1 unit(s)	

For tube Ø 13–16 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
40	4 x 10	200 x 83 x 59	blue	1LHY.1	1 unit(s)	
48	4 x 12	237 x 127 x 65	blue	1LK0.1	1 unit(s)	
72	6 x 12	241 x 127 x 65	blue	AXH2.1	1 unit(s)	
108	9 x 12	237 x 127 x 65	blue	1LK1.1	1 unit(s)	
108	9 x 12	237 x 127 x 65	red	AXH3.1	1 unit(s)	

For tube Ø 16–20 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
15	3 x 5	127 x 76 x 76	blue	AXH4.1	1 unit(s)	
40	4 x 10	240 x 90 x 80	blue	1LK5.1	1 unit(s)	
80	8 x 10	240 x 190 x 80	blue	1LK6.1	1 unit(s)	

For tube Ø 20–25 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
24	3 x 8	273 x 105 x 80	blue	1LK7.1	1 unit(s)	
40	4 x 10	290 x 118 x 80	blue	1LK8.1	1 unit(s)	

For tube Ø 25–30 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
8	2 x 4	147 x 80 x 89	blue	AXH5.1	1 unit(s)	
48	6 x 8	277 x 209 x 98	blue	1LK9.1	1 unit(s)	


For tube Ø 35–40 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
24	3 x 8	357 x 139 x 98	blue	1LKA.1	1 unit(s)	
48	6 x 8	350 x 260 x 98	blue	1LKC.1	1 unit(s)	

For tube Ø 55–60 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
6	2 x 3	2050 x 1352 x 104	blue	1LKE.1	1 unit(s)	

Racks and boxes



A 121 °C

Test tube stands ROTILABO® stainless steel

ROTH SELECTION. Material: Stainless steel 18/10 (1.4301), electrolytically polished. Autoclavable: Yes.

With four feet, 14 mm clearance from the floor.

For test tube Ø 16 mm


- Compartment size: 18 x 18 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
25	5 x 5	115 x 115 x 70	1C6K.1	1 unit(s)	
100	10 x 10	200 x 200 x 70	ENX9.1	1 unit(s)	

For test tube Ø 30 mm

- Compartment size: 30 x 30 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
12	2 x 6	200 x 70 x 70	1C6L.1	1 unit(s)	



Test tube stands ROTILABO® acrylic

ROTH SELECTION. Material: PMMA (Acrylic glass).

Hole Ø 18 mm

Base plate with holes (hole Ø 14 mm).

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
6	1 x 6	165 x 55 x 75	CPT9.1	1 unit(s)	
24	3 x 8	220 x 80 x 79	CPX0.1	1 unit(s)	
12	2 x 6	186 x 60 x 80	5119.1	1 unit(s)	
33	3 x 11	310 x 80 x 80	CPX1.1	1 unit(s)	

Hole Ø 36 mm

Base plate without holes.

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
12	2 x 6	275 x 100 x 80	5416.1	1 unit(s)	



Test tube stands ROTILABO®

ROTH SELECTION.

Stainless steel compartment size 14 x 14 mm

Material: Stainless steel. Autoclavable: Yes.

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
12	2 x 6	101 x 44 x 70	H439.1	1 unit(s)	
24	2 x 12	197 x 44 x 70	H440.1	1 unit(s)	
48	4 x 12	197 x 76 x 70	H441.1	1 unit(s)	

Stainless steel compartment size 18 x 18 mm

Material: Stainless steel. Autoclavable: Yes.

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
12	2 x 6	130 x 50 x 70	H442.1	1 unit(s)	
24	2 x 12	250 x 50 x 70	H443.1	1 unit(s)	
48	4 x 12	250 x 90 x 70	H444.1	1 unit(s)	

PA-coated compartment size 18 x 18 mm

Material: Steel PA-coated.

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
12	2 x 6	130 x 50 x 70	L385.1	1 unit(s)	
20	2 x 10	210 x 50 x 70	L386.1	1 unit(s)	
24	2 x 12	250 x 50 x 70	L387.1	1 unit(s)	
24	4 x 6	130 x 90 x 70	L389.1	1 unit(s)	
36	3 x 12	250 x 70 x 70	L388.1	1 unit(s)	
48	4 x 12	250 x 90 x 70	L390.1	1 unit(s)	
100	10 x 10	210 x 210 x 70	L391.1	1 unit(s)	



A 121 °C

Test tube stands ROTILABO® round

Bochem. Material: Stainless steel 18/10. Autoclavable: Yes.

- Compartment size (hole): Ø 20 mm
- Suitable for: Test tubes Ø 18 mm

No. of slots	Ø x H (mm)	Art. No.	Pack Qty.	
9	98 x 165	YT05.1	1 unit(s)	
12	116 x 165	YT06.1	1 unit(s)	
45	194 x 165	YT07.1	1 unit(s)	



Racks and boxes



A 121 °C

Test tube stands stainless steel

Bochem. Material: Stainless steel 18/10. Autoclavable: Yes.

Hole Ø 17 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
12	2 x 6	145 x 55 x 90	PH51.1	1 unit(s)	
24	2 x 12	278 x 55 x 90	PH52.1	1 unit(s)	
48	4 x 12	287 x 100 x 90	PH53.1	1 unit(s)	

Hole Ø 20 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
12	2 x 6	164 x 62 x 90	PH54.1	1 unit(s)	
24	2 x 12	310 x 60 x 90	PH55.1	1 unit(s)	



A 121 °C

Test tube stands ROTILABO® stainless steel Z shape

ROTH SELECTION. Material: Stainless steel 18/10 (Remanit 4301). Autoclavable: Yes. Stackable.

Hole Ø 12.5 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
50	155 x 78 x 50	X115.1	1 unit(s)	

Hole Ø 14.5 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
50	188 x 95 x 58	X116.1	1 unit(s)	

Hole Ø 17.5 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
25	110 x 110 x 40	PX26.1	1 unit(s)	
50	220 x 115 x 63	X117.1	1 unit(s)	

Hole Ø 30 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
15	220 x 115 x 63	PX27.1	1 unit(s)	



Test tube stands wood

Hecht. Material: Wood (Oak) varnished.

For glass Ø 18 mm

For test tubes with a height up of at least 100 mm. Optionally with drip pegs for drying test tubes after cleaning.

Version	No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
one tier, with 6 drip pegs	6	210 x 55 x 120	AEE1.1	1 unit(s)	
one tier, without drip pegs	6	210 x 55 x 120	AEE0.1	1 unit(s)	
two tiers, with 6 drip pegs	12	210 x 55 x 200	AEE5.1	1 unit(s)	
two tiers, without drip pegs	12	210 x 55 x 200	AEE4.1	1 unit(s)	
two tiers, without drip pegs	24	370 x 55 x 200	AEE8.1	1 unit(s)	

For glass Ø 22 mm

For test tubes with a height up of at least 100 mm. Optionally with drip pegs for drying test tubes after cleaning.

Version	No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
one tier, without drip pegs	6	230 x 55 x 120	AEE2.1	1 unit(s)	
two tiers, with 6 drip pegs	12	230 x 75 x 200	AEH0.1	1 unit(s)	

For glass Ø 30 mm

For test tubes with a height of at least 160 mm.

Version	No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
one tier, without drip pegs	6	270 x 55 x 180	AEE3.1	1 unit(s)	

Racks and boxes



Test tube stands extra large

ROTH SELECTION. Material: PP. Autoclavable: Yes.
For applications involving high sample throughput.
With alphanumeric labelling

For tube Ø 10–13 mm

No. of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
216	337 x 235 x 91	blue	TP50.1	1 unit(s)	
216	337 x 235 x 91	white	TP51.1	1 unit(s)	

For tube Ø 13–16 mm

No. of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
120	337 x 267 x 91	blue	TP55.1	1 unit(s)	
120	337 x 267 x 91	white	TP54.1	1 unit(s)	

For tube Ø 18–20 mm

No. of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
120	377 x 301 x 91	blue	1LKL.1	1 unit(s)	
120	377 x 301 x 91	white	1LKN.1	1 unit(s)	



Test tube stands slant rack

Nalgene. Material: PC. Autoclavable: Yes.

For simple and safe production of slant agar with two different inclination angles, for incubating liquid tissue cultures and for use in a water bath (not floating). Tubes containing culture media can be autoclaved in the rack in an upright position.

- Slanted rack enables tilting at an inclination angle of 5° or 20° so that culture media can harden at an angle
- 5° inclination angle: Larger surface, e.g. for aerobic cultures
- 20° inclination angle: Smaller surface, e.g. for anaerobic cultures; reduces the risk of culture media drying out during long-term storage

For glass Ø 16 mm

- Compartment size (hole): Ø 17 mm

No. of slots	Arrangement of slots	Art. No.	Pack Qty.	
40	4 x 10	E100.1	1 unit(s)	

For glass Ø 20 mm

- Compartment size (hole): Ø 20 mm

No. of slots	Arrangement of slots	Art. No.	Pack Qty.	
40	4 x 10	E101.1	1 unit(s)	

Алматы (7273)495-231	Иваново (4932)77-34-06	Магнитогорск (3519)55-03-13	Ростов-на-Дону (863)308-18-15	Тольятти (8482)63-91-07
Ангарск (3955)60-70-56	Ижевск (3412)26-03-58	Москва (495)268-04-70	Рязань (4912)46-61-64	Томск (3822)98-41-53
Архангельск (8182)63-90-72	Иркутск (395)279-98-46	Мурманск (8152)59-64-93	Самара (846)206-03-16	Тула (4872)74-02-29
Астрахань (8512)99-46-04	Казань (843)206-01-48	Набережные Челны (8552)20-53-41	Саранск (8342)22-96-24	Тюмень (3452)66-21-18
Барнаул (3852)73-04-60	Калининград (4012)72-03-81	Нижний Новгород (831)429-08-12	Санкт-Петербург (812)309-46-40	Ульяновск (8422)24-23-59
Белгород (4722)40-23-64	Калуга (4842)92-23-67	Новокузнецк (3843)20-46-81	Саратов (845)249-38-78	Улан-Удэ (3012)59-97-51
Благовещенск (4162)22-76-07	Кемерово (3842)65-04-62	Ноябрьск (3496)41-32-12	Севастополь (8692)22-31-93	Уфа (347)229-48-12
Брянск (4832)59-03-52	Киров (8332)68-02-04	Новосибирск (383)227-86-73	Симферополь (3652)67-13-56	Хабаровск (4212)92-98-04
Владивосток (423)249-28-31	Коломна (4966)23-41-49	Омск (3812)21-46-40	Ставрополь (8652)20-65-13	Чебоксары (8352)28-53-07
Владикавказ (8672)28-90-48	Кострома (4942)77-07-48	Орел (4862)44-53-42	Сочи (862)225-72-31	Челябинск (351)202-03-61
Владимир (4922)49-43-18	Краснодар (861)203-40-90	Оренбург (3532)37-68-04	Сыктывкар (8212)25-95-17	Череповец (8202)49-02-64
Волгоград (844)278-03-48	Красноярск (391)204-63-61	Пенза (8412)22-31-16	Тамбов (4752)50-40-97	Чита (3022)38-34-83
Вологда (8172)26-41-59	Курск (4712)77-13-04	Петрозаводск (8142)55-98-37	Сургут (3462)77-98-35	Якутск (4112)23-90-97
Воронеж (473)204-51-73	Курган (3522)50-90-47	Псков (8112)59-10-37	Тверь (4822)63-31-35	Ярославль (4852)69-52-93
Екатеринбург (343)384-55-89	Липецк (4742)52-20-81	Пермь (342)205-81-47		
Россия (495)268-04-70	Казахстан (772)734-952-31	Киргизия (996)312-96-26-47		