

Резервуары для реагентов CARL ROTH Bio-pure sterile, ROTILABO, Trifill Multi-channel

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47




Reagent reservoirs

Reagent reservoirs

ROTH SELECTION. Material: PS.

- Minimum reagent loss through conical collecting channel on the bottom
- Graduated volume markings on the inner wall
- Reinforced side walls and wider base for better stability



Volume 25 ml

Pic.	Colour	Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.
(1)	white	146	57	26	Non-sterile	10 x 10	CYN5.1	100 unit(s)
(1)	white	146	57	26	Sterile	10 x 10	HT65.1	100 unit(s)
(1)	white	146	57	26	Sterile	indiv.	HT66.1	50 unit(s)

Volume 25 ml partitioned

Reservoir partitioned in 2 chambers with capacity 5 or 10 ml, suitable for 8 and 4-channel pipettes.

Pic.	Colour	Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.
(2)	white	146	57	26	Non-sterile	10 x 10	CYN4.1	100 unit(s)
(2)	white	146	57	26	Sterile	10 x 10	HT67.1	100 unit(s)
(2)	white	146	57	26	Sterile	indiv.	HT68.1	50 unit(s)

Volume 100 ml

Pic.	Colour	Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.
(3)	white	146	83	32	Non-sterile	10 x 10	CYN6.1	100 unit(s)
(3)	white	146	83	32	Sterile	10 x 10	HT69.1	100 unit(s)
(3)	white	146	83	32	Sterile	indiv.	HT70.1	50 unit(s)

Volume 100 + 15 ml sterile

With additional reservoir at the side.

Pic.	Colour	Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.
(4)	white	146	102	32	Sterile	10 x 10	HT71.1	100 unit(s)

Reagent reservoirs Bio-pure™ sterile

ROTH SELECTION. Material: PET, biodegradable.

- Minimum reagent loss through conical collecting channel on the bottom
- Reinforced side walls and wider base for better stability
- Packed in biodegradable bags and recyclable boxes



Volume 25 ml

Pic.	Colour	Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.
(1)	white	146	57	26	Sterile	10 x 10	AXC7.1	100 unit(s)

Volume 25 ml (partitioned)

Reservoir partitioned in 2 chambers with capacity 5 or 10 ml, suitable for 8 and 4-channel pipettes.

Pic.	Colour	Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.
(2)	white	146	57	26	Sterile	10 x 10	AXC8.1	100 unit(s)

Volume 100 ml

Pic.	Colour	Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.
(3)	white	146	83	32	Sterile	10 x 10	AXC9.1	100 unit(s)

Reagent reservoirs

Reagent reservoirs

ROTH SELECTION. Material: PS.

- Compatible with 8 and 12-channel pipettes
- Robust design for safe handling
- V-shaped base for non-loss pipetting
- Graduations on the inside wall
- Spouts at the corners
- **RNase/DNase-free, pyrogen-free**



Volume 25 ml

Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.	
147	61	25	Non-sterile	1 x 100	EKT6.1	100 unit(s)	
147	61	25	Sterile	20 x 5	EKT7.1	100 unit(s)	
147	61	25	Sterile	indiv.	EKT8.1	50 unit(s)	

Volume 50 ml

Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.	
147	76	28	Non-sterile	1 x 100	EKT9.1	100 unit(s)	
147	76	28	Sterile	20 x 5	EKX0.1	100 unit(s)	
147	76	28	Sterile	indiv.	EKX1.1	50 unit(s)	

Volume 100 ml

Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.	
147	80	32	Non-sterile	1 x 100	EKX2.1	100 unit(s)	
147	80	32	Sterile	20 x 5	EKX3.1	100 unit(s)	
147	80	32	Sterile	indiv.	EKX4.1	50 unit(s)	

Reagent reservoirs ROTILABO®

ROTH SELECTION.

- Graduations on the inside wall
- Spout for pouring without loss of substance or dripping
- V-shaped channel in the base for pipetting without loss of substance
- Strengthened side walls and a wide base ensure stability
- Stackable
- **RNase-/DNase-free, pyrogen-free**



Volume 5 ml

Graduated at 5 ml.

Length (mm)	Inner length (mm)	Width (mm)	Inner width (mm)	Height (mm)	Inner height (mm)	Sterility	Pack.	Art. No.	Pack Qty.	
107	72	54	21	19	18	Sterile	50 x 1	KLE5.1	50 unit(s)	

Volume 10 ml

Graduated at 5 ml/10 ml.

Length (mm)	Inner length (mm)	Width (mm)	Inner width (mm)	Height (mm)	Inner height (mm)	Sterility	Pack.	Art. No.	Pack Qty.	
107	72	54	21	19	18	Non-sterile	1 x 100	KLE6.1	100 unit(s)	
107	72	54	21	19	18	Sterile	100 x 1	KLE7.1	100 unit(s)	



Reagent reservoirs

Reagent reservoirs ROTILABO® transparent

ROTH SELECTION. Material: PVC.

- Graduations on the inside wall
- Spout for pouring without loss of substance or dripping
- V-shaped channel in the base for pipetting without loss of substance
- Strengthened side walls and a wide base ensure stability
- Stackable
- **RNase-/DNase-free, pyrogen-free**



Volume 50 ml

Graduated at 25 ml/50 ml.

Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.
102	28	30	Sterile	20 x 5	KLE8.1	100 unit(s)
102	28	30	Non-sterile	1 x 100	KLE9.1	100 unit(s)

Volume 100 ml

Graduated at 50 ml/100 ml.

Length (mm)	Width (mm)	Height (mm)	Sterility	Pack.	Art. No.	Pack Qty.
109	48	32	Non-sterile	1 x 100	KLH0.1	100 unit(s)

Reagent reservoir cooler ROTILABO®

ROTH SELECTION. Material: PP.

Filled with coolant. Suitable for reagent reservoirs ROTILABO®, transparent. Cools media for up to three hours at a temperature of 5 °C. Stackable.



For 50 ml reservoirs

Colour	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.
violet	144	74	45	KLE1.1	2 unit(s)
blue	144	74	45	KLE2.1	2 unit(s)

For 100 ml reservoirs

Colour	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.
violet	150	90	43	KLE3.1	2 unit(s)
blue	150	90	43	KLE4.1	2 unit(s)

Reagent reservoirs Trifill Multi-channel

ROTH SELECTION. Material: PP. Autoclavable: Yes.

Trifill is a combination of three reservoir units. With V-shaped, graduated liquid container, max. volume 50 ml, as well as 8 and 12 graduated liquid reservoirs. Max. filling capacity of 5 ml per numbered, V-shaped chamber.



A 121 °C

Colour	Length (mm)	Width (mm)	Height (mm)	Sterility	Art. No.	Pack Qty.
white	120	111	40	Non-sterile	AA87.1	10 unit(s)

Reagent reservoirs

Reagent reservoirs 175 ml

ROTH SELECTION. Material: PP, transparent. Autoclavable: Yes.

- Volume 175 ml
- Graduated at: 25, 50, 75, 100, 125, 150 ml
- Stable
- Reusable
- Compatible with 8/12-channel pipettes as well as 96-well pipetting systems



A 121 °C

With lid

Removable lid, easy to open.

Inner length (mm)	Inner width (mm)	Inner height (mm)	Sterility	Pack.	Art. No.	Pack Qty.
113.6	75.5	28.7	Sterile	20 x 5	APX0.1	100 unit(s)
113.6	75.5	28.7	Non-sterile	10 x 10	APX1.1	100 unit(s)

Without lid

Inner length (mm)	Inner width (mm)	Inner height (mm)	Sterility	Pack.	Art. No.	Pack Qty.
113.6	75.5	28.7	Sterile	20 x 5	APX2.1	100 unit(s)
113.6	75.5	28.7	Non-sterile	10 x 10	APX3.1	100 unit(s)

Reagent reservoirs for 8-channel pipettes

ROTH SELECTION. Material: PS. Non-sterile.

Single compartments with round bottom for minimal reagent loss.
Capacity per well approx. 7 ml.



Colour	Length (mm)	Width (mm)	Sterility	Art. No.	Pack Qty.
white	126.5	84	Non-sterile	NA12.1	100 unit(s)

Reagent reservoirs for 8/12-channel pipettes

ROTH SELECTION. Material: PP. Autoclavable: Yes.

With 12 numbered, V-shaped wells. Capacity per well approx. 3 ml.
Temperature stable from -80 to +121 °C.



A 121 °C



Delivery incl. lid.

Sterility	Pack.	Art. No.	Pack Qty.
Non-sterile	indiv.	EKT5.1	50 unit(s)
Sterile	indiv.	EKT4.1	20 unit(s)

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Россия (495)268-04-70

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Казахстан (772)734-952-31

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Киргизия (996)312-96-26-47

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93