

Банки для образцов, проб CARL ROTH ROTILABO, SpecTainer, SecurTainer III

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47

Sample cans

Square boxes ROTILABO®

ROTH SELECTION. Material: HDPE, slip lid made of LDPE.

Rectangular, slightly tapered. Stackable. Food-safe, dishwasher-safe and microwave-safe. Temperature stability -20 to +70 °C.



Not a medical device / Not an IVD product

Flat form

Vol. (ml)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
500	103	103	64	1597.1	5 unit(s)	
1,000	208	103	64	1598.1	5 unit(s)	
2,000	208	208	64	A274.1	5 unit(s)	

High form

Vol. (ml)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
1,000	103	103	127	1599.1	5 unit(s)	
1,500	208	103	94	A272.1	5 unit(s)	
2,000	255	103	94	A273.1	5 unit(s)	
3,300	208	208	94	A275.1	5 unit(s)	

Square boxes ROTILABO®

ROTH SELECTION. Material: PP translucent, snap-on lid made of LDPE, blue.

For storing granulates, powders and paste-like substances.

Temperature stability -40 to +85 °C.

Dishwasher-safe and microwave-safe. Suitable for foodstuffs.



Not a medical device / Not an IVD product

Vol. (ml)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
300	110	110	44	HEK2.1	10 unit(s)	
450	150	100	50	HEK3.1	10 unit(s)	
500	230	140	32	HEK4.1	10 unit(s)	
800	180	120	65	HEK6.1	10 unit(s)	
1,000	235	140	52	HEK7.1	10 unit(s)	
1,500	230	150	65	HEK8.1	10 unit(s)	
2,000	280	180	65	HEK9.1	5 unit(s)	
2,500	250	158	92	HEL0.1	5 unit(s)	
3,000	280	180	80	HEL1.1	5 unit(s)	
4,000	280	180	115	HEL2.1	5 unit(s)	



Sample cans

Square boxes ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes.
With lid. For storing solid and liquid substances. Suitable for foodstuffs.
7.5 litre container with two folding handles on the sides.



A121 °C

Not a medical device / Not an IVD product

Vol. (ml)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.
460	80	80	80	9502.1	1 unit(s)
2,500	140	140	140	9506.1	1 unit(s)
7,500	200	200	200	9509.1	1 unit(s)

Storage tins ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10.
Odour- and moisture-proof. Suitable for foodstuffs.



Not a medical device / Not an IVD product

With hinged lid

Lid made of acrylic glass with silicone seal, leakproof.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.
700	125	95	YX66.1	1 unit(s)
1,000	125	125	YX67.1	1 unit(s)
1,250	125	155	YX68.1	1 unit(s)
1,700	125	195	YX69.1	1 unit(s)

With slip lid

Lid made of stainless steel with acrylic glass viewing window.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.
500	80	110	YX70.1	1 unit(s)
800	90	125	YX71.1	1 unit(s)
2,000	125	170	YX72.1	1 unit(s)
4,500	175	210	YX73.1	1 unit(s)

Sample cans

Square boxes ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10, matt polished.
Autoclavable: Yes.
Stackable. Suitable for foodstuffs.



A 121 °C



Delivery without lid, please order separately.

Not a medical device / Not an IVD product

Vol. (ml)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.
500	95	95	60	YH09.1	1 unit(s)
1,000	195	95	60	YH10.1	1 unit(s)
2,000	195	195	60	YH11.1	1 unit(s)

Accessories lids

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes.

A 121 °C

Not a medical device / Not an IVD product

Suitable for	Art. No.	Pack Qty.
Deep-freeze boxes, Art. No. YH09.1	YH12.1	1 unit(s)
Deep-freeze boxes, Art. No. YH10.1	YH13.1	1 unit(s)
Deep-freeze boxes, Art. No. YH11.1	YH14.1	1 unit(s)

Wide-neck container

Kartell. Material: HDPE.

With press-in plugs and screw cap. Ideal for storing powders or granulated substances. Easy to fill and clean thanks to large wide-neck opening. Suitable for foodstuffs.



In vitro diagnostic medical device according to IVDD

Vol. (ml)	Ø ext. (mm)	Ø Inner neck (mm)	Height (mm)	Art. No.	Pack Qty.
70	50	36	60	HT55.1	100 unit(s)
120	69	36	71	HT56.1	100 unit(s)
250	56	50	94	HT57.1	50 unit(s)
500	86	70	107	HT58.1	50 unit(s)
1,000	111	85	128	HT59.1	20 unit(s)
2,000	111	87	235	HT60.1	10 unit(s)



Sample cans

Sample tub

Nalgene. Autoclavable: Yes.

Easy to fill and clean thanks to large wide-neck opening.

Free of undercuts.

- 125 ml to 1000 ml tubs can be stacked thanks to recess in the lid
- 125 ml and 250 ml sizes suitable for mixers with Mason thread



A 121 °C



Note: Do not use to transport liquids.

Not a medical device / Not an IVD product

Type 2116

Material: PC transparent, screw closure made of PP. Shockproof, non-toxic. Particularly suitable for storage in refrigerators and cold storage rooms.

Temperature stable from -135 to +130 °C.

Vol. (ml)	Ø ext. (mm)	Outer lid Ø (mm)	Height (mm)	Art. No.	Pack Qty.
15	32	38	46	K682.1	4 unit(s)
30	36	43	48	K683.1	4 unit(s)
60	48	53	48	K684.1	4 unit(s)
125	69	70	74	K685.1	4 unit(s)
250	66	70	119	K686.1	4 unit(s)
500	117	120	86	K687.1	4 unit(s)
1,000	117	120	142	K688.1	4 unit(s)

Type 2117

Material: PMP transparent, screw closure made of PP. Good chemical resistance. Particularly suitable as an archive container: Formaldehyde loss only approx. 2% per year.

Vol. (ml)	Ø ext. (mm)	Outer lid Ø (mm)	Height (mm)	Art. No.	Pack Qty.
60	48	53	48	K693.1	4 unit(s)
125	69	70	74	K694.1	4 unit(s)
250	66	70	119	K695.1	4 unit(s)
500	117	120	86	K696.1	4 unit(s)
1,000	117	120	145	K697.1	4 unit(s)

Type 2118

Material: PP translucent, screw closure made of PP. Good chemical resistance. Temperature stable from 0 to +130 °C.

Vol. (ml)	Ø ext. (mm)	Outer lid Ø (mm)	Height (mm)	Art. No.	Pack Qty.
15	32	38	46	K698.1	12 unit(s)
30	36	43	48	K699.1	12 unit(s)
60	48	53	48	K700.1	12 unit(s)
125	64	70	74	K701.1	12 unit(s)
250	66	70	119	K702.1	6 unit(s)
500	117	120	86	K703.1	6 unit(s)
1,000	117	120	145	K704.1	6 unit(s)

Sample cans

Wide-neck container ROTILABO®

ROTH SELECTION.

Suitable for packaging and storing liquids and solids. Easy to fill thanks to wide-neck opening. Suitable for foodstuffs. Height without cap.



Not a medical device / Not an IVD product

Clear PVC

Material: PVC clear, screw closure made of PP with PE foam insert.

Vol. (ml)	Ø Inner neck (mm)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.
50	23	38	38	60.5	H312.1	25 unit(s)
100	31	46	46	71	H313.1	24 unit(s)
200	41	59.5	59.5	88.5	H314.1	14 unit(s)
300	41	67.5	67.5	95	H315.1	15 unit(s)
500	54	80.5	80.5	108	H316.1	10 unit(s)
750	69	91.5	91.5	120	H317.1	10 unit(s)
1,000	69	97.5	97.5	144	H318.1	8 unit(s)
2,000	69	122	122	187	H319.1	4 unit(s)

Brown PVC

Material: PVC brown-transparent, screw closure made of PP with PE foam insert.

Vol. (ml)	Ø Inner neck (mm)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.
50	23	38	38	60.5	H320.1	25 unit(s)
100	31	46	46	71	H321.1	24 unit(s)
200	41	59.5	59.5	88.5	H322.1	14 unit(s)
300	41	67.5	67.5	95	H323.1	15 unit(s)
500	54	80.5	80.5	108	H324.1	10 unit(s)
750	69	91.5	91.5	120	H325.1	10 unit(s)
1,000	69	97.5	97.5	144	H326.1	8 unit(s)
2,000	69	122	122	187	H327.1	4 unit(s)

PETG

Material: PETG transparent, screw closure made of PP with PE foam insert.

Vol. (ml)	Ø Inner neck (mm)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.
50	23	38	38	60.5	H328.1	25 unit(s)
100	31	46	46	71	H329.1	24 unit(s)
200	41	59.5	59.5	88.5	H330.1	14 unit(s)
500	54	80.5	80.5	108	H331.1	10 unit(s)
1,000	69	97.5	97.5	144	H332.1	8 unit(s)
2,000	69	122	122	187	T408.1	4 unit(s)



Sample cans

Sample tub ROTILABO® with screw closure

ROTH SELECTION. Material: HDPE.

Suitable for storing liquids or solids. Thick-walled.



Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Ø Internal (mm)	Height (mm)	Art. No.	Pack Qty.	
5	21	17	30	K760.1	10 unit(s)	
10	21	17	52	K761.1	10 unit(s)	
30	35	31	53	K762.1	10 unit(s)	
60	35	31	92	K763.1	10 unit(s)	
90	55	50	62	K764.1	5 unit(s)	
180	55	50	108	K765.1	5 unit(s)	

Sample tub with screw closure

Vitlab. Material: PP.

For sampling, transporting and storing samples. Suitable for foodstuffs.



Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Thread	Art. No.	Pack Qty.	
30	38	54	GL 40	YA17.1	10 unit(s)	
60	38	90	GL 40	YA18.1	10 unit(s)	
90	54	62	GL 56	YA19.1	10 unit(s)	
180	54	112	GL 56	YA20.1	10 unit(s)	

Wide-neck jar Packco

CurTec. Material: HDPE.

White jar, smooth inside, dark red screw-on lid and semi-transparent stopper. Screw-on lid can be easily opened and closed with a quarter turn. For storing and dispatching solids and pastes.

- cGMP-certified/ISO 15378 (primary packaging materials for pharmaceutical products)
- With UN approval (1H2) for solids
- Compliant with food regulations
- Resistant to steam
- With tamper-evident cap



Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Ø Internal (mm)	Height (mm)	UN approval	Art. No.	Pack Qty.	
300	99	96.8	55.4	1H2/X0,4/Y0,6/Z0,9/S	1A5K.1	10 unit(s)	
500	99	96.8	82.4	1H2/X0,7/Y1,0/Z1,5/S	1A5L.1	10 unit(s)	
650	99	96.8	104.4	1H2/X0,9/Y1,3/Z1,9/S	1A5N.1	10 unit(s)	
1,000	99	96.8	157.4	1H2/X1,3/Y2,0/Z3,0/S	1A5P.1	10 unit(s)	
1,300	99	96.8	197.4	1H2/X1,8/Y2,0/Z3,0/S	1A5T.1	10 unit(s)	

Sample cans

Sample tub ROTILABO®

ROTH SELECTION.

With tamper-evident cap. For powders, pastes, granulates and other bulk goods. Suitable for foodstuffs.



Not a medical device / Not an IVD product

With snap-on lid

Material: HDPE lid made of PE.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
30	35	54	5553.1	10 unit(s)	
75	45	67	5602.1	10 unit(s)	
150	64	67	5623.1	10 unit(s)	
300	64	110	5642.1	10 unit(s)	

With slip lid

Material: PP lid made of PE.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
750	108	96	5657.1	10 unit(s)	
1,250	108	155	5694.1	5 unit(s)	

Sample tub with snap-on lid

Nalgene. Material: LDPE, translucent.

With friction-fit, snap-on lid. Versatile use as a sample container, beaker for pH measurements, weighing container, etc. Temperature stability -100 to +80 °C



Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
5	23.6	24.9	ATK3.1	12 unit(s)	
12	27.1	35.6	ATK4.1	12 unit(s)	
18	27.1	52.6	ATK5.1	12 unit(s)	
28	27.1	84.6	ATK6.1	12 unit(s)	
50	37.3	70.9	ATK7.1	12 unit(s)	
75	37.3	107.3	ATK8.1	12 unit(s)	



Sample cans



A 121 °C



Sample tub ROTILABO® with snap-on lid

ROTH SELECTION. Material: PP. Autoclavable: Yes.

With lid attached, can be attached vertically. Sterile (γ-sterilised) or non-sterile. With labelling area and graduation (except Art. No. EP22.1). 10 ml graduation. Leakproof. Temperature stability -20 to +121 °C. Autoclavable with open lid.

Not a medical device / Not an IVD product

Vol. (ml)	Sterility	Ø Internal (mm)	Height (mm)	Art. No.	Pack Qty.
45	Non-sterile	29	85	EP22.1	650 unit(s)
50	Non-sterile	31	78	EP24.1	650 unit(s)
50	Sterile	31	78	EP25.1	650 unit(s)
90	Non-sterile	43	80	EP26.1	350 unit(s)
90	Sterile	43	80	EP27.1	350 unit(s)
300	Non-sterile	75	88	EP28.1	240 unit(s)
300	Sterile	75	88	EP29.1	240 unit(s)



Sample tub with screw closure

Hünorsdorff. Material: PP, Lid made of PE.

For storing solid samples. Temperature stability 0 to +100 °C. Free of undercuts, fully drainable.

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.
250	104	50	XY70.1	1 unit(s)
500	104	100	XY71.1	1 unit(s)
1,000	126	126	XY72.1	1 unit(s)



A 121 °C

Wide-neck jar

Bürkle. Material: PP. Autoclavable: Yes.

Well-designed packaging box with screw cap. Free of undercuts, total emptying possible.

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.
12	25	31	1CX4.1	10 unit(s)
25	31	40	1CX5.1	10 unit(s)
60	41	56	1CX6.1	10 unit(s)
90	52	52	1CX7.1	10 unit(s)
125	52	67	1CX8.1	10 unit(s)
250	68	79	1CX9.1	10 unit(s)
625	90	113	1CXA.1	5 unit(s)
1,000	102	129	1CXC.1	5 unit(s)



Sample tub SpecTainer™

Simport. Material: PP, highly transparent, screw cap made of PE.

Shatterproof, stackable and leakproof.

The tamper-evident cap guarantees the safety of the contents: The security tab breaks when unscrewing the cap, indicating that the container has been opened. Temperature stability -90 to +100 °C. Contains no plasticisers, release agents and latex.

Purpose: For use as urine beaker.

In vitro diagnostic medical device according to IVDD

Vol. (ml)	Sterility	Ø ext. (mm)	Height (mm)	Lid colour	Art. No.	Pack Qty.
60	Non-sterile	48	62	yellow	AYT6.1	500 unit(s)
60	Sterile	48	62	turquoise	AYT9.1	500 unit(s)
90	Non-sterile	58	62	yellow	AYT7.1	400 unit(s)
90	Sterile	58	62	turquoise	AYX0.1	400 unit(s)
120	Non-sterile	58	78	yellow	AYT8.1	300 unit(s)
120	Sterile	58	78	turquoise	AYX1.1	300 unit(s)

Sample cans



Sample tub SecurTainer™ III

Simport. Material: PP transparent, screw cap made of PE.

Shatterproof, stackable and leakproof.

Purpose: For storage or transport of histological samples, urine, sputum and other liquids or solids. Temperature stability -90 to +100 °C.

In vitro diagnostic medical device according to IVDD

SecurTainer™

Vol. (ml)	Ø ext. (mm)	Outer lid Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
20	28	35	54	HYE7.1	500 unit(s)	
40	43	48	54	HYE8.1	500 unit(s)	
60	43	48	64	HYE9.1	500 unit(s)	
90	52	56	63	HYH0.1	400 unit(s)	
120	52	56	78	HYH1.1	300 unit(s)	



Accessories lid inserts

ROTH SELECTION. Material: PP.

For colour-coded, easy identification of:

- Reaction vials (Art. No. CK47.1 ff, CK57.1 ff and ATC4.1)
- Slide boxes with tamper-evident cap (Art. No. KL89.1)
- Sample boxes (Art. No. AYT6.1 ff)

Ø ext. (mm)	Colour	Designation	Art. No.	Pack Qty.	
11	blue	Lid inserts, blue	CK72.1	500 unit(s)	
11	yellow	Lid inserts, yellow	CK73.1	500 unit(s)	
11	red	Lid inserts, red	CK74.1	500 unit(s)	
11	white	Lid inserts, white	CK75.1	500 unit(s)	
11	green	Lid inserts, green	CK76.1	500 unit(s)	
11	assorted colours	Lid inserts, assorted colours	CK77.1	500 unit(s)	



S



Cryogenic boxes with screw closure

Nalgene. Material: PP, lid made of HDPE. **Sterile.**

Type 5005. For storing solid samples. For use in freezers and/or deep freezers. Pyrogen-free and non-cytotoxic.

In vitro diagnostic medical device according to IVDD

Vol. (ml)	Sterility	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
15	Sterile	33	47	PP19.1	75 unit(s)	114,15



S



Sample tub

Greiner Bio-One. Material: PS, transparent with white screw cap made of HDPE. **Sterile.**

Impervious to fluids. Temperature stability -20 to +60 °C. Suitable for universal use in the laboratory for non-human and non-medical samples.

Not a medical device / Not an IVD product

Vol. (ml)	Nominal volume (ml)	Sterility	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
50	60	Sterile	40	60	HYH7.1	300 unit(s)	



Sample cans



Sample tub

Greiner Bio-One. Material: PS transparent, screw cap made of metal. Impervious to fluids. Temperature stability -20 to +60 °C. Suitable for universal use in the laboratory for non-human and non-medical samples.

Not a medical device / Not an IVD product

Vol. (ml)	Nominal volume (ml)	Sterility	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
80	100	Non-sterile	45	77	HYK0.1	160 unit(s)	
150	175	Non-sterile	49	107	HYK1.1	120 unit(s)	
150	175	Sterile	49	107	HYK2.1	120 unit(s)	
260	290	Non-sterile	58	121	HYK3.1	48 unit(s)	



Glass jar DURAN®

DWK Life Sciences. Material: DURAN® glass. Autoclavable: Yes. High temperature and chemical resistance.

Delivery incl. lid.

Not a medical device / Not an IVD product

Vol. (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
75	60	40	CCL2.1	1 unit(s)	
175	80	50	CCL3.1	1 unit(s)	
325	100	60	CCL4.1	1 unit(s)	
1,000	150	80	CCL5.1	1 unit(s)	



Cryogenic boxes ROTILABO®

ROTH SELECTION. Material: Aluminium 99.5%, non-coated. Autoclavable: Yes.

With embedded rubber seal in the screw-on lid. Temperature stable up to -196 °C.

Not a medical device / Not an IVD product

Vol. (ml)	Nominal volume (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
2	3	14.5	23	1083.1	100 unit(s)	
10	15	30	29	1084.1	100 unit(s)	



Sample tub ROTILABO® fluoroplastics

ROTH SELECTION. Material: PTFE. Autoclavable: Yes.

Chemically inert. Temperature stability -200 to +260 °C. Thick-walled version, stackable.

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Ø Inner neck (mm)	Height (mm)	Art. No.	Pack Qty.	
15	34	26	34	LT85.1	1 unit(s)	
30	34	26	62	LT86.1	1 unit(s)	
60	60	47	46	LT87.1	1 unit(s)	
120	72	54	62	LT88.1	1 unit(s)	
240	72	60	100	LT89.1	1 unit(s)	
360	90	73	95	LT90.1	1 unit(s)	
480	90	73	125	LT91.1	1 unit(s)	
1,000	110	92	160	LT92.1	1 unit(s)	
2,200	134	110.6	260	LT93.1	1 unit(s)	

Sample cans



Sample tub ROTILABO® fluoroplastics

ROTH SELECTION. Material: PTFE. Autoclavable: Yes.

Chemically inert. Thick-walled type with conical base. Temperature stability -200 to +260 °C.

Delivery incl. screw closure.

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Ø Inner neck (mm)	Height (mm)	Art. No.	Pack Qty.	
25	33	25	64	LT82.1	1 unit(s)	
50	43	31	78	LT83.1	1 unit(s)	



Sample tub fluoroplastics

Saint-Gobain. Material: PFA, highly transparent. Autoclavable: Yes.

With labelling area and seal, leak-proof. Cryo- and microwave-safe.

- Temperature stability -196 to +250 °C
- Inert to virtually all chemicals, including hydrofluoric and nitric acid
- High purity, low metal ion content, ideal for use in trace analysis and for the storage of highly pure substances
- No memory effect
- Hydrophobic and anti-adhesive, extremely smooth surfaces allow easy cleaning

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
60	50.8	46.8	YY50.1	1 unit(s)	
120	66	53.2	YY51.1	1 unit(s)	
240	75.4	96	YY52.1	1 unit(s)	
360	95	90.4	YY53.1	1 unit(s)	
480	65.2	94.3	YY54.1	1 unit(s)	
1,000	106.7	150	YY55.1	1 unit(s)	



Sample tub fluoroplastics

Vitlab. Material: PFA, highly transparent. Autoclavable: Yes.

With screw cap.

- For sampling, transportation and storage
- Chemically inert to nearly all chemicals
- Hydrophobic and anti-adhesive, extremely smooth surfaces allow easy cleaning
- High purity, ideal for use in trace analysis and for storage of reference solutions
- No memory effect
- Temperature stability -200 to +250 °C

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Thread	Art. No.	Pack Qty.
30	38	54	GL 40	YA21.1	1 unit(s)
60	38	90	GL 40	YA22.1	1 unit(s)
90	54	62	GL 56	YA23.1	1 unit(s)
180	54	112	GL 56	YA24.1	1 unit(s)

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47