

Чаши, лотки, поддоны для проб, образцов CARL ROTH ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47



Sampling bowls



A 121 °C

Sample tray ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10, high-gloss. Autoclavable: Yes.

Suitable for foodstuffs.

Delivery without lid, please order separately.

Not a medical device / Not an IVD product

Vol. (ml)	Ø Internal (mm)	Height (mm)	Art. No.	Pack Qty.
250	120	35	YH20.1	1 unit(s)
400	130	40	YH21.1	1 unit(s)
650	150	45	YH22.1	1 unit(s)
850	170	45	YH23.1	1 unit(s)
1,250	190	50	YH24.1	1 unit(s)
1,400	190	55	YH25.1	1 unit(s)
1,600	210	55	YH26.1	1 unit(s)
2,000	225	55	YH27.1	1 unit(s)

A 121 °C

Accessories lid for ROTILABO® sample trays

ROTH SELECTION. Material: Stainless steel 18/10, high-gloss. Autoclavable: Yes.

Not a medical device / Not an IVD product

Suitable for	Art. No.	Pack Qty.
Sample tray Order No. YH20.1	YL03.1	1 unit(s)
Sample tray Order No. YH21.1	YL04.1	1 unit(s)
Sample tray Order No. YH22.1	YL05.1	1 unit(s)
Sample tray Order No. YH23.1	YL06.1	1 unit(s)
Sample tray Order No. YH24.1	YL07.1	1 unit(s)
Sample tray Order No. YH25.1	YL08.1	1 unit(s)
Sample tray Order No. YH26.1	YL09.1	1 unit(s)
Sample tray Order No. YH27.1	YL10.1	1 unit(s)



A 121 °C

Tablet pans

Contacto Bander. Material: Stainless steel 18/10, high-gloss. Autoclavable: Yes.

Suitable for foodstuffs.

Not a medical device / Not an IVD product

Flat shape

Ø (mm)	Height (mm)	Art. No.	Pack Qty.
120	20	YX74.1	1 unit(s)
140	25	YX75.1	1 unit(s)

Tall

Ø (mm)	Height (mm)	Art. No.	Pack Qty.
65	20	YX76.1	1 unit(s)
100	30	YX77.1	1 unit(s)



Sample tray ROTILABO® with lid

ROTH SELECTION. Material: Aluminium.

Pressed from one piece. Non-permeable. The cardboard lids that are aluminium-coated on one side are held in place by bending the tray edge. The white top side can be marked or labelled.

Not a medical device / Not an IVD product

Vol. (ml)	Upper length (mm)	Upper width (mm)	Height (mm)	Art. No.	Pack Qty.
250	120	95	35	0522.1	100 unit(s)
500	140	115	40	0906.1	100 unit(s)

Sampling bowls



Sample tray ROTILABO®

ROTH SELECTION. Material: Aluminium.

Free of oil residues and coatings. Highly versatile as weighing or drying dishes, as dispensing or storage containers or for covering vessels. Stackable.

Not a medical device / Not an IVD product

Vol. (ml)	Ø Internal (mm)	Height (mm)	Art. No.	Pack Qty.	
20	44	12	PP65.1	100 unit(s)	
60	63	17	PP66.1	100 unit(s)	
75	62	30	PP67.1	100 unit(s)	
200	110	35	PP68.1	50 unit(s)	
500	181	25	PP69.1	50 unit(s)	



Sample tray ROTILABO®

ROTH SELECTION. Material: Aluminium.

Not a medical device / Not an IVD product

Vol. (ml)	Ø Outer upper (mm)	Ø Internal (mm)	Height (mm)	Art. No.	Pack Qty.	
28	64	55	13	5494.1	100 unit(s)	
125	96	85	26	5496.1	100 unit(s)	



Dappenglas

Hecht. Material: Soda-lime glass.

For mixing small quantities of substances. Usable on both sides, volume approx. 1.5 ml or 3 ml.

Not a medical device / Not an IVD product

Ø (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
31	29	colourless	CKY4.1	1 unit(s)	
31	29	green	CKY5.1	1 unit(s)	
31	29	pink	CKY6.1	1 unit(s)	
31	29	brown	CKY7.1	1 unit(s)	



①



②

Specimen jars ROTILABO®

ROTH SELECTION. Material: Borosilicate glass 3.3.

Not a medical device / Not an IVD product

(1) With knob

Ø ext. (mm)	Outer lid Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
85	95	85	NT38.1	1 unit(s)	
100	110	100	NT39.1	1 unit(s)	
120	130	120	NT40.1	1 unit(s)	
150	162	150	NT41.1	1 unit(s)	
200	214	200	NT42.1	1 unit(s)	
250	264	250	NT43.1	1 unit(s)	
300	314	300	NT44.1	1 unit(s)	

(2) With recessed knob

Ø ext. (mm)	Outer lid Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
65	60	65	NT46.1	1 unit(s)	
85	79	85	NT47.1	1 unit(s)	
100	94	100	NT48.1	1 unit(s)	
100	94	500	NT49.1	1 unit(s)	
120	103	120	NT50.1	1 unit(s)	
120	103	500	NT51.1	1 unit(s)	
150	130	150	NT52.1	1 unit(s)	
150	130	500	NT53.1	1 unit(s)	



Glass jar DURAN®

DWK Life Sciences. Material: DURAN® glass. Autoclavable: Yes. High temperature and chemical resistance.

Delivery incl. lid.

Not a medical device / Not an IVD product

Vol. (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
75	60	40	CCL2.1	1 unit(s)	
175	80	50	CCL3.1	1 unit(s)	
325	100	60	CCL4.1	1 unit(s)	
1,000	150	80	CCL5.1	1 unit(s)	



Sampling bowls



A 121 °C

Cotton and swab containers

Contacto Bander. Material: Stainless steel 18/10. Autoclavable: Yes. Suitable for foodstuffs.

Delivery without lid, please order separately.

Not a medical device / Not an IVD product

Vol. (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.
150	60	60	YT80.1	1 unit(s)
250	85	50	YT81.1	1 unit(s)
350	100	50	YT82.1	1 unit(s)
400	85	85	YT83.1	1 unit(s)
600	85	125	YT84.1	1 unit(s)
750	100	100	YT85.1	1 unit(s)
1,750	150	100	YT86.1	1 unit(s)
2,500	150	150	YT87.1	1 unit(s)
3,000	200	100	YT88.1	1 unit(s)

A 121 °C

Accessories lids for cotton wool and swab containers

Contacto Bander. Material: Stainless steel 18/10. Autoclavable: Yes.

Not a medical device / Not an IVD product

Suitable for	Art. No.	Pack Qty.
Container-Ø 60 mm	YT89.1	1 unit(s)
Container-Ø 85 mm	YT90.1	1 unit(s)
Container-Ø 100 mm	YT91.1	1 unit(s)
Container-Ø 150 mm	YT92.1	1 unit(s)
Container-Ø 200 mm	YT93.1	1 unit(s)



A 121 °C

Crystallization dishes ROTILABO®

ROTH SELECTION. Material: Borosilicate glass 3.3. Autoclavable: Yes.

With spout

Acc. to DIN 12338.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.
20	40	25	AYE1.1	10 unit(s)
40	50	30	AYE2.1	10 unit(s)
60	60	35	AYE3.1	10 unit(s)
100	70	40	AYE4.1	10 unit(s)
150	80	45	AYE5.1	10 unit(s)
300	95	55	AYE6.1	10 unit(s)
500	115	65	AYE7.1	10 unit(s)
900	140	75	AYE8.1	10 unit(s)
2,000	190	90	AYE9.1	1 unit(s)
3,500	230	100	AYH0.1	1 unit(s)

Without spout

Acc. to DIN 12337.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.
20	40	25	AYH1.1	10 unit(s)
40	50	30	AYH2.1	10 unit(s)
60	60	35	AYH3.1	10 unit(s)
100	70	40	AYH4.1	10 unit(s)
150	80	45	AYH5.1	10 unit(s)
300	95	55	AYH6.1	10 unit(s)
500	115	65	AYH7.1	10 unit(s)
900	140	75	AYH8.1	10 unit(s)
2,000	190	90	AYH9.1	1 unit(s)
3,500	230	100	AYK0.1	1 unit(s)

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47