

Камеры, штативы для окрашивания мазков CARL ROTH StainTray, Schiefferdecker, Wash-N-Dry, Columbia, ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тургут (3462)77-98-35
Тверь (4822)63-31-35


Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47

Staining racks and chambers



Staining rack

Marienfeld.

For custom adaptation to a variety of dishes and bowls. With 2 fixing clips and 2 handles, each extendible by 115 mm. For tilting and fixing in a horizontal position. When used in combination with staining tray CEX6.1, a locking bolt holds the tray securely in the tilted position. The tray can also be fixed in the horizontal position using the second locking bolt. Total length: Approx. 380 mm.

Delivery without staining tray.

Version	Art. No.	Pack Qty.	
Staining rack	CEX5.1	1 unit(s)	

Accessories staining tray

Marienfeld. Material: PVC.

Staining tray inner dimensions: L 380 x W 145 x H 75 mm.

Designation	Version	Art. No.	Pack Qty.	
Staining tray, inner dimensions: L 380 x W 145 x H 75 mm	Staining tray	CEX6.1	1 unit(s)	



Staining chamber

Nalgene. Material: PMP, lid made of PP copolymer, drain cap made of LDPE.

For staining, fixing, bleaching membranes, electrophoresis gels, etc. Chamber with handles on both sides for easy transportation. Capacity with gel added under shaking bath conditions: 500 to 750 ml. Temperature stability: -70 to +70 °C.

Not UV-transparent.

L x W x H	Art. No.	Pack Qty.	
12.5 x 12.5 x 50	H591.1	1 unit(s)	



Staining chamber StainTray™

Simport. Material: ABS.

Suitable as a staining chamber or a humidity chamber in immune histology, haematology, cytology and microbiology.

- Slides securely adhere to the rubber strip
- Drain with plug for simple emptying of the chamber
- Stackable

Black lid

Material: Lid made of ABS. Temperature stability: -20 to +80 °C.

Pic.	Suitable for	L x W x H (mm)	Art. No.	Pack Qty.	
(1)	20 microscope slides	380 x 240 x 45	HA51.1	1 unit(s)	
(2)	10 microscope slides	240 x 240 x 45	HL98.1	1 unit(s)	

Transparent lid

Material: Lid made of PETG. Temperature stability: -20 to +60 °C.

Pic.	Suitable for	L x W x H (mm)	Art. No.	Pack Qty.	
(1)	20 microscope slides	380 x 240 x 45	HA50.1	1 unit(s)	
(2)	10 microscope slides	240 x 240 x 45	HL97.1	1 unit(s)	



Disposable staining chambers

ROTH SELECTION. Material: PP, PE.

For staining, washing and drying microscope slides with dimensions of 76 x 26 mm.

- Chamber volume: 38 ml
- Recessed handles for easy transportation
- Two drain openings
- Dark lid for light-sensitive applications

Delivery incl. four staining chambers and one lid.

Suitable for	L x W x H (mm)	Art. No.	Pack Qty.	
8 microscope slides	305 x 127 x 30	TP31.1	1 set	



Staining racks and chambers



Staining trough ROTILABO®

Simport. Material: POM.

Resistant to most staining reagents, but irreversible discolouration may occur.

Delivery incl. connection clips, without microscope slide holder.

Vol. (ml)	Colour	L x W x H (mm)	Art. No.	Pack Qty.	
80	blue	64 x 76 x 92	HA44.1	1 unit(s)	
80	yellow	64 x 76 x 92	HA45.1	1 unit(s)	
80	green	64 x 76 x 92	HA46.1	1 unit(s)	
80	pink	64 x 76 x 92	HA47.1	1 unit(s)	
80	white	64 x 76 x 92	HA48.1	1 unit(s)	



Accessories slide holder ROTILABO®

Simport. Material: POM.

For 76 x 26 mm microscope slides. Vertical storage of the microscope slides means the labelling areas do not come into contact with staining solution.

Suitable for	Colour	L x W x H (mm)	Art. No.	Pack Qty.	
12 microscope slides	grey	60 x 64 x 97	HA49.1	1 unit(s)	



Staining trough

Kartell. Material: PMP.

With two lids: One which seals tightly and the other with a small opening for attaching the insert during staining.

Microwave safe.

Delivery without inserts.

L x W x H (mm)	Art. No.	Pack Qty.	
101 x 81 x 70	2290.2	4 unit(s)	



Accessories insert for staining trough

Kartell. Material: PP.

For microscope slides 76 x 26 mm.

Designation	Suitable for	L x W x H (mm)	Art. No.	Pack Qty.	
Inserts for 20 microscope slides with dimensions of 76 x 26 mm,	20 microscope slides	86 x 70 x 21	2291.2	2 unit(s)	

Staining racks and chambers



Staining dish

Kartell. Material: PMP.

For microscope slides with dimensions of 76 x 26 mm. Microwave-safe.

Delivery incl. lid.

Suitable for	L x W x H (mm)	Art. No.	Pack Qty.	
8 microscope slides	58 x 54 x 86	2292.2	4 unit(s)	



Staining trough acc. to Schiefferdecker

Kartell. Material: PMP.

For microscope slides 76 x 26 mm. Microwave safe.

Suitable for	L x W x H (mm)	Art. No.	Pack Qty.	
20 microscope slides	89 x 70 x 50	XC73.1	4 unit(s)	



A 121 °C

Staining trough acc. to Coplin

ROTH SELECTION. Material: PP. Autoclavable: Yes.

For microscope slides with dimensions of 76 x 26 mm.

Opening Ø of 53 mm for simple placement.

Delivery incl. screw-on lid.

Suitable for	Ø x H (mm)	Art. No.	Pack Qty.	
10 microscope slides	53 x 114	YC09.1	1 unit(s)	



Staining trough acc. to Coplin

DWK Life Sciences. Material: Soda-lime glass.

For microscope slides 76 x 26 mm.

Suitable for	Ø x H (mm)	Art. No.	Pack Qty.	
10 microscope slides	66 x 108	H548.1	1 unit(s)	



Staining dish acc. to Coplin with screw cap

Mariefeld. Material: Soda-lime glass, screw cap made of plastic (urea 131.5), with PE foam insert.

For microscope slides 76 x 26 mm.

Delivery incl. screw cap.

Suitable for	Height without closure (mm)	Art. No.	Pack Qty.	
10 microscope slides	108	ETE2.1	6 unit(s)	
10 microscope slides	85	ETE3.1	6 unit(s)	



Staining racks and chambers



Staining dish acc. to Coplin

DWK Life Sciences. Material: Soda-lime glass, screw top made of PP. For microscope slides 76 x 26 mm. Dimensions: L 26 x W 26 x D 90 mm.

Delivery incl. screw top.

Suitable for	Vol. (ml)	L x W x H (mm)	Art. No.	Pack Qty.	
10 microscope slides	60	26 x 26 x 90	HLA7.1	1 unit(s)	



A 121 °C

Microscope slide holder Wash-N-Dry™

ROTH SELECTION. Material: PP. Autoclavable: Yes.

For staining, washing and drying microscope slides with dimensions of 76 x 26 mm. Ideal for tall 150 ml beakers. Distance to floor, 10 mm below the microscope slide for the use of magnetic stir bars.

Suitable for	L x W x H (mm)	Art. No.	Pack Qty.	
12 microscope slides	31.5 x 38.4 x 110	AXC6.1	1 unit(s)	



Staining dish acc. to Columbia

DWK Life Sciences. Material: Soda-lime glass, screw top made of PP, with PTFE-coated PE insert.

For coverslips, W 17–23 mm, H 30 mm Longer coverslips can be inserted with the screw top removed.

Delivery incl. screw top.

Suitable for	Art. No.	Pack Qty.	
4 coverslips	HLA6.1	1 unit(s)	



A 121 °C

Coverslip rack Wash-N-Dry™

ROTH SELECTION. Material: PP. Autoclavable: Yes.

For square or round coverslips with dimensions of 18 to 25 mm. With notches on the underside for fixing the coverslips, e.g. with elastic bands.

Colour	Slots	L x W x H (mm)	Art. No.	Pack Qty.	
white	10	39 x 31 x 24	TH51.1	1 unit(s)	

Staining racks and chambers



Staining trough acc. to Hellendahl

ROTH SELECTION. Material: Soda-lime glass.
For microscope slides 76 x 26 mm.

Without expansion

Pic.	Suitable for	L x W x H (mm)	Art. No.	Pack Qty.
(1)	8 microscope slides	93 x 43 x 97	H549.1	1 unit(s)

With expansion

Autoclavable: Yes.

Pic.	Suitable for	L x W x H (mm)	Art. No.	Pack Qty.
(2)	8 microscope slides	63 x 58 x 98	H550.1	1 unit(s)



A 121 °C

Staining trough porcelain

ROTH SELECTION. Material: Porcelain, glazed. Autoclavable: Yes.
For the protection of light-sensitive staining solutions.

- For microscope slides with dimensions of 76 x 26 mm
- Lid dimensions: L 75 x W 48 x H 5 mm
- Temperature stability: to +1000 °C
- Microwave-safe

Delivery incl. lid.

Suitable for	L x W x H (mm)	Art. No.	Pack Qty.
6 microscope slides	70 x 40 x 86	XX63.1	2 unit(s)



Staining trough acc. to Schiefferdecker

ROTH SELECTION. Material: Soda-lime glass.
For microscope slides 76 x 26 mm.

Suitable for	L x W x H (mm)	Art. No.	Pack Qty.
20 microscope slides	89 x 71 x 46	H551.1	1 unit(s)



Staining trough

ROTH SELECTION.

Material: Soda-lime glass.

Delivery incl. lid.

L x W x H (mm)	Art. No.	Pack Qty.
105 x 85 x 70	H554.1	1 unit(s)



Accessories glass insert for staining trough

ROTH SELECTION. Material: Soda-lime glass.

For microscope slides with dimensions of 76 x 26 mm or any width up to 52 mm.


Delivery without wire hanger.

Suitable for	L x W x H (mm)	Art. No.	Pack Qty.
10 microscope slides	85 x 60 x 45	H552.1	1 unit(s)

Accessories wire hanger for glass staining rack

ROTH SELECTION. Material: Stainless steel.

Version	Art. No.	Pack Qty.
Stainless steel	H553.1	1 unit(s)



Accessories metal insert for staining trough

Mariefeld. Material: Stainless steel.

For microscope slides with dimensions of 76 x 26 mm or any width up to 52 mm.

Suitable for	L x W x H (mm)	Art. No.	Pack Qty.
10 microscope slides	80 x 55 x 53	CET8.1	1 unit(s)



Staining racks and chambers



Staining trough

DWK Life Sciences. Material: Soda-lime glass, stainless steel rack.
Staining system for microscope slides with dimensions of 76 x 26 mm, 76 x 38 mm and 76 x 52 mm.
Rack with folding handle for closing the dish and inserting and removing the microscope slides.
Inner dimensions at base: L 185 x W 88 x D 78 mm.

Delivery incl. lid and rack.

Suitable for	Art. No.	Pack Qty.	
50 microscope slides	HLA2.1	1 set	



Staining trough acc. to Hausser/Gedigk

ROTH SELECTION. Material: Soda-lime glass.

Delivery without lid and staining insert.

L x W x H (mm)	Art. No.	Pack Qty.	
130 x 93 x 121	T725.1	1 unit(s)	

Accessories for Hausser/Gedigk staining trough

ROTH SELECTION.

For microscope slides 76 x 26 mm.

Lid

Designation	L x W x H (mm)	Art. No.	Pack Qty.	
Lid, L 130 x W 93 x H 3 mm	130 x 93 x 3	T726.1	1 unit(s)	

Staining insert

Designation	Suitable for	L x W x H (mm)	Art. No.	Pack Qty.	
Staining insert made of stainless steel, for microscope slides with dimensions of 76 x 26 mm, L 110 x W 75 x H 95 mm	60 microscope slides	110 x 75 x 95	T727.1	1 unit(s)	



Preparation box

DWK Life Sciences. Material: DURAN® glass.

With ground glass plate. Versatile for use as a slide storage box, TLC separation chamber, etc.

Vol. (ml)	L x W x H (mm)	Art. No.	Pack Qty.	
250	60 x 50 x 100	T067.1	1 unit(s)	
550	100 x 50 x 120	T068.1	1 unit(s)	
800	130 x 50 x 130	T069.1	1 unit(s)	
1,000	150 x 50 x 150	T070.1	1 unit(s)	
1,200	120 x 60 x 180	T071.1	1 unit(s)	
4,000	210 x 100 x 210	T072.1	1 unit(s)	
8,500	250 x 140 x 250	T073.1	1 unit(s)	



Glass boxes

ROTH SELECTION. Material: Soda-lime glass.

Ground rim.

Vol. (l)	L x W x H (mm)	Art. No.	Pack Qty.	
0.8	100 x 100 x 100	NY95.1	1 unit(s)	
1.8	100 x 150 x 150	NY96.1	1 unit(s)	
5.7	220 x 300 x 150	NY97.1	1 unit(s)	
8.7	180 x 250 x 220	NY98.1	1 unit(s)	
19.9	230 x 360 x 260	NY99.1	1 unit(s)	



A121 °C

Staining dish silicone

Hammacher. Material: Silicone, soft. Autoclavable: Yes.

Temperature stability: -55 to +200 °C.

Vol. (ml)	Art. No.	Pack Qty.	
160	HCY2.1	1 unit(s)	

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47