

Подставки, штативы для флаконов, пробирок CARL ROTH ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тургут (3462)77-98-35
Тверь (4822)63-31-35


Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47

Floating racks



Floating stands ROTILABO®

ROTH SELECTION. Material: HDPE (foamed).

For bringing reaction vials and centrifuge tubes up to the right temperature in water baths. With handles for insertion and removal. When loaded, the stands float on the surface of the water, while the inserted tubes do not float.

For 0.2–2 ml reaction vials

Pic.	No. of slots	Ø x H (mm)	Colour	Art. No.	Pack Qty.	
(1)	18	82 x 60	blue	LP92.1	5 unit(s)	

For 1.5–2 ml reaction vials

Pic.	No. of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
(2)	24	150 x 82 x 69	yellow	P199.1	5 unit(s)	

For 15 ml centrifuge tubes

Pic.	No. of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
(3)	8	130 x 59 x 65	blue	P200.1	5 unit(s)	

For 50 ml centrifuge tubes

Pic.	No. of slots	L x W x H (mm)	Colour	Art. No.	Pack Qty.	
(4)	4	105 x 105 x 65	green	P201.1	5 unit(s)	



Floating stands

Nalgene. Material: PP. Autoclavable: Yes.

For bringing reaction vials up to the right temperature in water bath vessels. With handle for insertion and removal.

- Racks float on the surface of the water when fully loaded
- Reaction vials are fixed in position in the rack and therefore cannot float
- Suitable for 0.5–2 ml reaction vials
- Temperature stability -70 to +100 °C

Version	No. of slots	Art. No.	Pack Qty.	
round	8	E072.1	1 unit(s)	
round	20	E073.1	1 unit(s)	
square	16	E075.1	1 unit(s)	



Floating stands Cryogenic

ROTH SELECTION. Material: PE (foamed).

Ideal for bringing samples in PCR vessels to the right temperature. The loaded stands are still fully capable of floating; the tubes are fitted securely, meaning that it is not possible for them to float.

- Suitable for: 0.2 ml PCR vessels
- Temperature stability -80 to +80 °C

Round

No. of slots	Ø x H (mm)	Art. No.	Pack Qty.	
10	70 x 8.5	LH39.1	10 unit(s)	

Square

Suitable for PCR single vessels and PCR 8-strips.

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
32	4 x 8	80 x 80 x 8.5	LH40.1	10 unit(s)	




Stands for vials, tubes

Storage box ROTILABO®

ROTH SELECTION. Material: PP. Autoclavable: Yes.

- With transparent lid
- Stackable
- Temperature stability -80 to +100 °C



A121 °C

(1) For 1.5 ml autosampler vials (ND8–ND11)

With alphanumeric labelling on all four edges and the base cavities.

Version	Compartment size (mm)	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	12 x 12	130	130	45	pink	TY99.1	1 unit(s)	
Slip lid	12 x 12	130	130	45	yellow	XA00.1	1 unit(s)	
Slip lid	12 x 12	130	130	45	orange	XA01.1	1 unit(s)	
Slip lid	12 x 12	130	130	45	blue	XA02.1	1 unit(s)	
Slip lid	12 x 12	130	130	45	green	XA03.1	1 unit(s)	
Slip lid	12 x 12	130	130	45	colourless	XA04.1	1 unit(s)	

(2) For 4 ml autosampler vials (ND13)

With alphanumeric labelling on all four edges and the base cavities.

Version	Compartment size (mm)	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	Ø 15	130	130	52	red	XA05.1	1 unit(s)	

(3) For 20 ml EPA screw-top bottles

Version	Compartment size (mm)	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	30 x 30	130	130	70	violet	XA07.1	1 unit(s)	

(4) For 30 and 40 ml EPA screw-top bottles

Version	Compartment size (mm)	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	Ø 30	130	130	105	violet	XA08.1	1 unit(s)	

(3) For 5, 10 and 20 ml headspace bottles

Version	Compartment size (mm)	Length (mm)	Width (mm)	Height (mm)	Colour	Art. No.	Pack Qty.	
Slip lid	25 x 25	130	130	102	blue	XA06.1	1 unit(s)	

Storage box for sample vials

ROTH SELECTION. Material: PP. Autoclavable: Yes.

- Suitable for ND8–ND11 autosampler vials



A121 °C

Version	Compartment size (mm)	Slots	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
Slip lid	12 x 12	96 (8 x 12)	172.5	116.5	40.5	A618.1	1 unit(s)	

Stands for vials, tubes



Sample stands ROTILABO®

ROTH SELECTION. Material: PMMA (Acrylic), transparent.
With silicon feet for stability and stackability.

For autosampler vial micro inserts

Also suitable for tubes with a conical base.

- Compartment size (hole): Ø 8,5 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
50	5 x 10	173 x 95 x 20	TY95.1	1 unit(s)	

For ND8 to ND11 autosampler vials

Also suitable for tubes with a conical base.

- Compartment size (hole): Ø 12 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
50	5 x 10	173 x 95 x 20	TY96.1	1 unit(s)	

For ND13 autosampler vials

- Compartment size (hole): Ø 15,1 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
40	5 x 8	176 x 116 x 20	TY97.1	1 unit(s)	

For ND18, ND20 autosampler vials (up to 20 ml)

- Compartment size (hole): Ø 24 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
25	5 x 5	160 x 160 x 30	TY98.1	1 unit(s)	



Sample stands wood

ROTH SELECTION. Material: Wood.

- Compartment size (hole): Ø 11 mm
- Suitable for: Test tubes Ø 10 mm

L x W x H (mm)	Art. No.	Pack Qty.	
257 x 47 x 24	0403.1	1 unit(s)	



A 121 °C



Sample stands

ROTH SELECTION. Material: PP. Autoclavable: Yes.

For incubation of small test tubes, serum / RIA test tubes and T3/T4 reaction vials.

- Suitable for most water baths
- Height of second level: 41 mm
- Temperature stability -80 to +121 °C

For tube Ø 10–13 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
72	238 x 111 x 89	YC06.1	1 unit(s)	

For tube Ø 13–16 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
40	254 x 111 x 89	YC07.1	1 unit(s)	

For tube Ø 16–20 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
40	254 x 111 x 89	YC08.1	1 unit(s)	



A 121 °C



Sample stands ROTILABO® PTFE

ROTH SELECTION. Material: PTFE. Autoclavable: Yes.

- Chemically inert
- Temperature stability -200 to +250 °C

For tube Ø 6–8 mm

- Compartment size (hole): Ø 8 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
27	3 x 9	180 x 60 x 45	NY05.1	1 unit(s)	

For tube Ø 16–18 mm

- Compartment size (hole): Ø 19 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
10	2 x 5	180 x 60 x 70	NY07.1	1 unit(s)	



Stands for vials, tubes

For tube Ø 25–30 mm

- Compartment size (hole): Ø 30 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
4	1 x 4	180 x 60 x 80	NY08.1	1 unit(s)	



A 121 °C

Sample stands for scintillation flasks

ROTH SELECTION. Material: PP. Autoclavable: Yes.

Perfect for storage, transportation or as a rack on the laboratory bench.

- Suitable for: Scintillation vials up to Ø 28 mm
- Compartment size (hole): Ø 30 mm
- Height to intermediate floor 44 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
24	3 x 8	298 x 127 x 95	K993.1	1 unit(s)	



Sample stands acrylic

ROTH SELECTION. Material: PMMA (Acrylic), transparent.

- Suitable for: Sample vials Ø 20–22 mm
- Compartment size (hole): Ø 24 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
30	3 x 10	330 x 110 x 22	YC05.1	1 unit(s)	



A 121 °C

Sample stands

ROTH SELECTION. Material: Steel, epoxy-coated. Autoclavable: Yes.
Temperature stability to +121 °C.

For vessel Ø 23 mm

- Compartment size: 24 x 24 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
12	2 x 6	164 x 64 x 32	AXE5.1	1 unit(s)	
50	5 x 10	271 x 145 x 32	AXE6.1	1 unit(s)	

For vessel Ø 32 mm

- Compartment size: 33 x 33 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
12	2 x 6	212 x 80 x 60	AXE7.1	1 unit(s)	
50	5 x 10	355 x 189 x 60	AXE8.1	1 unit(s)	

For vessel Ø 56 mm

- Compartment size: 58 x 58 mm

No. of slots	Arrangement of slots	L x W x H (mm)	Art. No.	Pack Qty.	
24	4 x 6	393 x 262 x 50	AXE9.1	1 unit(s)	



A 121 °C

Sample stands for NMR tubes

ROTH SELECTION. Material: PP. Autoclavable: Yes.

With two handles on the side for easy transportation. Stability is provided by hemispherical holes in the lower level.

- Suitable for NMR tubes Ø 5 mm

No. of slots	L x W x H (mm)	Art. No.	Pack Qty.	
72	210 x 110 x 220	YA94.1	1 unit(s)	

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Россия (495)268-04-70

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Казахстан (772)734-952-31

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Киргизия (996)312-96-26-47

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93