

Пробирки, виалы для испытаний CARL ROTH ROTILABO

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47

Test vials



A 121 °C



Test tubes

Boettger. Material: PP, transparent. Autoclavable: Yes.

Delivery without screw caps, please order separately.

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
40	30	85	PT32.1	500 unit(s)	
60	30	126	PT33.1	300 unit(s)	
120	44	115	PT34.1	250 unit(s)	

A 121 °C

Accessories screw caps

Boettger. Material: PP. Autoclavable: Yes.

Not a medical device / Not an IVD product

Version	Suitable for	Art. No.	Pack Qty.	
Standard	PT29.1, PT30.1, PT32.1 and PT33.1	PT35.1	500 unit(s)	
With stand for tube-Ø 16 mm	PT29.1, PT30.1, PT32.1 and PT33.1	PT36.1	300 unit(s)	
With stand for tube-Ø 25 mm	PT31.1 and PT34.1	PT38.1	250 unit(s)	
With attachment Ø inside 16 mm	PT29.1, PT30.1, PT32.1 and PT33.1	PT37.1	300 unit(s)	



Culture vials with screw cap

DWK Life Sciences.

With DIN thread.

Delivery incl. screw caps.

Not a medical device / Not an IVD product

Soda-lime glass

Material: Soda-lime glass screw cap made of PP with white rubber seal. Autoclavable: Limited. Note: Can cause glass corrosion when autoclaved several times. Heat resistant up to approx. 100 °C.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Thread	Art. No.	Pack Qty.	
6	12	100	GL 14	K251.1	100 unit(s)	
12	16	100	GL 18	K252.1	100 unit(s)	
20	16	160	GL 18	K253.1	100 unit(s)	
30	18	180	GL 18	K254.1	100 unit(s)	

DURAN® glass

Material: DURAN® glass screw caps made of PBT with PTFE-coated seals. Autoclavable: Yes. Hot air sterilization to max. 180 °C.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Thread	Art. No.	Pack Qty.	
6	12	100	GL 14	K248.1	50 unit(s)	
9	13	100	GL 14	PY94.1	50 unit(s)	
20	16	160	GL 18	K249.1	50 unit(s)	
20	16	150	GL 18	PY95.1	50 unit(s)	
30	18	180	GL 18	K250.1	50 unit(s)	
34	20	150	GL 18	PY96.1	50 unit(s)	



Test vials

Test tubes made of fluoroplastic

Vitlab.

Highly transparent. Optionally with screw cap GL 25, available with or without stopper.



Not a medical device / Not an IVD product

With screw cap made of PFA

Material: PFA. Autoclavable: Yes.

Version	Vol. (ml)	Ø ext. (mm)	Height (mm)	Thread	Art. No.	Pack Qty.	
Without ring marking	15	22	110	GL 25	YA28.1	1 unit(s)	
With ring marking at 10 ml	15	22	110	GL 25	YA27.1	1 unit(s)	

With stopper made of PE

Material: PFA, PE.

Version	Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
With ring marking at 10 ml	12	16	110	YA25.1	1 unit(s)	

Without stopper

Material: PFA. Autoclavable: Yes.

Version	Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
Without ring marking	12	16	110	YA26.1	1 unit(s)	

Test tubes made of fluoroplastic

Saint-Gobain. Material: PFA, Highly transparent. Autoclavable: Yes.

- Temperature-stable from -196 to +250 °C
- Inert to virtually all chemicals, including hydrofluoric and nitric acid
- High purity, low metal ion content, ideal for use in trace analysis and for the storage of highly pure substances
- No memory effect
- Hydrophobic and anti-adhesive, extremely smooth surfaces allow easy cleaning



Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Ø Internal (mm)	Height (mm)	Art. No.	Pack Qty.	
5	10	8	132	YY79.1	1 unit(s)	
15	16	14	162	YY81.1	1 unit(s)	
20	16	14	175	YY82.1	1 unit(s)	
25	19	16	185	YY83.1	1 unit(s)	
50	22	19	235	YY84.1	1 unit(s)	
60	22	19	245	YY85.1	1 unit(s)	

Test vials



Test tubes with tamper-evident seal

Camlab. Material: PP. Autoclavable: Yes.

Temperature-stable from -90 to +121 °C. Cap dimensions: Ø 24 mm, height 2 mm.

- Leak-proof tamper-evident seal
- Free-standing, conical base for optimal sample recovery
- Frosted label field
- Graduated in 1 ml increments

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Lid colour	Art. No.	Pack Qty.
5	16	57	pink	X327.1	250 unit(s)
5	16	57	blue	X328.1	250 unit(s)
5	16	57	white	X329.1	250 unit(s)
10	15	94	pink	X330.1	250 unit(s)
10	15	94	blue	X331.1	250 unit(s)
10	15	94	white	X332.1	250 unit(s)



Culture vials with screw cap

DWK Life Sciences. Material: Borosilicate glass acc. to USP type I and ASTM E 438 type I class B; seals made of black phenol resin with seals made of PTFE-coated rubber.

With round bottom, suitable for use in shakers and tilt/roller mixers.

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Thread	Art. No.	Pack Qty.
11	16	100	15-415	ETK4.1	144 unit(s)
14	16	125	15-415	ETK5.1	144 unit(s)
27	20	150	18-415	ETK6.1	144 unit(s)



Test tube ROTILABO®

ROTH SELECTION. Autoclavable: Yes.

Good temperature stability.

Not a medical device / Not an IVD product

With beaded rim, thick-walled

Material: Borosilicate glass 3.3.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Art. No.	Pack Qty.
4	10	75	0,8-1,0	HA87.1	100 unit(s)
5	10	100	0,8-1,0	HA88.1	100 unit(s)
8	12	100	0,8-1,0	HA90.1	100 unit(s)
16	14	130	0,8-1,0	HA91.1	100 unit(s)
17	16	130	1,0-1,2	HA92.1	100 unit(s)
21	16	160	1,0-1,2	HA93.1	100 unit(s)
32	18	180	1,0-1,2	HA94.1	100 unit(s)
34	20	150	1,0-1,2	HA95.1	100 unit(s)
40	20	180	1,0-1,2	HA96.1	100 unit(s)
55	25	150	1,0-1,2	HA97.1	50 unit(s)
70	25	200	1,0-1,2	HA98.1	50 unit(s)
100	30	200	1,0-1,4	HA99.1	50 unit(s)

With beaded rim, thin-walled

Material: Borosilicate glass 5.1.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Art. No.	Pack Qty.
2	8	70	0,4-0,5	HC00.1	100 unit(s)
4	10	75	0,4-0,5	HC01.1	100 unit(s)
6	10	100	0,4-0,5	HC02.1	100 unit(s)
6,5	12	75	0,4-0,5	HC03.1	100 unit(s)
9	12	100	0,4-0,5	HC04.1	100 unit(s)
16	14	130	0,4-0,5	HC05.1	100 unit(s)
20	16	130	0,5-0,6	HC06.1	100 unit(s)
25	16	160	0,5-0,6	HC07.1	100 unit(s)
35	18	180	0,5-0,6	HC08.1	100 unit(s)
45	20	180	0,5-0,6	HC10.1	100 unit(s)
60	25	150	0,6-0,7	HC11.1	50 unit(s)
80	25	200	0,6-0,7	HC12.1	50 unit(s)
110	30	200	0,7-0,8	HC13.1	50 unit(s)



Test tubes

With straight edge, thin-walled

Material: Borosilicate glass 3.3.

Vol. (ml)	Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Art. No.	Pack Qty.	
4	10	75	0,8-1,0	HA73.1	100 unit(s)	
5	10	100	0,8-1,0	HA74.1	100 unit(s)	
6	12	75	0,8-1,0	HA75.1	100 unit(s)	
8	12	100	0,8-1,0	HA76.1	100 unit(s)	
16	14	130	1,0-1,2	HA77.1	100 unit(s)	
17	16	130	1,0-1,2	HA78.1	100 unit(s)	
21	16	160	1,0-1,2	HA79.1	100 unit(s)	
32	18	180	1,0-1,2	HA80.1	100 unit(s)	
34	20	150	1,0-1,2	HA81.1	100 unit(s)	
40	20	180	1,0-1,2	HA82.1	100 unit(s)	
55	25	150	1,0-1,2	HA83.1	50 unit(s)	
70	25	200	1,0-1,2	HA84.1	50 unit(s)	
100	30	200	1,0-1,4	HA85.1	50 unit(s)	



A121°C



Test tube DURAN®

DWK Life Sciences. Material: DURAN® glass. Autoclavable: Yes.

Thick-walled, therefore particularly resistant to mechanical loads. Good temperature stability.

Not a medical device / Not an IVD product

With beaded rim

Vol. (ml)	Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Art. No.	Pack Qty.	
2	8	70	0,8-1,0	C196.1	100 unit(s)	
4	10	75	0,8-1,0	C197.1	100 unit(s)	
5	10	100	0,8-1,0	C198.1	100 unit(s)	
6	12	75	0,8-1,0	C199.1	100 unit(s)	
8	12	100	0,8-1,0	C205.1	100 unit(s)	
9	13	100	0,8-1,0	PY92.1	100 unit(s)	
16	14	130	0,8-1,0	C206.1	100 unit(s)	
17	16	130	1,0-1,2	C207.1	100 unit(s)	
21	16	160	1,0-1,2	C208.1	100 unit(s)	
32	18	180	1,0-1,2	C209.1	100 unit(s)	
34	20	150	1,0-1,2	C210.1	100 unit(s)	
40	20	180	1,0-1,2	C211.1	100 unit(s)	
55	25	150	1,0-1,2	C212.1	50 unit(s)	
70	25	200	1,0-1,2	C213.1	50 unit(s)	
100	30	200	1,0-1,4	C214.1	50 unit(s)	

With straight edge

Vol. (ml)	Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Art. No.	Pack Qty.	
2	8	70	0,8-1,0	C182.1	100 unit(s)	
4	10	75	0,8-1,0	C183.1	100 unit(s)	
5	10	100	0,8-1,0	C184.1	100 unit(s)	
6	12	75	0,8-1,0	C185.1	100 unit(s)	
8	12	100	0,8-1,0	C186.1	100 unit(s)	
9	13	100	0,8-1,0	PY93.1	100 unit(s)	
16	14	130	0,8-1,0	C187.1	100 unit(s)	
17	16	130	1,0-1,2	C188.1	100 unit(s)	
21	16	160	1,0-1,2	C189.1	100 unit(s)	
32	18	180	1,0-1,2	C190.1	100 unit(s)	
34	20	150	1,0-1,2	C191.1	100 unit(s)	
40	20	180	1,0-1,2	C192.1	100 unit(s)	
55	25	150	1,0-1,2	C193.1	50 unit(s)	
70	25	200	1,0-1,2	C194.1	50 unit(s)	
100	30	200	1,0-1,4	C195.1	50 unit(s)	

Test tube

Mariefeld. Material: Soda-lime glass.

With round bottom and straight edge.

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Art. No.	Pack Qty.	
4	10	75	0,6	TX18.1	1,000 unit(s)	
6	12	75	0,6	TX19.1	500 unit(s)	
9	12	100	0,6	TX20.1	500 unit(s)	
14	16	100	0,8	TX21.1	310 unit(s)	
20	16	130	0,8	TX22.1	100 unit(s)	
23	16	150	0,8	TX23.2	100 unit(s)	
25	16	160	0,8	TX24.2	100 unit(s)	
35	18	180	1,0	TX25.1	100 unit(s)	
75	25	200	1,2	TX26.1	63 unit(s)	



Test tube ROTILABO®

ROTH SELECTION. Material: Soda-lime glass.

With round bottom and smooth rim, slightly narrowed at the top.

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Art. No.	Pack Qty.	
20	16	160	0,6	0486.2	250 unit(s)	

Test tubes



A 121 °C



Test tube Fiolax®

DWK Life Sciences. Material: Borosilicate glass 4.9. Autoclavable: Yes. With beaded rim. Thin-walled and relatively resistant to sudden temperature changes and localised heating.

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Art. No.	Pack Qty.
2	8	70	0,4-0,5	K223.1	100 unit(s)
4	10	75	0,4-0,5	K224.1	100 unit(s)
6	10	100	0,4-0,5	K225.1	100 unit(s)
6.5	12	75	0,4-0,5	K226.1	100 unit(s)
9	12	100	0,4-0,5	K227.1	100 unit(s)
16	14	130	0,4-0,5	K228.1	100 unit(s)
20	16	130	0,5-0,6	K229.1	100 unit(s)
25	16	160	0,4-0,5	K230.1	100 unit(s)
35	18	180	0,4-0,5	K231.1	100 unit(s)
39	20	150	0,4-0,5	K232.1	100 unit(s)
45	20	180	0,4-0,5	K233.1	100 unit(s)
60	25	150	0,6-0,7	K234.1	50 unit(s)
80	25	200	0,6-0,7	K235.1	50 unit(s)
110	30	200	0,7-0,8	K236.1	50 unit(s)



Test tube with flat bottom

DWK Life Sciences. Material: Soda-lime glass.

Colourless, without rim. For use as reaction vials for chemical reactions and for storing small amounts of liquids.

Not a medical device / Not an IVD product

Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Art. No.	Pack Qty.
11	43	0,7-0,8	1H3A.1	200 unit(s)
15	65	0,7-0,8	1H3H.1	300 unit(s)
16	40	0,7-0,8	1H39.1	200 unit(s)
16	50	0,9-1,0	1H3C.1	100 unit(s)
19	50	0,7-0,8	1H3E.1	100 unit(s)
19	70	0,7-0,8	1H3K.1	100 unit(s)
24	100	0,9-1,0	1H3L.1	50 unit(s)
30	100	1,1-1,2	1H3N.1	81 unit(s)



A 121 °C



Test tube with ground glass stopper

BRAND. Material: Borosilicate glass 3.3, stopper made of PP. Autoclavable: Yes.

With standard ground sleeve and PP stopper. White graduated.

Not a medical device / Not an IVD product

Vol. (ml)	Graduations (ml)	Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Standard	Art. No.	Pack Qty.
10	0.1	15	165	1,2	12/21	Y752.1	1 unit(s)
20	0.2	17	200	1,2	14/23	Y753.1	1 unit(s)
25	0.2	17	220	1,2	14/23	Y754.1	1 unit(s)
30	0.5	19	220	1,2	14/23	Y755.1	1 unit(s)



Test tubes



Test tube with ground glass stopper

Bohemia. Material: Borosilicate glass 3.3. Autoclavable: Yes.

With ground glass sleeve, ground glass stopper and white labelling space.

Not a medical device / Not an IVD product

Vol. (ml)	Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Standard ground glass joint	Art. No.	Pack Qty.	
10	16	150	1.8	12/21	NY90.1	1 unit(s)	
20	18	180	1.8	14/23	NY91.1	1 unit(s)	
30	22	150	1.8	19/26	NY92.1	1 unit(s)	
50	22	200	1.8	19/26	NY93.1	1 unit(s)	
100	34	200	2.0	29/32	NY94.1	1 unit(s)	



Test tubes

Kartell. Material: PS.

In vitro diagnostic medical device according to IVDD

Version	Vol. (ml)	Ø ext. (mm)	Height (mm)	Wall thickness (mm)	Art. No.	Pack Qty.	
Round base	1	11.5	55	1	K940.1	1,000 unit(s)	
Round base	2	11	70	1	K939.1	1,000 unit(s)	
Round base	4	12	75	1	K937.1	1,000 unit(s)	
Conical	10	16	105	1	K935.1	1,000 unit(s)	
Round base	10	16	100	1	K936.2	500 unit(s)	
Round base	15	16	150	1	K938.2	850 unit(s)	



Test tubes with rim

Kartell.

In vitro diagnostic medical device according to IVDD

Version	Vol. (ml)	Ø ext. (mm)	Height (mm)	Art. No.	Pack Qty.	
Conical, material: PS	10	16	100	TT74.1	2,000 unit(s)	
Cylindrical material: PP	10	16	100	TT73.1	2,000 unit(s)	



Test tube clips ROTILABO®

ROTH SELECTION. Material: Wood.

Length 180 mm.

Version	Art. No.	Pack Qty.	
For test tubes Ø 12-20 mm	0645.1	12 unit(s)	



Test tube holder ROTILABO® with finger grip

ROTH SELECTION. Material: Steel nickel-plated.

Length 130 mm.

Version	Art. No.	Pack Qty.	
For test tubes up to Ø 30 mm	X114.1	1 unit(s)	

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47