

Пинцеты CARL ROTH ROTILABO, Wafer, Aspirette

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47

Tweezers



A 121 °C

Hooked tweezers

ROTH SELECTION. Material: Stainless steel X15Cr13 (1.4024). Autoclavable: Yes.

Surgical. With teeth 1:2 for a secure hold on tissue.

Tweezers for grasping and preparing tissue samples and vessels.

Medical device class I MDD

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	teeth 1:2	115	2860.1	1 unit(s)	
straight	teeth 1:2	130	2861.1	1 unit(s)	
straight	teeth 1:2	145	2862.1	1 unit(s)	



A 121 °C

Tweezers anatomical

Hammacher. Material: Stainless steel Remanit 4024. Autoclavable: Yes. With fine tips, grooved on the inside.

Tweezers for gentle collection of tissue samples.

Medical device class I MDD

Adson fine

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	blunt	120	YK36.1	1 unit(s)	
straight	blunt	150	YK37.1	1 unit(s)	

Micro Adson

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	blunt	120	YK38.1	1 unit(s)	

①



②



A 121 °C

Tweezers special shape (splinter tweezers)

Hammacher. Material: Stainless steel Remanit 4024. Autoclavable: Yes.

Tweezers for removing very small foreign bodies from the skin or from wounds.

Medical device class I MDD

Pic.	Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
(1)	straight	pointed	80	YK48.1	1 unit(s)	
(2)	straight	pointed	105	YK49.1	1 unit(s)	

①



②



A 121 °C

Tweezers with guide pin

ROTH SELECTION. Material: Stainless steel 18/10 (Remanit 4301). Autoclavable: Yes.

Fine type. Inside serrated tips.

Not a medical device / Not an IVD product

(1) Straight

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	pointed	105	LT35.1	1 unit(s)	
straight	pointed	130	LT36.1	1 unit(s)	

(2) Curved

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
curved	pointed	105	LT37.1	1 unit(s)	
curved	pointed	130	LT38.1	1 unit(s)	



Tweezers



A 121 °C

Tweezers with guide pin

Hammacher. Material: Stainless steel Wironit. Autoclavable: Yes.

Straight, tips ribbed on the inside

Not a medical device / Not an IVD product

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	pointed	155	HCL3.1	1 unit(s)	

Angled, tips ribbed on the inside

Tweezers for gentle collection of tissue samples.

Medical device class I MDD


Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
angled	pointed	155	HCL2.1	1 unit(s)	

Angled, tips smooth

Tweezers for gentle collection of tissue samples.

Medical device class I MDD

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
angled	pointed	155	HCL1.1	1 unit(s)	



A 121 °C

Tweezers with insulated tips

Hammacher. Material: Stainless steel Wironit. Autoclavable: Yes.

Blunt tips, 45° angled, for secure gripping of small parts.

Not a medical device / Not an IVD product

Pic.	Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
(1)	angled	blunt	165	ETC8.1	1 unit(s)	
(2)	angled	blunt	300	ETC9.1	1 unit(s)	



A 121 °C

Tweezers

ROTH SELECTION. Material: Stainless steel 18/10 (Remanit 4301). Autoclavable: Yes.

Inside of tip is fluted for secure grip.

Not a medical device / Not an IVD product

(1) Straight blunt

Anatomical.

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	blunt	105	2687.1	1 unit(s)	
straight	blunt	115	2688.1	1 unit(s)	
straight	blunt	130	2689.1	1 unit(s)	
straight	blunt	145	2690.1	1 unit(s)	
straight	blunt	160	2691.1	1 unit(s)	
straight	blunt	200	2692.1	1 unit(s)	
straight	blunt	250	2693.1	1 unit(s)	

(2) Straight pointed

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	pointed	105	2695.1	1 unit(s)	
straight	pointed	115	2801.1	1 unit(s)	
straight	pointed	130	2802.1	1 unit(s)	
straight	pointed	145	2851.1	1 unit(s)	
straight	pointed	160	2852.1	1 unit(s)	
straight	pointed	200	2853.1	1 unit(s)	

(3) Curved pointed

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
curved	pointed	105	2854.1	1 unit(s)	
curved	pointed	115	2855.1	1 unit(s)	
curved	pointed	130	2856.1	1 unit(s)	
curved	pointed	145	2857.1	1 unit(s)	
curved	pointed	160	2858.1	1 unit(s)	
curved	pointed	200	2859.1	1 unit(s)	

Curved blunt

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
curved	blunt	115	PX28.1	1 unit(s)	
curved	blunt	130	PX29.1	1 unit(s)	
curved	blunt	145	PX30.1	1 unit(s)	
curved	blunt	160	PX31.1	1 unit(s)	
curved	blunt	200	PX32.1	1 unit(s)	

Tweezers



Tweezer set ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10 (Remanit® 4301), Case made of plastic.

Contents:

- 1 tweezers with fine tip, curved, L 120 mm
- 1 tweezers straight, blunt, 130 mm
- 1 tweezers with guide pin, straight, 130 mm
- 1 double-glass tweezers, curved, 115 mm

Not a medical device / Not an IVD product

Length (mm)	Width (mm)	Art. No.	Pack Qty.
150	90	PX40.1	1 unit(s)



A 121 °C

Forceps long

ROTH SELECTION. Material: Stainless steel 18/10 (Remanit 4301). Autoclavable: Yes.

Blunt. Ribbed on the inside for a secure hold.

Not a medical device / Not an IVD product

Straight

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.
straight	blunt	300	0951.1	1 unit(s)

Curved

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.
curved	blunt	300	HHA0.1	1 unit(s)



Tweezers round

ROTH SELECTION. Material: Stainless steel.

With rounded tips, inside grooved for secure hold.

Not a medical device / Not an IVD product

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.
straight	blunt	145	HEC7.1	1 unit(s)



A 121 °C

Tweezers self-tensioning

Usbeck. Material: Stainless steel. Autoclavable: Yes.

Not a medical device / Not an IVD product

Angled

Pic.	Version	Tip shape	Length (mm)	Art. No.	Pack Qty.
(1)	angled	flat	105	HL63.1	1 unit(s)

Straight, pointed

Pic.	Version	Tip shape	Length (mm)	Art. No.	Pack Qty.
(2)	straight	pointed	120	HL64.1	1 unit(s)



A 121 °C

Fixation forceps Graefe straight

Hammacher. Material: Stainless steel Wironit (1.4301).

Autoclavable: Yes.

With ratchet for a secure hold.

Universal gripping tweezers for grasping tissues.

Not a medical device / Not an IVD product

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.
straight	flat	110	KAP3.1	1 unit(s)



Tweezers



A 121 °C

Cover slip tweezers

ROTH SELECTION. Material: Stainless steel 18/10 (Remanit 4301). Autoclavable: Yes.

According to Kühne.

Not a medical device / Not an IVD product

Straight

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	flat	115	K717.1	1 unit(s)	

Curved

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
curved	flat	115	K718.1	1 unit(s)	
curved	flat	130	K719.1	1 unit(s)	
curved	flat	145	K720.1	1 unit(s)	



A 121 °C

Cover slip tweezers

Hammacher. Material: Stainless steel Wironit. Autoclavable: Yes.

Not a medical device / Not an IVD product

(1) Straight

Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
straight	flat	5	145	EXC6.1	1 unit(s)	

(2) Curved

Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
curved	flat	5	120	EXC7.1	1 unit(s)	



A 121 °C

Special tweezers Wafer

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes.

Not a medical device / Not an IVD product

PTFE-coated

Wafer tweezers with two stop angles and narrow upper part.

Pic.	Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
(1)	straight	flat	7	120	LX39.1	1 unit(s)	

3 teeth

Pic.	Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
(2)	straight	flat	9	125	LX40.1	1 unit(s)	

4 teeth

Pic.	Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
(3)	straight	flat	12	125	LX41.1	1 unit(s)	



Special tweezers for filter paper

Bochem. Material: Stainless steel 18/10.

For picking up filter papers, cover slips and other flat objects.

Not a medical device / Not an IVD product

Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
straight	flat	12	125	C828.1	1 unit(s)	

Tweezers



A 121 °C

Special tweezers with ring-shaped pick-up

ROTH SELECTION. Material: Stainless steel 18/10. Autoclavable: Yes.
Tip dimensions: W 6.5 x Ø 3.5 mm.

Not a medical device / Not an IVD product

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	ring-shaped	120	LX48.1	1 unit(s)	



A 121 °C

Special tweezers ring-shaped

Hammacher. Material: Stainless steel 18/10 (Remanit 4301).
Autoclavable: Yes.

For delicate tissue and specimens.

Not a medical device / Not an IVD product

Version	Tip shape	Length (mm)	Ø ext. (mm)	Ø Internal (mm)	Art. No.	Pack Qty.	
straight	ring-shaped	90	3	2.2	LL00.1	1 unit(s)	
straight	ring-shaped	90	6	4.8	LL01.1	1 unit(s)	



Spring steel tweezers

ROTH SELECTION. Material: Stainless steel 1.4301.
Lightweight, thin and flexible. Very gentle on delicate items, such as insects.

Not a medical device / Not an IVD product

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	pointed	105	AL99.1	1 unit(s)	
straight	blunt	115	AN00.1	1 unit(s)	



Weighing tweezers ROTILABO®

Usbeck. Material: Stainless steel, PVC tips.

Not a medical device / Not an IVD product

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
curved	blunt	90	YT63.1	1 unit(s)	



A 121 °C

Tweezers PTFE-coated

ROTH SELECTION. Material: Stainless steel 18/10 (1.4301), PTFE-coated. Autoclavable: Yes.

Suitable for liquids such as acids, bases, solvents, as long as coating is intact.

- Non-stick
- Temperature stability: -200 to +260 °C

Not a medical device / Not an IVD product

Straight pointed

Version	Tip shape	Colour	Length (mm)	Art. No.	Pack Qty.	
straight	pointed	black	115	C833.2	1 unit(s)	
straight	pointed	black	130	C834.2	1 unit(s)	
straight	pointed	black	145	C835.2	1 unit(s)	

Straight blunt

Version	Tip shape	Colour	Length (mm)	Art. No.	Pack Qty.	
straight	blunt	black	115	C830.1	1 unit(s)	
straight	blunt	black	130	C831.1	1 unit(s)	
straight	blunt	black	145	C832.1	1 unit(s)	

Curved pointed

Version	Tip shape	Colour	Length (mm)	Art. No.	Pack Qty.	
curved	pointed	black	115	C836.2	1 unit(s)	
curved	pointed	black	130	C837.2	1 unit(s)	
curved	pointed	black	145	C838.2	1 unit(s)	



Tweezers



Tweezer set ROTILABO® PTFE-coated

ROTH SELECTION. Material: Stainless steel 18/10 (material 1.4301) PTFE-coated, case made of plastic.

Red PTFE coating with good anti-adhesive properties forms a separating layer between the tweezers material and the sample. Dull colour of the coating enables glare-free working under the bright lighting in microscopy. The chemical-resistant coating also protects against contact with the nickel contained in the stainless steel.

Contents:

- 1 tweezers with fine tip, L 120 mm
- 1 pointed tweezers with guide pin, L 130 mm
- 1 pointed tweezers, curved, L 130 mm
- 1 blunt tweezers, L 145 mm

Not a medical device / Not an IVD product

Length (mm)	Width (mm)	Art. No.	Pack Qty.
150	90	ENX5.1	1 unit(s)



Special tweezers aspirette®

Hirschmann.

With vacuum pump. Suitable for picking up and positioning small and micro parts weighing up to 70 g. Can be used to position membrane filters, cover slips, microscope slides, optical lenses, electronic chips and much more.

- With single-button operation
- Ergonomically designed
- Battery can be charged approx. 1000 times

Delivery incl. 230 V charging unit, three lifting needles and three suction cups, Ø 3.2/6 and 10 mm.

Version	Art. No.	Pack Qty.
Battery powered vacuum tweezer	A247.1	1 unit(s)



Tweezers long nylon

ROTH SELECTION. Material: Nylon, Handle made of stainless steel 18/10.

Temperature resistant from -20 to +150 °C.

Not a medical device / Not an IVD product

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.
straight	blunt	300	PT95.1	1 unit(s)



Tweezers ceramic

ROTH SELECTION. Material: Ceramic.

- With non-stick surface
- Low abrasion
- Antimagnetic
- Low susceptibility to electrostatic charging
- Resists acids and bases (except for HF)
- Temperature stability: to +700 °C

Not a medical device / Not an IVD product

(1) Straight

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.
straight	pointed	120	5327.1	1 unit(s)


(2) Curved

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.
curved	pointed	125	5328.1	1 unit(s)

Tweezers

Tweezers ceramic

ROTH SELECTION. Material: Ceramic, PAA handle.
Hardened. Metal-free and antimagnetic.



Not a medical device / Not an IVD product

Pic.	Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
(1)	straight	pointed	0.5	125	PE65.1	1 unit(s)	
(2)	straight	pointed	1	125	PE66.1	1 unit(s)	
(3)	straight	blunt	2	125	PE67.1	1 unit(s)	

Tweezers ROTILABO® disposable

ROTH SELECTION.
Anatomical form.



(1) ABS, sterile

Material: ABS. **Sterile**. Intended use: For handling clinical cotton wool, gauze, plasters and other items in contact with human skin. Suitable for human samples. For use by trained laboratory employees only.
NOT suitable for taking samples directly from the human body.

In vitro diagnostic medical device according to IVDD

Version	Tip shape	Colour	Tip width (mm)	Length (mm)	Pack.	Art. No.	Pack Qty.	
straight	blunt	white	4	120	100 x 1	KL05.1	100 unit(s)	

Tweezers disposable

ROTH SELECTION. Material: ABS. Non-sterile.



Not a medical device / Not an IVD product

Version	Tip shape	Colour	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
straight	blunt	blue	3.7	111	1L90.1	6 unit(s)	



Tweezers

Tweezers disposable

Bürkle. Material: PS.

- Lightweight plastic tweezers for fatigue-free work
- Secure handling thanks to grooved grips
- Manufactured in a clean room (clean room class 7)
- Compliant with the EU food directives and FDA guidelines



Standard

Non-sterile.

Version	Tip shape	Colour	Length (mm)	Pack.	Art. No.	Pack Qty.	
angled	pointed	white	130	1 x 100	1A54.1	100 unit(s)	
straight	pointed	white	130	1 x 100	1A53.1	100 unit(s)	
straight	flat	white	130	1 x 100	1A55.1	100 unit(s)	

SteriPlast® sterile

Sterile.

Version	Tip shape	Colour	Length (mm)	Pack.	Art. No.	Pack Qty.	
angled	pointed	white	130	100 x 1	1A57.1	100 unit(s)	
straight	pointed	white	130	100 x 1	1A56.1	100 unit(s)	
straight	flat	white	130	100 x 1	1A58.1	100 unit(s)	

SteriPlast® can be detected, sterile

Sterile. Can be detected using a colour scanner, metal detector or X-rays.

Version	Tip shape	Colour	Length (mm)	Pack.	Art. No.	Pack Qty.	
angled	pointed	blue	130	100 x 1	1A5A.1	100 unit(s)	
straight	pointed	blue	130	100 x 1	1A59.1	100 unit(s)	
straight	flat	blue	130	100 x 1	1A5C.1	100 unit(s)	

Tweezers plastic

ROTH SELECTION. Material: POM (glass fibre reinforced), with carbon fibres.

- Good chemical resistance
- Anti-static
- High strength due to glass fibre reinforcement
- Temperature stability: to +90 °C



Not a medical device / Not an IVD product

Straight flat tip

Pic.	Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
(1)	straight	flat	3.0	110	C317.1	3 unit(s)	

Straight pointed

Pic.	Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
(2)	straight	pointed	0.3	110	C318.1	3 unit(s)	
(3)	straight	pointed	1.0	110	C319.1	3 unit(s)	

Curved flat tip

Pic.	Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
(4)	curved	flat	2.0	110	C320.1	3 unit(s)	

Tweezers

Tweezers plastic

ROTH SELECTION. Material: CTFE.

High chemical resistance, acid-resistance. Suitable for use in chemical laboratories and in electronics. Temperature stability: to +200 °C.



Not a medical device / Not an IVD product

Ribbed tips

Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
straight	blunt	3	112	C322.1	1 unit(s)	

Tips smooth

Version	Tip shape	Tip width (mm)	Length (mm)	Art. No.	Pack Qty.	
straight	blunt	3	112	C321.1	1 unit(s)	

Tweezers plastic

Vitlab. Material: POM.

Flexible, with rounded tips. Very good resilience. With stria on the outside for optimal handling and grip. Tweezer TT37.1 additionally with stria on the inside of the tip.



Not a medical device / Not an IVD product

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	blunt	118	TT35.1	5 unit(s)	
straight	blunt	147	TT36.1	5 unit(s)	
straight	blunt	253	TT37.1	5 unit(s)	

Tweezers plastic PTFE

ROTH SELECTION. Material: PTFE. Autoclavable: Yes.

- Flexible, with smooth tips
- High chemical resistance
- Temperature stability: -200 to +250 °C



A 121 °C

Not a medical device / Not an IVD product

Fine tips

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	blunt	100	1A0C.1	1 unit(s)	
straight	blunt	150	0504.1	1 unit(s)	
straight	blunt	200	1A0L.1	1 unit(s)	

Flat tips

Version	Tip shape	Length (mm)	Art. No.	Pack Qty.	
straight	flat	100	1A0E.1	1 unit(s)	
straight	flat	150	1A0H.1	1 unit(s)	
straight	flat	200	1A0K.1	1 unit(s)	



Tweezers



Tweezers long wood

ROTH SELECTION. Material: Wood (Beech, untreated).
Suitable for removing samples from liquid nitrogen.

Not a medical device / Not an IVD product

Version	Tip shape	Length (mm)	Width (mm)	Art. No.	Pack Qty.	
straight	flat	300	20	EET0.1	1 unit(s)	
straight	flat	400	20	EET1.1	1 unit(s)	



Tweezers stands ROTILABO®

ROTH SELECTION. Material: Stainless steel 18/10 (1.4301), ESD-coated.

For safe storage of sensitive tweezers as well as for drying tweezers after cleaning. Tweezers holder can be screwed upside down under a shelf for space-saving storage.

Outer length (mm)	Outer width (mm)	Outer height (mm)	Art. No.	Pack Qty.	
185	60	135	PX20.1	1 unit(s)	



A 121 °C

Tweezers stands

Hammacher. Material: Stainless steel, Wironit (1.4301).
Autoclavable: Yes.

For safe storage of sensitive tweezers as well as for drying tweezers after cleaning.

Outer length (mm)	Outer width (mm)	Outer height (mm)	Art. No.	Pack Qty.	
140	120	140	KAP8.1	1 unit(s)	



A 121 °C

Scalpel solid

ROTH SELECTION. Material: Stainless steel. Autoclavable: Yes.

Not a medical device / Not an IVD product

Chromium steel 1.4034

Material: 1.4034.

Pic.	Type	Length (mm)	Art. No.	Pack Qty.	
(1)	anatomic	137	PX21.1	1 unit(s)	

Chromium steel 1.4021

Material: 1.4021.

Pic.	Type	Length (mm)	Art. No.	Pack Qty.	
(2)	pointed	125	PX23.1	1 unit(s)	
(3)	spherical	130	PX22.1	1 unit(s)	



A 121 °C

Scalpel handles (hinge insert)

Hammacher. Material: Stainless steel Remanit 4021. Autoclavable: Yes.

Handle for holding and guiding scalpel blades. Hinge mechanism for inserting and removing the blades reduces risk of injury. Suitable for scalpel blades types 10 to 16 (CE20.1ff and 3608.1ff).

Medical device class I MDD

Version	Art. No.	Pack Qty.	
No. 3	HCY7.1	1 unit(s)	

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Россия (495)268-04-70

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Казахстан (772)734-952-31

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Киргизия (996)312-96-26-47

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93