

Бутыли промывочные CARL ROTH ROTILABO, GL 18, GL 25, GL 32, GL 40, GL 45, GL 52, GL 56, GL 63

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47

Wash bottles



A 121 °C



Wash bottle ROTILABO®

ROTH SELECTION. Material: Borosilicate glass 3.3. Autoclavable: Yes. Spray bottle consisting of flat-bottomed flask with sleeve NS 29/32 and spray bottle attachment.

Delivery without suction aid, please order separately.

Vol. (ml)	Standard ground glass joint	Art. No.	Pack Qty.	
250	29/32	ACK8.1	1 unit(s)	
500	29/32	ACK9.1	1 unit(s)	
1,000	29/32	ACL0.1	1 unit(s)	



Wash bottle ROTILABO® wide neck

ROTH SELECTION. Material: LDPE.

Comes with trilingual contents label and hazard identification symbol to GHS. Spray cap, colour-coded. Empties completely due to shiftable riser.

Note: The chemical resistance of the spray bottles is limited when aggressive filling media are used. The bottles should therefore be replaced regularly and not used for analytical or critical applications.

Version	Vol. (ml)	Art. No.	Pack Qty.	
Acetone	250	T398.1	1 unit(s)	
Acetone	500	T399.1	1 unit(s)	
Ethanol	250	T400.1	1 unit(s)	
Ethanol	500	T401.1	1 unit(s)	
Methanol	250	T404.1	1 unit(s)	
Methanol	500	T405.1	1 unit(s)	
Isopropanol	250	T402.1	1 unit(s)	
Isopropanol	500	T403.1	1 unit(s)	
Distilled water	250	T406.1	1 unit(s)	
Distilled water	500	T407.1	1 unit(s)	



Wash bottle ROTILABO® with overpressure valve

ROTH SELECTION.

Material: LDPE (bottle types Acetone and MEK made of PP), screw cap and spray pipe made of PP.

Overpressure valve in the screw cap extensively prevents leakage, dripping or escape of the liquid from the spray pie, e.g. in the case of temperature fluctuations.

With permanent four-language content imprint (German, English, French, Spanish) and hazard warning acc. to GHS.

Note: The chemical resistance of the spray bottles is limited when aggressive filling media are used. The bottles should therefore be replaced regularly and not used for analytical or critical applications.

Volume 500 ml

Version	Vol. (ml)	Thread	Art. No.	Pack Qty.	
Dist. water	500	GL 25	HL52.1	1 unit(s)	
Acetone	500	GL 25	HL51.1	1 unit(s)	
Ethanol	500	GL 25	HL54.1	1 unit(s)	
Isopropanol	500	GL 25	HL55.1	1 unit(s)	
Methyl ethyl ketone (MEK)	500	GL 25	PY44.1	1 unit(s)	

Volume 1000 ml

Version	Vol. (ml)	Thread	Art. No.	Pack Qty.	
Dist. water	1,000	GL 32	HL58.1	1 unit(s)	
Acetone	1,000	GL 32	HL57.1	1 unit(s)	
Ethanol	1,000	GL 32	HL59.1	1 unit(s)	
Isopropanol	1,000	GL 32	HL61.1	1 unit(s)	



Wash bottle ROTILABO®

ROTH SELECTION. Material: LDPE transparent.

Height data without riser.

Vol. (ml)	Ø (mm)	Height (mm)	Art. No.	Pack Qty.	
250	60	145	5330.1	1 unit(s)	
500	75	180	5331.1	1 unit(s)	
1,000	100	230	5332.1	1 unit(s)	

Spray bottles set ROTILABO®

ROTH SELECTION.

Consisting of 3 wash bottles, volumes 250, 500 and 1000 ml.

Version	Art. No.	Pack Qty.	
Spray bottle set ROTILABO®	5333.1	1 set	



Wash bottles

Wash bottle

Vitelab. Material: LDPE, transparent, screw cap and spray tube with elongated tip made of PP.



Narrow neck

Vol. (ml)	Ø (mm)	Height (mm)	Thread	Art. No.	Pack Qty.	
50	37	73	GL 18	XY89.1	1 unit(s)	
100	43	95	GL 18	XY90.1	1 unit(s)	
250	58	135	GL 25	XY91.1	1 unit(s)	
500	74	180	GL 25	XY93.1	1 unit(s)	
1,000	92	221	GL 32	XY95.1	1 unit(s)	

Wide neck

Vol. (ml)	Ø (mm)	Height (mm)	Thread	Art. No.	Pack Qty.	
250	58	146	GL 45	XY92.1	1 unit(s)	
500	76	166	GL 45	XY94.1	1 unit(s)	
1,000	91	226	GL 63	XY96.1	1 unit(s)	

Screw closures

Vitelab. Autoclavable: Yes.

With integrated lip seal, leak-proof.



A121 °C

Polypropylene

Material: PP.

Thread	Art. No.	Pack Qty.	
GL 18	YA32.1	24 unit(s)	
GL 25	YA33.1	12 unit(s)	
GL 32	YA34.1	12 unit(s)	
GL 40	YA35.1	12 unit(s)	
GL 45	YA36.1	12 unit(s)	
GL 52	YA37.1	12 unit(s)	
GL 56	YA38.1	12 unit(s)	
GL 63	YA39.1	12 unit(s)	

Perfluoralkoxy

Material: PFA.

Thread	Art. No.	Pack Qty.	
GL 18	YA40.1	1 unit(s)	
GL 25	YA41.1	1 unit(s)	
S 28	YA42.1	1 unit(s)	
S 40	YA43.1	1 unit(s)	

Wash bottles

Wash bottle ROTILABO®

ROTH SELECTION. Material: LDPE, screw cap and spray tube made of PP.

Precision jet and optimised medium backflow through fine tip of spray tube. Easy identification thanks to selection of differently coloured bottled. Height data without riser.



Volume 250 ml

Vol. (ml)	Colour	Ø (mm)	Height (mm)	Thread	Art. No.	Pack Qty.	
250	red	58	135	GL 25	APK1.1	1 unit(s)	
250	green	58	135	GL 25	APK2.1	1 unit(s)	
250	yellow	58	135	GL 25	APK3.1	1 unit(s)	
250	blue	58	135	GL 25	APK4.1	1 unit(s)	

Volume 500 ml

Vol. (ml)	Colour	Ø (mm)	Height (mm)	Thread	Art. No.	Pack Qty.	
500	red	75	180	GL 25	6530.1	1 unit(s)	
500	green	75	180	GL 25	6533.1	1 unit(s)	
500	yellow	75	180	GL 25	6536.1	1 unit(s)	
500	blue	75	180	GL 25	6539.1	1 unit(s)	

Spray bottles set ROTILABO®

ROTH SELECTION.

Consisting of 4 wash bottles volume 500 ml in red, green, yellow and blue.

Version	Art. No.	Pack Qty.	
Spray bottle set ROTILABO®	APK5.1	1 set	

Wash bottle ROTILABO®

ROTH SELECTION. Material: PE.

With integrated stable riser and removable sealing cap.



Vol. (ml)	Length (mm)	Width (mm)	Height (mm)	Art. No.	Pack Qty.	
250	90	60	140	Y801.1	1 unit(s)	
500	96	65	178	Y802.1	1 unit(s)	



Wash bottles



Wash bottle

Kartell. Material: PE.

With curved riser and removable spray nozzle.

Vol. (ml)	Ø (mm)	Art. No.	Pack Qty.	
50	35	Y767.1	1 unit(s)	
100	43	Y768.1	1 unit(s)	
250	60	Y769.1	1 unit(s)	
500	75	Y770.1	1 unit(s)	
1,000	95	Y771.1	1 unit(s)	



Wash bottle wide neck

Nalgene. Material: LDPE, screw cap made of PP.

With wide-neck opening for fast and easy filling. Bottle and riser are made from one piece, which allows the bottle to be emptied completely. With leak-proof cap.

Type 2402

Vol. (ml)	Colour	Thread	Art. No.	Pack Qty.	
125	colourless	24 mm	C464.1	6 unit(s)	
250	colourless	38 mm	C465.1	4 unit(s)	
500	colourless	38 mm	C466.1	4 unit(s)	
750	colourless	38 mm	C467.1	4 unit(s)	
1,000	colourless	43 mm	C468.1	2 unit(s)	

Type DS 2408

Vol. (ml)	Colour	Thread	Art. No.	Pack Qty.	
250	red	38 mm	K669.1	1 unit(s)	
500	red	38 mm	K670.1	1 unit(s)	



Wash bottle

Nalgene. Material: LDPE transparent, screw cap and riser made of coloured PP.

Type 2422. With coloured caps for easy identification of the contents.

Leak-proof, since cap and riser are made from one piece.

Vol. (ml)	Colour	Art. No.	Pack Qty.	
500	yellow	K665.1	1 unit(s)	
500	blue	K666.1	1 unit(s)	
500	green	K667.1	1 unit(s)	
500	red	K668.1	1 unit(s)	



A121 °C

Wash bottle made of fluoroplastic

Nalgene. Material: FEP, riser and screw cap made of ETFE.

Autoclavable: Yes.

Type 2403. High chemical resistance. For contamination-free working with trace element analyses. **Leak-proof**, since riser holder and cap are made from one piece. Phthalat-free.

Vol. (ml)	Art. No.	Pack Qty.	
125	K671.1	1 unit(s)	
250	K672.1	1 unit(s)	
500	K673.1	1 unit(s)	
1,000	K674.1	1 unit(s)	

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Россия (495)268-04-70

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Казахстан (772)734-952-31

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Киргизия (996)312-96-26-47

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93