

Бани водяные лабораторные CARL ROTH

Технические характеристики

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тургут (3462)77-98-35
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия (495)268-04-70

Казахстан (772)734-952-31

Киргизия (996)312-96-26-47



Water baths



Water bath WBT series

For a working temperature range from 25 to 100 °C.

- Stainless steel inner tank
- Microprocessor-controlled temperature controller
- Digital temperature display
- Adjustable upper and lower setpoint value limit
- Visual and/or audible limit alarm (adjustable)
- With dry run protection
- Memory function saves three combinations of setpoint temperature and time parameters
- Bath is emptied via a floor drain with a locking screw. A drain cock is also available as an option (please order separately)
- WBT 1 water bath supplied with round, flat cover, all other models have a slanting cover
- Heater exposed on floor

Technical specifications:

Art. No.	EA02.1	EA03.1	EA04.1	EA05.1
Type	WBT1	WBT6	WBT12	WBT22
Bath volume	1–2 l	2–6 l	4–12 l	6–22 l
Inner length		215 mm	390 mm	455 mm
Inner width		230 mm		
Inner height	120 mm			170 mm
Inner Ø x H	130 x 120 mm			
Outer length	230 mm	320 mm	490 mm	570 mm
Outer width	240 mm	340 mm		
Outer height	220 mm	260 mm		310 mm
Operating temperature	+25 to +100 °C			
Temperature constancy	±0,25 K			
Heating output	700 W	1000 W	1600 W	2000 W
Weight	5.0 kg	7.5 kg	10.5 kg	12.5 kg
Mains connection	230 V, 50/60 Hz			

Delivery incl. removable bottom.

Round bath opening

Type	Bath volume (l)	Art. No.	Pack Qty.	
WBT1	1–2	EA02.1	1 unit(s)	

Square bath opening

Type	Bath volume (l)	Art. No.	Pack Qty.	
WBT6	2–6	EA03.1	1 unit(s)	
WBT12	4–12	EA04.1	1 unit(s)	
WBT22	6–22	EA05.1	1 unit(s)	

Accessories for WBT series water baths

Round test tube rack for WBT 1 water bath

Made of stainless steel. Water bath holding capacity: 1 piece

Suitable for	Ø Hole (mm)	Art. No.	Pack Qty.	
36 test tubes Ø 10.5 mm	10.5	EA07.1	1 unit(s)	
36 test tubes Ø 11.5 mm	11.5	EA08.1	1 unit(s)	
24 test tubes Ø 14.5 mm	14.5	EA09.1	1 unit(s)	
16 test tubes Ø 18.5 mm	18.5	EA10.1	1 unit(s)	
6 test tubes Ø 32 mm	32.0	EA11.1	1 unit(s)	

Test tube racks for water baths WBT 6 to WBT 22

Made of stainless steel. Water bath holding capacity:
WBT 6 = 3 pieces, WBT 12 = 5 pieces, WBT 22 = 6 pieces.

Suitable for	Ø Hole (mm)	Art. No.	Pack Qty.	
56 test tubes Ø 10.5 mm	10.5	EA12.1	1 unit(s)	
56 test tubes Ø 11.5 mm	11.5	EA13.1	1 unit(s)	
33 test tubes Ø 14.5 mm	14.5	EA14.1	1 unit(s)	
18 test tubes Ø 18.5 mm	18.5	EA15.1	1 unit(s)	
6 test tubes Ø 32 mm	32.0	EA16.1	1 unit(s)	

Drain cock for water baths WBT 1 to WBT 22

Designation	Art. No.	Pack Qty.	
Drain cock, complete	EA17.1	1 unit(s)	

Water baths

Water bath Hydro-series

Lauda.

For an operating temperature range from 25 to 100 °C.

- TFT display showing all process data
- Adjustment and display of temperature values at increments of 0.1 °C
- With timer function, max. 999 hrs, 59 mins
- Adjustable trigger points for overheating and insufficient temperature alarm up to max. ± 10 K
- Visual and audible alarm in the event of lack of water, overheating, insufficient temperature and gauge break or short circuit
- Three application-specific timer functions
- Bath unit, cover frame, sieve plate, lid and heater made of stainless steel
- With a drain valve



Technical specifications:

Art. No.	1A81.1	1A82.1	1A87.1	1A83.1	1A88.1	1A84.1	1A85.1	1A86.1
Type	H 4	H 8	H 8 A	H 16	H 16 A	H 22	H 24	H 41
Bath volume	4 l	8 l		16 l		22 l	24 l	41 l
Inner length	245 mm			400 mm			600 mm	410 mm
Inner width	100 mm	200 mm		245 mm				296 mm
Inner height	115 mm					180 mm	115 mm	285 mm
Outer length	340 mm			500 mm			700 mm	510 mm
Outer width	290 mm	395 mm		440 mm				490 mm
Outer height	275 mm		345 mm	275 mm	345 mm		275 mm	476 mm
Operating temperature	25 to 100 °C							
Temperature constancy	0,1 K							
Heating output	500 W	1000 W		1500 W				
Mains connection	230 V, 50/60 Hz							

Delivery incl. stainless steel lid.

Standard version

Type	Bath volume (l)	Art. No.	Pack Qty.	
H 4	4	1A81.1	1 unit(s)	
H 8	8	1A82.1	1 unit(s)	
H 16	16	1A83.1	1 unit(s)	
H 22	22	1A84.1	1 unit(s)	
H 24	24	1A85.1	1 unit(s)	
H 41	41	1A86.1	1 unit(s)	

Design with a circulation system

Type	Bath volume (l)	Art. No.	Pack Qty.	
H 8 A	8	1A87.1	1 unit(s)	
H 16 A	16	1A88.1	1 unit(s)	



Water baths

Water bath WTB series

Memmert.

For temperature control of ointments, emulsions, samples, plates and nutrient solutions in the laboratory, as well as warm storage and calibration in industrial applications.

- For an operating temperature range from 5 °C above ambient temperature to +100 °C
- Intuitive 3.5" colour graphic display with touchscreen
- Integrated digital clock from 1 min to 99.59 h
- Independent temperature limiter for heater switch-off in the event of a fault
- Corrosion-protected heating outside of the working area
- Stainless steel tray with no corners or edges for ease of cleaning, incl. central drain (no obstructing fittings)
- Bath drainage by drain valve



With slanting cover

Technical specifications:

Art. No.	1EYA.1	1EYC.1	1EYE.1	1EYH.1	1EYK.1	1EYL.1
Type	WTB 6	WTB 11	WTB 15	WTB 24	WTB 35	WTB 50
Bath volume	7 l	10 l	17 l	23 l	37 l	51 l
Inner length	277 mm	267 mm	482 mm	472 mm	602 mm	592 mm
Inner width	217 mm	207 mm	277 mm	267 mm	482 mm	472 mm
Inner height	150 mm	200 mm	150 mm	200 mm	150 mm	200 mm
Outer length	349 mm		554 mm		674 mm	
Outer width	289 mm		349 mm		554 mm	
Outer height	245 mm	295 mm	245 mm	295 mm	245 mm	295 mm
Height with slanting cover	348 mm	398 mm	348 mm	398 mm	348 mm	398 mm
Operating temperature	From 5 °C above room temperature to 100 °C					
Set temperature range	+10 to +100 °C					
Setting accuracy	0,1 °C					
Temporal temp. variance	±0,1 K					
Spatial temp. deviation	±0,3 K					
Heating output	1000 W		1420 W		2800 W	
Weight	10 kg	12 kg	14 kg	16 kg	20 kg	23 kg
Mains connection	230 V, 50/60 Hz					

Delivery incl. stainless steel slanting cover and drain valve.

Type	Bath volume (l)	Art. No.	Pack Qty.	
WTB 6	7	1EYA.1	1 unit(s)	
WTB 11	10	1EYC.1	1 unit(s)	
WTB 15	17	1EYE.1	1 unit(s)	
WTB 24	23	1EYH.1	1 unit(s)	
WTB 35	37	1EYK.1	1 unit(s)	
WTB 50	51	1EYL.1	1 unit(s)	

Water baths

With flat cover and sets of rings

Technical specifications:

Art. No.	1EXX.1	1EXY.1	1EY0.1	1EY1.1	1EY2.1	1EY3.1
Type	WTB 6	WTB 11	WTB 15	WTB 24	WTB 35	WTB 50
Bath volume	7 l	10 l	17 l	23 l	37 l	51 l
Inner length	277 mm	267 mm	482 mm	472 mm	602 mm	592 mm
Inner width	217 mm	207 mm	277 mm	267 mm	482 mm	472 mm
Inner height	150 mm	200 mm	150 mm	200 mm	150 mm	200 mm
Outer length	349 mm		554 mm		674 mm	
Outer width	289 mm		349 mm		554 mm	
Outer height	245 mm	295 mm	245 mm	295 mm	245 mm	295 mm
Height with flat cover	245 mm	295 mm	245 mm	295 mm	245 mm	295 mm
Ring set openings	3		6		8	
Ø Flat cover opening	87 mm		107 mm			
Operating temperature	From 5 °C above room temperature to 100 °C					
Set temperature range	+10 to +100 °C					
Setting accuracy	0,1 °C					
Temporal temp. variance	±0,1 K					
Spatial temp. deviation	±0,3 K					
Heating output	1000 W		1420 W		2800 W	
Weight	10 kg	12 kg	14 kg	16 kg	20 kg	23 kg
Mains connection	230 V, 50/60 Hz					


Delivery incl. flat cover with sets of rings and drain valve.

Type	Bath volume (l)	Art. No.	Pack Qty.	
WTB 6	7	1EXX.1	1 unit(s)	
WTB 11	10	1EXY.1	1 unit(s)	
WTB 15	17	1EY0.1	1 unit(s)	
WTB 24	23	1EY1.1	1 unit(s)	
WTB 35	37	1EY2.1	1 unit(s)	
WTB 50	51	1EY3.1	1 unit(s)	

Accessories floor grids for water baths WTB series

Memmert.

Usable on both sides. Height: 30 mm or 60 mm.



Suitable for	Art. No.	Pack Qty.	
WTB 6/11	1EYN.1	1 unit(s)	
WTB 15/24	1EYP.1	1 unit(s)	
WTB 35/50	1EYT.1	1 unit(s)	



Water baths

Water bath Precision series

Thermo.

- **For a working temperature range of 30 to 100 °C** (except for the GP 02 model which can only be used up to 90 °C)
- Microprocessor-controlled temperature controller
- Digital temperature display
- Integrated timer for programming on/off times, timer function (1 to 99 mins) or continuous operation-
- Memory function can save up to four setpoint temperatures
- Equipped with protection against overheating
- Heating elements under the bath tank (no fittings in bath tank)
- Internal water drainage tube with integrated quick-release coupling (only with models GP 10, GP 20, GP 28 and GP 15 Dual)



Technical specifications:

Art. No.	HAC0.1	HAC1.1	HAC2.1	HAC3.1	HAC4.1	HAC5.1
Type	GP 02	GP 05	GP 10	GP 20	GP 28	GP 15 Dual
Bath volume	2 l	5 l	10 l	20 l	28 l	5/10 l
Inner length	155 mm	300 mm	330 mm	500 mm		
Inner width	138 mm	154 mm	301 mm	297 mm		
Inner height	150 mm				200 mm	
Outer length	199 mm	355 mm	383 mm	555 mm		587 mm
Outer width	230 mm	246 mm	393 mm	392 mm		
Outer height	344 mm	363 mm			412 mm	363 mm
Operating temperature	30 to 90 °C	30 to 100 °C				
Temperature constancy	±0,2 K					
Heating output	200 W	300 W	800 W	1200 W		300/800 W
Mains connection	100–155 / 200–230 V, 50/60 Hz					

Delivery incl. transparent cover.

Standard

Pic.	Type	Bath volume (l)	Art. No.	Pack Qty.
	GP 02	2	HAC0.1	1 unit(s)
(1)	GP 05	5	HAC1.1	1 unit(s)
	GP 10	10	HAC2.1	1 unit(s)
(2)	GP 20	20	HAC3.1	1 unit(s)
	GP 28	28	HAC4.1	1 unit(s)

Dual bath

Two independently controllable baths in one device. A different temperature can be set for each bath.

Bath 1: Volume 5 l; L x W x H internal: 300 x 154 x 150 mm

Bath 2: Volume 10 l; L x W x H internal: 330 x 301 x 150 mm

Pic.	Type	Bath volume (l)	Art. No.	Pack Qty.
(3)	GP 15 Dual	5/10	HAC5.1	1 unit(s)

Water baths

Water bath PURA™ PURA™ series

Julabo.

For a temperature range up to 99.9 °C.

- Digital temperature display
- Internal timer function (max. 99 h, 59 mins in 1-min increments), with signal after expiry of the time, optionally with or without heating switch-off
- Integrated drain screw (except PURA™ 4)
- Easy to clean as there are no fittings in the bath
- Integrated run-dry protection and recessed grips
- **Foldable bath cover available as accessory, please order separately**



Technical specifications:

Art. No.	HCH1.1	HCH2.1	HCH3.1	HCH4.1	HCH5.1
Type	Pura™ 4	Pura™ 10	Pura™ 14	Pura™ 22	Pura™ 30
Bath volume	0,8–4,5 l	1,4–9,5 l	2–14 l	3,4–25,5 l	4,8–36 l
Inner length	120 mm	220 mm	330 mm	550 mm	770 mm
Inner width	270 mm				
Inner height	140 mm			180 mm	
Outer length	210 mm	310 mm	420 mm	640 mm	860 mm
Outer width	350 mm				
Outer height	220 mm			260 mm	
Height with lid	360 mm			400 mm	
Display resolution	0,1 °C				
Temperature constancy	0,15 K				
Heating output	500 W	1200 W	1800 W	2000 W	
Weight	4,9 kg	7 kg	8,5 kg	11,5 kg	14,5 kg
Mains connection	230 V, 50/60 Hz				

Pic.	Type	Bath volume (l)	Art. No.	Pack Qty.	
(1)	Pura™ 4	0,8–4,5	HCH1.1	1 unit(s)	
(2)	Pura™ 10	1,4–9,5	HCH2.1	1 unit(s)	
	Pura™ 14	2–14	HCH3.1	1 unit(s)	
(3)	Pura™ 22	3,4–25,5	HCH4.1	1 unit(s)	
	Pura™ 30	4,8–36	HCH5.1	1 unit(s)	

Accessories bath cover folds back, transparent for PURA™ series

Julabo.

Suitable for	Art. No.	Pack Qty.	
for Pura™ 4	HCH6.1	1 unit(s)	
for Pura™ 10	HCH7.1	1 unit(s)	
for Pura™ 14	HCH8.1	1 unit(s)	
for Pura™ 22	HCH9.1	1 unit(s)	
for Pura™ 30	HCK0.1	1 unit(s)	



Water baths

Water bath with circulation pump WTB series

Memmert.

With pump device for continuous circulation of the temperature control medium.

- For an operating temperature range from 5 °C above ambient temperature to +100 °C
- Pump can be switched on and off
- Intuitive 3.5" colour graphic display with touchscreen
- Integrated digital clock from 1 min to 99.59 h
- Independent temperature limiter for heater switch-off in the event of a fault
- Corrosion-protected heating outside of the working area
- Stainless steel tray with no corners or edges for ease of cleaning, incl. central drain (no obstructing fittings)
- Bath drainage by drain valve



Technical specifications:

Art. No.	1EY4.1	1EY5.1	1EY6.1	1EY7.1	1EY8.1	1EY9.1
Type	WTB 6	WTB 11	WTB 15	WTB 24	WTB 35	WTB 50
Bath volume	7 l	10 l	17 l	23 l	37 l	51 l
Inner length	277 mm	267 mm	482 mm	472 mm	602 mm	592 mm
Inner width	217 mm	207 mm	277 mm	267 mm	482 mm	472 mm
Inner height	150 mm	200 mm	150 mm	200 mm	150 mm	200 mm
Outer length	349 mm		554 mm		674 mm	
Outer width	289 mm		349 mm		554 mm	
Outer height	245 mm	295 mm	245 mm	295 mm	245 mm	295 mm
Height with slanting cover	348 mm	398 mm	348 mm	398 mm	348 mm	398 mm
Operating temperature	From 5 °C above room temperature to 100 °C					
Set temperature range	+10 to +100 °C					
Setting accuracy	0,1 °C					
Temporal temp. variance	±0,1 K					
Spatial temp. deviation	±0,3 K					
Heating output	1000 W		1420 W		2800 W	
Weight	10 kg	12 kg	14 kg	16 kg	20 kg	23 kg
Mains connection	230 V, 50/60 Hz					


Delivery incl. stainless steel slanting cover and pump device.

Type	Bath volume (l)	Art. No.	Pack Qty.
WTB 6	7	1EY4.1	1 unit(s)
WTB 11	10	1EY5.1	1 unit(s)
WTB 15	17	1EY6.1	1 unit(s)
WTB 24	23	1EY7.1	1 unit(s)
WTB 35	37	1EY8.1	1 unit(s)
WTB 50	51	1EY9.1	1 unit(s)

Accessories floor grids for water baths WTB series

Memmert.

Usable on both sides. Height: 30 mm or 60 mm.



Suitable for	Art. No.	Pack Qty.
WTB 6/11	1EYN.1	1 unit(s)
WTB 15/24	1EYP.1	1 unit(s)
WTB 35/50	1EYT.1	1 unit(s)

Water baths



Heating bath with magnetic drive, HBR 4 control

IKA.

Effective volume 4 l. Interior free of heating coils (no components in the bath tub). Made of stainless steel. Available with low viscosity oil (50 mPas) or water as a heat conductor. Digital display for setpoint, actual and safety temperature as well as the speed. **Integrated magnetic drive for rotating the heating medium**, creating a better heat distribution in the bath. Heating element under the bath tank (no fittings in bath tank). With RS 232 port. Two handles.

Technical specifications:

Art. No.	LN86.1
Type	HBR 4
Bath volume	4 l
Speed range	100 to 800 rpm
Inner height	160 mm
Inner Ø x H	200 x 160 mm
Outer height	250 mm
Operating temperature	Room temperature up to +200 °C
Setting accuracy	±1 K
Control deviation	±1 K
Heating output	1000 W
Weight	4.4 kg
Mains connection	230 V, 50/60 Hz

Type	Bath volume (l)	Art. No.	Pack Qty.
HBR 4	4	LN86.1	1 unit(s)



Accessories set of rings for HBR 4 control heating bath

IKA.

For covering, comprises 10 rings with a Ø of 25 to 185 mm. Prevents dust gathering, heat radiation, intake of water and formation of oil mist when working with oil.

Designation	Art. No.	Pack Qty.
Set of rings	9700.1	1 unit(s)



Accessories intermediate bottom HBR 4 control

IKA.

Made of stainless steel. For achieving a better heat distribution without disrupting the movement of the rotating stirring bar.

Type	Art. No.	Pack Qty.
Intermediate bottom	9689.1	1 unit(s)



Serial water baths Hydro H-V series

Lauda.

For an operating temperature range from 5 °C above room temperature to 100 °C. Special baths for gentle exhaust steam applications from columns, Erlenmeyer flasks or beakers. For unsupervised continuous operation based on water level controllers and a low-water safety device. Bath unit, cover frame with openings and split sieve plate over the tubular heating element, made of stainless steel. Bath clearance 150 mm.

Technical specifications:

Art. No.	1H1P.1	1H1N.1
Type	Hydro H 9 V	Hydro H 6 V
Bath volume	9 l	5,3 l
Inner height	90 mm	
Outer length	982 mm	682 mm
Outer width	232 mm	
Outer height	190 mm	
Ring set openings	6	4
Ø Flat cover opening	131 mm	
Operating temperature	25 to 100 °C	
Temperature constancy	±3 K	
Heating output	1500 W	1000 W
Weight	20 kg	15 kg
Mains connection	230 V, 50/60 Hz	

Delivery incl. removable flat/perforated lids with multi-part set of rings and stand rods (L 600 x Ø 12 mm).

Type	Bath volume (l)	Art. No.	Pack Qty.
Hydro H 9 V	9	1H1P.1	1 unit(s)
Hydro H 6 V	5,3	1H1N.1	1 unit(s)

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Россия (495)268-04-70

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Казахстан (772)734-952-31

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Сургут (3462)77-98-35
Тверь (4822)63-31-35

Киргизия (996)312-96-26-47

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93